



PERCENT SOLID

Supervisor: Iwona
Analyst: RUBINA
Date: 2/10/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 13:15
In Date: 02/09/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 108
Time OUT: 08:50
Out Date: 02/10/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB124002

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O1420-01	01A-01B-01C	1	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O1420-02	02A-02B-02C	2	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O1432-01	DOMINICK-STREET-UTLIITY-RELOCATION-EXCAVATION	3	1.15	8.78	9.93	9.06	90.1	
O1432-02	DOMINICK-STREET-UTLIITY-RELOCATION-EXCAVATION	4	1.13	8.47	9.6	8.81	90.7	
O1470-01	NWB-1823	5	1.13	8.51	9.64	9.22	95.1	
O1485-01	EA-FILTERBAGS-0208	6	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1487-01	340-341	7	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O1487-02	342-343	8	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O1489-01	BARN-EAST-END-DRAIN	9	1.15	8.84	9.99	6.82	64.1	
O1490-01	CLI-2023-1A	10	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1490-02	CLI-2023-1B	11	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1490-03	KZZ169K-1-1	12	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1490-04	KZZ169K-2-1	13	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1490-05	KMA614E-1-1	14	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1490-06	KMA614E-2-1	15	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1491-01	2742	16	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O1502-01	KMA652M-1-1	17	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1502-02	KMA652M-1-2	18	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1502-03	KMA652M-2-1	19	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1502-04	KMA652M-2-2	20	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1502-05	KMA891I-1-1	21	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1502-06	KMA891I-1-2	22	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1502-07	JEH655I-1-1	23	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1502-08	JEH655I-1-2	24	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1502-09	BC248186-1-1	25	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1502-10	BC248186-1-2	26	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: Iwona
Analyst: RUBINA
Date: 2/10/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 13:15
In Date: 02/09/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 108
Time OUT: 08:50
Out Date: 02/10/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB124002

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O1503-01	D2880	27	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1503-02	D2895	28	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1503-03	D2896	29	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1504-01	COMP-1	30	1.15	8.74	9.89	9.38	94.2	
O1504-03	EPH-1	31	1.15	8.74	9.89	9.38	94.2	
O1504-04	EPH-2	32	1.15	8.74	9.89	9.38	94.2	
O1507-01	ETGI-330	33	1.14	8.76	9.9	9.31	93.3	
O1507-02	ETGI-330-EPH-2	34	1.14	8.76	9.9	9.31	93.3	
O1508-01	RS-1	35	1.14	8.49	9.63	9.31	96.2	
O1508-02	UST-1	36	1.14	8.85	9.99	9.28	92.0	
O1508-03	UST-2	37	1.15	8.52	9.67	9.29	95.5	
O1508-04	UST-3	38	1.15	8.68	9.83	9.35	94.5	
O1508-05	UST-4	39	1.14	8.57	9.71	8.81	89.5	
O1508-06	RS-2	40	1.15	8.81	9.96	8.66	85.2	
O1508-07	PA-1	41	1.14	8.78	9.92	9.02	89.7	
O1508-08	PA-2	42	1.15	8.84	9.99	9.5	94.5	
O1508-09	RS-3	43	1.12	8.87	9.99	9.2	91.1	
O1508-10	CV-1	44	1.15	8.83	9.98	8.65	84.9	
O1510-01	TWM-TP14-020823	45	1.14	8.75	9.89	9.15	91.5	
O1510-03	TWM-TP14-020823	46	1.15	8.65	9.8	8.82	88.7	
O1510-05	JPP-51B-020823	47	1.15	8.61	9.76	7.46	73.3	
O1510-07	JPP-51B-020823	48	1.14	8.65	9.79	8.38	83.7	
O1510-09	JPP-55-020823	49	1.13	8.86	9.99	8.84	87.0	
O1510-11	JPP-55-020823	50	1.15	8.84	9.99	8.73	85.7	
O1510-13	JPP-59-020823	51	1.14	8.81	9.95	9.58	95.8	
O1510-15	JPP-59-020823	52	1.14	8.85	9.99	9.39	93.2	
O1513-01	VNJ-208	53	1.13	8.84	9.97	8.91	88.0	
O1513-02	VNJ-208	54	1.16	8.72	9.88	8.83	88.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$