



PERCENT SOLID

Supervisor: Iwona

Analyst: Al-Terek

Date: 2/22/2023

OVENTEMP IN Celsius(°C): 108

Time IN: 12:30

In Date: 02/16/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 106

Time OUT: 09:30

Out Date: 02/17/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB124122

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O1534-01	C0B39	1	1.11	8.83	9.94	8.11	79.3	
O1561-02	BUR-1527	2	1.00	1.00	2.00	2.00	100.0	
O1575-01	NWB-1887	3	1.16	8.81	9.97	9.72	97.2	
O1576-01	RE137-BAGFILTER-20220215	4	1.00	1.00	2.00	2.00	100.0	
O1579-01	186	5	1.00	1.00	2.00	2.00	100.0	
O1581-01	MH-A	6	1.14	8.66	9.8	8.42	84.1	
O1581-03	MH-A-EPH-2	7	1.13	8.82	9.95	8.72	86.1	
O1581-04	MH-C	8	1.14	8.84	9.98	8.62	84.6	
O1581-06	MH-C-EPH-2	9	1.12	8.77	9.89	8.58	85.1	
O1582-01	LEONIA-TRENCH-LINE	10	1.13	8.78	9.91	9.33	93.4	
O1582-03	LEONIA-TRENCH-LINE-EPH-2	11	1.13	8.78	9.91	9.33	93.4	
O1582-04	LEONIA-BARRIER	12	1.13	8.78	9.91	9.33	93.4	
O1583-01	KMH303C-1-1	13	1.00	1.00	2.00	2.00	100.0	
O1583-02	KMH303C-1-2	14	1.00	1.00	2.00	2.00	100.0	
O1584-01	IRV2023-9A	15	1.00	1.00	2.00	2.00	100.0	
O1584-02	IRV2023-9B	16	1.00	1.00	2.00	2.00	100.0	
O1585-01	EDA810Q-1-1	17	1.00	1.00	2.00	2.00	100.0	
O1585-02	EDA810Q-1-2	18	1.00	1.00	2.00	2.00	100.0	
O1586-01	COMP-1	19	1.13	8.76	9.89	8.9	88.7	
O1586-03	COMP-2	20	1.11	8.86	9.97	8.16	79.6	
O1586-05	COMP-3	21	1.12	8.84	9.96	8.52	83.7	
O1586-07	COMP-4	22	1.13	8.84	9.97	8.99	88.9	
O1586-09	COMP-5	23	1.12	8.82	9.94	9.06	90.0	
O1587-01	OK-02-02162023	24	1.12	8.82	9.94	9.14	90.9	
O1587-02	OK-02-02162023-EPH-2	25	1.12	8.82	9.94	9.14	90.9	
O1589-01	MIXED-DEMO	26	1.00	1.00	2.00	2.00	100.0	
O1590-01	BC260967LOZ-1-1	27	1.00	1.00	2.00	2.00	100.0	



PERCENT SOLID

Supervisor: Iwona
Analyst: Al-Terek
Date: 2/22/2023

OVENTEMP IN Celsius(°C): 108
Time IN: 12:30
In Date: 02/16/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 106
Time OUT: 09:30
Out Date: 02/17/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB124122

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O1590-02	BC260967LOZ-1-2	28	1.00	1.00	2.00	2.00	100.0	
O1590-03	BC271315LOZ-1-1	29	1.00	1.00	2.00	2.00	100.0	
O1590-04	BC271315LOZ-1-2	30	1.00	1.00	2.00	2.00	100.0	
O1590-05	BC271242LOZ-1-1	31	1.00	1.00	2.00	2.00	100.0	
O1590-06	BC271242LOZ-1-2	32	1.00	1.00	2.00	2.00	100.0	
O1592-01	BC260952-1-1	33	1.00	1.00	2.00	2.00	100.0	
O1592-02	BC260952-1-2	34	1.00	1.00	2.00	2.00	100.0	
O1592-03	BC226467-1-1	35	1.00	1.00	2.00	2.00	100.0	
O1592-04	BC226467-1-2	36	1.00	1.00	2.00	2.00	100.0	
O1592-05	BC254741-1-1	37	1.00	1.00	2.00	2.00	100.0	
O1592-06	BC254741-1-2	38	1.00	1.00	2.00	2.00	100.0	
O1592-07	HEF499K-1-1	39	1.00	1.00	2.00	2.00	100.0	
O1592-08	HEF499K-1-2	40	1.00	1.00	2.00	2.00	100.0	
O1592-09	HEF500T-1-1	41	1.00	1.00	2.00	2.00	100.0	
O1592-10	HEF500T-1-2	42	1.00	1.00	2.00	2.00	100.0	
O1592-11	KMA772A-1-1	43	1.00	1.00	2.00	2.00	100.0	
O1592-12	KMA772A-1-2	44	1.00	1.00	2.00	2.00	100.0	
O1592-13	BC260897-1-1	45	1.00	1.00	2.00	2.00	100.0	
O1592-14	BC260897-1-2	46	1.00	1.00	2.00	2.00	100.0	
O1593-01	B-13-SB01	47	1.13	8.69	9.82	8.71	87.2	
O1593-02	B-13-SB03	48	1.11	8.81	9.92	8.67	85.8	
O1593-03	B-26-SB01	49	1.12	8.86	9.98	8.66	85.1	
O1593-04	B-26-SB03	50	1.12	8.68	9.8	8.55	85.6	
O1593-05	SED-2	51	1.11	8.87	9.98	5.29	47.1	
O1593-08	SED-3	52	1.12	8.83	9.95	5.21	46.3	
O1595-01	S-322-1	53	1.12	8.82	9.94	9.77	98.1	
O1595-02	S-322-2	54	1.14	8.67	9.81	9.11	91.9	
O1595-03	S-322-3	55	1.12	8.79	9.91	9.18	91.7	



PERCENT SOLID

Supervisor: Iwona
Analyst: Al-Terek
Date: 2/22/2023

OVENTEMP IN Celsius(°C): 108
Time IN: 12:30
In Date: 02/16/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 106
Time OUT: 09:30
Out Date: 02/17/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB124122

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O1595-04	S-655-1	56	1.11	8.77	9.88	9.41	94.6	
O1595-05	S-655-2	57	1.12	8.79	9.91	9.08	90.6	
O1595-06	S-655-3	58	1.1	8.83	9.93	9.06	90.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$