

Clear table


| Instrument ID: IC-2 |                      |                       |                       |                       |                       |                        |                       |             |               |                          | Analyst : | IZ | Method: | 300.0 / 9056A |  |
|---------------------|----------------------|-----------------------|-----------------------|-----------------------|-----------------------|------------------------|-----------------------|-------------|---------------|--------------------------|-----------|----|---------|---------------|--|
| ident               | concentra<br>tion F- | concentratio<br>n CL- | concentrati<br>on NO2 | concentrati<br>on BR- | concentrati<br>on NO3 | concentrati<br>on HPO4 | concentrati<br>on SO4 | file name   | date time     | Initial wt/<br>Final Vol | Analyst   |    |         |               |  |
| STD1                | 0.0000               | 0.0000                | 0.0000                | 0.0000                | 0.0000                | 0.0000                 | 0.0000                | 2023-02-21_ | 2/21/23 10:44 |                          | IZ/MA     |    |         |               |  |
| STD2                | 0.4920               | 0.6780                | 0.7170                | 2.3360                | 0.5860                | 1.2450                 | 3.5390                | 2023-02-21_ | 2/21/23 11:16 |                          | IZ/MA     |    |         |               |  |
| STD3                | 0.8220               | 1.2480                | 1.2550                | 4.1190                | 1.0260                | 2.0990                 | 6.2150                | 2023-02-21_ | 2/21/23 11:48 |                          | IZ/MA     |    |         |               |  |
| STD4                | 0.9980               | 1.4880                | 1.4540                | 4.9470                | 1.2430                | 2.2960                 | 7.4400                | 2023-02-21_ | 2/21/23 13:37 |                          | IZ/MA     |    |         |               |  |
| STD5                | 1.8430               | 2.8780                | 2.8510                | 9.5750                | 2.3920                | 4.8400                 | 14.2760               | 2023-02-21_ | 2/21/23 14:09 |                          | IZ/MA     |    |         |               |  |
| STD6                | 3.8900               | 5.8880                | 5.9030                | 19.5800               | 4.8850                | 9.8330                 | 29.2840               | 2023-02-21_ | 2/21/23 14:40 |                          | IZ/MA     |    |         |               |  |
| STD7                | 5.1240               | 7.6200                | 7.6200                | 25.4430               | 6.3680                | 12.6870                | 38.2450               | 2023-02-21_ | 2/21/23 15:12 |                          | IZ/MA     |    |         |               |  |
| ICV                 | 1.8340               | 2.9800                | 3.0180                | 9.9890                | 2.4670                | 4.8900                 | 14.7310               | 2023-02-21_ | 2/21/23 17:02 |                          | IZ/MA     |    |         |               |  |
| ICB                 | 0.0000               | 0.0000                | 0.0000                | 0.0000                | 0.0000                | 0.0000                 | 0.0000                | 2023-02-21_ | 2/21/23 18:03 |                          | IZ/MA     |    |         |               |  |
| CCV                 | 1.9810               | 3.1480                | 3.0690                | 10.4550               | 2.5620                | 4.9560                 | 15.3390               | 2023-02-21_ | 2/21/23 10:28 |                          | IZ/MA     |    |         |               |  |
| CCB                 | 0.0000               | 0.0000                | 0.0000                | 0.0000                | 0.0000                | 0.0000                 | 0.0000                | 2023-02-22_ | 2/22/23 11:01 |                          | IZ/MA     |    |         |               |  |
| LB124209BLW         | 0.0000               | 0.0000                | 0.0000                | 0.0000                | 0.0000                | 0.0000                 | 0.0000                | 2023-02-22_ | 2/22/23 11:32 |                          | IZ/MA     |    |         |               |  |
| LB124209BSW         | 2.0060               | 3.1580                | 3.0990                | 10.5660               | 2.6080                | 5.1780                 | 15.6290               | 2023-02-22_ | 2/22/23 12:03 |                          | IZ/MA     |    |         |               |  |
| O1634-01            | 0.0000               | 0.0000                | 0.0000                | 0.0000                | 0.0000                | 0.0000                 | 0.0000                | 2023-02-22_ | 2/22/23 12:34 |                          | IZ/MA     |    |         |               |  |
| O1634-01MS          | 1.7970               | 2.8540                | 2.8230                | 9.5590                | 2.4170                | 4.6520                 | 14.2060               | 2023-02-22_ | 2/22/23 13:04 |                          | IZ/MA     |    |         |               |  |
| O1634-01MSD         | 1.7660               | 2.8230                | 2.8220                | 9.4160                | 2.3750                | 4.5090                 | 13.9490               | 2023-02-22_ | 2/22/23 13:35 |                          | IZ/MA     |    |         |               |  |
| O1646-01            | 0.0000               | 0.0000                | 0.0000                | 0.0000                | 0.2010                | 0.0000                 | 0.0000                | 2023-02-22_ | 2/22/23 14:05 |                          | IZ/MA     |    |         |               |  |
| O1646-02            | 0.0000               | 0.0000                | 0.0000                | 0.0000                | 0.1520                | 0.0000                 | 0.0000                | 2023-02-22_ | 2/22/23 14:36 |                          | IZ/MA     |    |         |               |  |
| O1646-03            | 0.0000               | 0.0000                | 0.0000                | 0.0000                | 0.1690                | 0.0000                 | 0.0000                | 2023-02-22_ | 2/22/23 15:07 |                          | IZ/MA     |    |         |               |  |
| CCV                 | 1.9720               | 3.1810                | 3.1440                | 10.6010               | 2.6360                | 5.2700                 | 15.6110               | 2023-02-22_ | 2/22/23 16:10 |                          | IZ/MA     |    |         |               |  |
| CCB                 | 0.0000               | 0.0000                | 0.0000                | 0.0000                | 0.0000                | 0.0000                 | 0.0000                | 2023-02-22_ | 2/22/23 16:41 |                          | IZ/MA     |    |         |               |  |
| O1627-08            | 0.4020               | 0.5160                | 0.5790                | 1.9470                | 0.4600                | 0.9440                 | 3.6690                | 2023-02-22_ | 2/22/23 17:12 |                          | IZ/MA     |    |         |               |  |
| O1627-07            | 0.2290               | 0.2180                | 0.3310                | 1.0180                | 0.2250                | 0.6320                 | 1.8060                | 2023-02-22_ | 2/22/23 17:42 |                          | IZ/MA     |    |         |               |  |
| O1627-07RE 0.25M    | 0.1730               | 0.0880                | 0.2160                | 0.5840                | 0.1100                | 0.4920                 | 1.3250                | 2023-02-22_ | 2/22/23 18:13 |                          | IZ/MA     |    |         |               |  |
| CCV                 | 2.0180               | 3.1530                | 3.0640                | 10.4120               | 2.5900                | 5.3920                 | 15.4310               | 2023-02-22_ | 2/22/23 18:44 |                          | IZ/MA     |    |         |               |  |
| CCB                 | 0.0000               | 0.0000                | 0.0000                | 0.0000                | 0.0000                | 0.0000                 | 0.0000                | 2023-02-22_ | 2/22/23 19:14 |                          | IZ/MA     |    |         |               |  |

Clear table

| Instrument ID: IC-2 |                      |                       |                       |                       |                       |                        |                       |             |               |                          |         |
|---------------------|----------------------|-----------------------|-----------------------|-----------------------|-----------------------|------------------------|-----------------------|-------------|---------------|--------------------------|---------|
| Analyst :           |                      |                       |                       | IZ                    |                       | Method: 300.0 / 9056A  |                       |             |               |                          |         |
| ident               | concentra<br>tion F- | concentratio<br>n CL- | concentrati<br>on NO2 | concentrati<br>on BR- | concentrati<br>on NO3 | concentrati<br>on HPO4 | concentrati<br>on SO4 | file name   | date time     | Initial wt/<br>Final Vol | Analyst |
| STD1                | 0.0000               | 0.0000                | 0.0000                | 0.0000                | 0.0000                | 0.0000                 | 0.0000                | 2023-02-21_ | 2/21/23 10:44 |                          | IZ/MA   |
| STD2                | 0.4920               | 0.6780                | 0.7170                | 2.3360                | 0.5860                | 1.2450                 | 3.5390                | 2023-02-21_ | 2/21/23 11:16 |                          | IZ/MA   |
| STD3                | 0.8220               | 1.2480                | 1.2550                | 4.1190                | 1.0260                | 2.0990                 | 6.2150                | 2023-02-21_ | 2/21/23 11:48 |                          | IZ/MA   |
| STD4                | 0.9980               | 1.4880                | 1.4540                | 4.9470                | 1.2430                | 2.2960                 | 7.4400                | 2023-02-21_ | 2/21/23 13:37 |                          | IZ/MA   |
| STD5                | 1.8430               | 2.8780                | 2.8510                | 9.5750                | 2.3920                | 4.8400                 | 14.2760               | 2023-02-21_ | 2/21/23 14:09 |                          | IZ/MA   |
| STD6                | 3.8900               | 5.8880                | 5.9030                | 19.5800               | 4.8850                | 9.8330                 | 29.2840               | 2023-02-21_ | 2/21/23 14:40 |                          | IZ/MA   |
| STD7                | 5.1240               | 7.6200                | 7.6200                | 25.4430               | 6.3680                | 12.6870                | 38.2450               | 2023-02-21_ | 2/21/23 15:12 |                          | IZ/MA   |
| ident               | True<br>Value F-     | True Value<br>CL-     | True Value<br>NO2     | True Value<br>BR-     | True Value<br>NO3     | True Value<br>HPO4     | True Value<br>SO4     |             |               |                          |         |
| STD1                | 0.0000               | 0.0000                | 0.0000                | 0.0000                | 0.0000                | 0.0000                 | 0.0000                |             |               |                          |         |
| STD2                | 0.4000               | 0.6000                | 0.6000                | 2.0000                | 0.5000                | 1.0000                 | 3.0000                |             |               |                          |         |
| STD3                | 0.8000               | 1.2000                | 1.2000                | 4.0000                | 1.0000                | 2.0000                 | 6.0000                |             |               |                          |         |
| STD4                | 1.0000               | 1.5000                | 1.5000                | 5.0000                | 1.2500                | 2.5000                 | 7.5000                |             |               |                          |         |
| STD5                | 2.0000               | 3.0000                | 3.0000                | 10.0000               | 2.5000                | 5.0000                 | 15.0000               |             |               |                          |         |
| STD6                | 4.0000               | 6.0000                | 6.0000                | 20.0000               | 5.0000                | 10.0000                | 30.0000               |             |               |                          |         |
| STD7                | 5.0000               | 7.5000                | 7.5000                | 25.0000               | 6.2500                | 12.5000                | 37.0000               |             |               |                          |         |
| ident               | Relative<br>Error F- | Relative<br>Error CL- | Relative<br>Error NO2 | Relative<br>Error BR- | Relative<br>Error NO3 | Relative<br>Error HPO4 | Relative<br>Error SO4 |             |               |                          |         |
| STD1                |                      |                       |                       |                       |                       |                        |                       |             |               |                          |         |
| STD2                | 23.0000              | 13.0000               | 19.5000               | 16.8000               | 17.2000               | 24.5000                | 17.9667               |             |               |                          |         |
| STD3                | 2.7500               | 4.0000                | 4.5833                | 2.9750                | 2.6000                | 4.9500                 | 3.5833                |             |               |                          |         |
| STD4                | -0.2000              | -0.8000               | -3.0667               | -1.0600               | -0.5600               | -8.1600                | -0.8000               |             |               |                          |         |
| STD5                | -7.8500              | -4.0667               | -4.9667               | -4.2500               | -4.3200               | -3.2000                | -4.8267               |             |               |                          |         |
| STD6                | -2.7500              | -1.8667               | -1.6167               | -2.1000               | -2.3000               | -1.6700                | -2.3867               |             |               |                          |         |
| STD7                | 2.4800               | 1.6000                | 1.6000                | 1.7720                | 1.8880                | 1.4960                 | 3.3649                |             |               |                          |         |

# CALIBRATION OF COMPONENT F-

Method: AnionsIC2-022123.mtw  
 Equation:  $Q = 0.0749918 \cdot A + 4.34204$   
 RSD: 5.920 %  
 Correlation coefficient: 0.998088




K3 = 0      K2 = 0      K1 = 0.0749918      K0 = 4.34204  
 Base: Area  
 Ref.channel: ch1  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height | Area  | Conc. | Vol/Dil | Retention | Used | File |
|-------|--------|-------|-------|---------|-----------|------|------|
| 1     | 0      | 0     | 0     | 0       | 0         |      | 0    |
| 2     | 1.369  | 81.6  | 0.4   | 20      | 6.113     |      |      |
| 3     | 3.148  | 161.3 | 0.8   | 20      | 6.113     |      |      |
| 4     | 4.373  | 208.2 | 1     | 20      | 6.113     |      |      |
| 5     | 8.63   | 433.6 | 2     | 20      | 6.113     |      |      |
| 6     | 20.89  | 979.5 | 4     | 20      | 6.113     |      |      |
| 7     | 27.03  | 1309  | 5     | 20      | 6.113     |      |      |

CALIBRATION OF COMPONENT CL-

Method: AnionsIC2-022123.mtw  
 Equation:  $Q = 0.126762 \cdot A + 3.40502$   
 RSD: 3.398 %  
 Correlation coefficient: 0.999370




K3 = 0      K2 = 0      K1 = 0.126762      K0 = 3.40502  
 Base: Area  
 Ref.channel: ch1  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height | Area  | Conc. | Vol/Dil | Retention | Used | File |
|-------|--------|-------|-------|---------|-----------|------|------|
| 1     | 0      | 0     | 0     | 0       | 0         |      | 0    |
| 2     | 4.987  | 80.06 | 0.6   | 20      | 9.089     |      |      |
| 3     | 10.35  | 170   | 1.2   | 20      | 9.089     |      |      |
| 4     | 12.75  | 207.9 | 1.5   | 20      | 9.089     |      |      |
| 5     | 26.5   | 427.3 | 3     | 20      | 9.089     |      |      |
| 6     | 56.44  | 902.1 | 6     | 20      | 9.089     |      |      |
| 7     | 73.41  | 1175  | 7.5   | 20      | 9.089     |      |      |

1.□

# CALIBRATION OF COMPONENT NO2

Method: AnionsIC2-022123.mtw  
 Equation:  $Q = 0.0557567 \cdot A + 5.76076$   
 RSD: 3.856 %  
 Correlation coefficient: 0.999189




K3 = 0      K2 = 0      K1 = 0.0557567      K0 = 5.76076  
 Base: Area  
 Ref.channel: ch1  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height | Area  | Conc. | Vol/Dil | Retention | Used | File |
|-------|--------|-------|-------|---------|-----------|------|------|
| 1     | 0      | 0     | 0     | 0       | 0         |      | 0    |
| 2     | 7.543  | 153.9 | 0.6   | 20      | 10.92     |      |      |
| 3     | 16.86  | 346.9 | 1.2   | 20      | 10.92     |      |      |
| 4     | 20.85  | 418.2 | 1.5   | 20      | 10.92     |      |      |
| 5     | 45.79  | 919.5 | 3     | 20      | 10.92     |      |      |
| 6     | 98.56  | 2014  | 6     | 20      | 10.92     |      |      |
| 7     | 128    | 2630  | 7.5   | 20      | 10.92     |      |      |

CALIBRATION OF COMPONENT BR-

Method: AnionsIC2-022123.mtw  
 Equation:  $Q = 0.281656 \cdot A + 14.8485$   
 RSD: 3.758 %  
 Correlation coefficient: 0.999230




K3 = 0      K2 = 0      K1 = 0.281656      K0 = 14.8485  
 Base: Area  
 Ref.channel: ch1  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height | Area  | Conc. | Vol/Dil | Retention | Used | File |
|-------|--------|-------|-------|---------|-----------|------|------|
| 1     | 0      | 0     | 0     | 0       | 0         |      | 0    |
| 2     | 5.403  | 113.2 | 2     | 20      | 14.05     |      |      |
| 3     | 11.55  | 239.8 | 4     | 20      | 14.05     |      |      |
| 4     | 14.44  | 298.5 | 5     | 20      | 14.05     |      |      |
| 5     | 30.61  | 627.2 | 10    | 20      | 14.05     |      |      |
| 6     | 65.85  | 1338  | 20    | 20      | 14.05     |      |      |
| 7     | 86.08  | 1754  | 25    | 20      | 14.05     |      |      |

1. ☐

# CALIBRATION OF COMPONENT NO3

Method: AnionsIC2-022123.mtw  
 Equation:  $Q = 0.0491063 \cdot A + 3.24923$   
 RSD: 3.939 %  
 Correlation coefficient: 0.999154




K3 = 0      K2 = 0      K1 = 0.0491063      K0 = 3.24923  
 Base: Area  
 Ref.channel: ch1  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height | Area  | Conc. | Vol/Dil | Retention | Used | File |
|-------|--------|-------|-------|---------|-----------|------|------|
| 1     | 0      | 0     | 0     | 0       | 0         |      | 0    |
| 2     | 6.522  | 172.4 | 0.5   | 20      | 15.97     |      |      |
| 3     | 13.59  | 351.7 | 1     | 20      | 15.97     |      |      |
| 4     | 17.08  | 440.1 | 1.25  | 20      | 15.97     |      |      |
| 5     | 35.76  | 908.2 | 2.5   | 20      | 15.97     |      |      |
| 6     | 77.05  | 1923  | 5     | 20      | 15.97     |      |      |
| 7     | 101.6  | 2527  | 6.25  | 20      | 15.97     |      |      |

CALIBRATION OF COMPONENT HPO4

Method: AnionsIC2-022123.mtw  
 Equation:  $Q = 0.0553071 \cdot A + 13.4923$   
 RSD: 4.068 %  
 Correlation coefficient: 0.999097



K3 = 0      K2 = 0      K1 = 0.0553071      K0 = 13.4923

Base: Area

Ref.channel: ch1

ISTD:


Formula: Linear

Weight: 1

| Level | Height | Area  | Conc. | Vol/Dil | Retention | Used | File |
|-------|--------|-------|-------|---------|-----------|------|------|
| 1     | 0      | 0     | 0     | 0       | 0         |      |      |
| 2     | 3.872  | 206.4 | 1     | 20      | 22.61     |      |      |
| 3     | 9.929  | 515   | 2     | 20      | 22.61     |      |      |
| 4     | 10.64  | 586.4 | 2.5   | 20      | 22.61     |      |      |
| 5     | 30.45  | 1506  | 5     | 20      | 22.61     |      |      |
| 6     | 69.97  | 3312  | 10    | 20      | 22.61     |      |      |
| 7     | 93.58  | 4344  | 12.5  | 20      | 22.61     |      |      |

CALIBRATION OF COMPONENT SO4

Method: AnionsIC2-022123.mtw  
 Equation:  $Q = 0.162666 \cdot A + 18.2166$   
 RSD: 4.213 %  
 Correlation coefficient: 0.999032



K3 = 0      K2 = 0      K1 = 0.162666      K0 = 18.2166  
 Base: Area  
 Ref.channel: ch1  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height | Area  | Conc. | Vol/Dil | Retention | Used | File |
|-------|--------|-------|-------|---------|-----------|------|------|
| 1     | 0      | 0     | 0     | 0       | 0         |      | 0    |
| 2     | 8.242  | 323.2 | 3     | 20      | 26.01     |      |      |
| 3     | 16.92  | 652.2 | 6     | 20      | 26.01     |      |      |
| 4     | 20.99  | 802.8 | 7.5   | 20      | 26.01     |      |      |
| 5     | 43.78  | 1643  | 15    | 20      | 26.01     |      |      |
| 6     | 94.3   | 3489  | 30    | 20      | 26.01     |      |      |
| 7     | 124.6  | 4590  | 37.5  | 20      | 26.01     |      |      |

Report date: 2/22/2023 11:53:29 AM  
Printed by: wet


Ident: STD1  
Analysis from: 2/21/2023 10:44:54 AM  
File: \_2023-02-21\_

Last save: 2/21/2023 4:00:42 PM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33358

Last save: 2/21/2023 2:38:

SAMPLE:  
:  
Vial number: 2'  
Volume: 20.0  $\mu$ L  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net  
METROHM LTD

Report date: 2/22/2023 11:53:38 AM  
Printed by: wet


Ident: STD2  
Analysis from: 2/21/2023 11:16:32 AM  
File: \_2023-02-21\_

Last save: 2/21/2023 4:05:53 PM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33359

Last save: 2/21/2023 2:38:

SAMPLE:  
: IZ/MA  
Vial number: 3  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.18             | 0.732          | 1.30         | 3.42        | 73.347         | 6.53      |
| 2  | 9.16             | 0.243          | 4.99         | 13.17       | 80.059         | 7.13      |
| 3  | 11.01            | 0.293          | 7.54         | 19.92       | 153.905        | 13.71     |
| 4  | 14.17            | 0.324          | 5.40         | 14.26       | 113.152        | 10.08     |
| 5  | 16.13            | 0.401          | 6.52         | 17.22       | 172.360        | 15.36     |
| 6  | 22.65            | 0.715          | 3.87         | 10.22       | 206.396        | 18.39     |
| 7  | 26.04            | 0.608          | 8.24         | 21.76       | 323.163        | 28.79     |
| 7  | 29.00            | 0.474          | 37.86        | 99.98       | 1122.383       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 2/22/2023 11:53:46 AM  
Printed by: wet


Ident: STD3  
Analysis from: 2/21/2023 11:48:10 AM  
File: \_2023-02-21\_

Last save: 2/21/2023 4:00:42 PM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33360

Last save: 2/21/2023 2:38:

SAMPLE:  
: IZ/MA  
Vial number: 4  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.11             | 0.668          | 3.15         | 3.82        | 161.313        | 6.62      |
| 2  | 9.09             | 0.242          | 10.35        | 12.57       | 170.040        | 6.98      |
| 3  | 10.92            | 0.292          | 16.86        | 20.47       | 346.928        | 14.24     |
| 4  | 14.05            | 0.320          | 11.55        | 14.02       | 239.776        | 9.84      |
| 5  | 15.97            | 0.393          | 13.60        | 16.51       | 351.707        | 14.43     |
| 6  | 22.61            | 0.698          | 9.93         | 12.06       | 515.049        | 21.13     |
| 7  | 26.01            | 0.598          | 16.92        | 20.55       | 652.177        | 26.76     |
| 7  | 29.00            | 0.459          | 82.35        | 99.99       | 2436.989       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 2/22/2023 11:53:54 AM  
Printed by: wet


Ident: STD4  
Analysis from: 2/21/2023 1:37:38 PM  
File: \_2023-02-21\_

Last save: 2/21/2023 4:00:42 PM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33364

Last save: 2/21/2023 2:38:

SAMPLE:  
: IZ/MA  
Vial number: 5  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.19             | 0.647          | 4.37         | 4.32        | 208.153        | 7.03      |
| 2  | 9.15             | 0.241          | 12.75        | 12.61       | 207.890        | 7.02      |
| 3  | 10.96            | 0.289          | 20.85        | 20.62       | 418.231        | 14.12     |
| 4  | 14.07            | 0.318          | 14.44        | 14.28       | 298.531        | 10.08     |
| 5  | 15.98            | 0.390          | 17.08        | 16.89       | 440.103        | 14.86     |
| 6  | 22.75            | 0.745          | 10.64        | 10.52       | 586.352        | 19.80     |
| 7  | 26.14            | 0.594          | 20.98        | 20.75       | 802.818        | 27.10     |
| 7  | 29.00            | 0.461          | 101.11       | 100.00      | 2962.078       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 2/22/2023 11:54:01 AM  
Printed by: wet


Ident: STD5  
Analysis from: 2/21/2023 2:09:16 PM  
File: \_2023-02-21\_

Last save: 2/21/2023 4:00:42 PM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33365

Last save: 2/21/2023 2:38:

SAMPLE:  
: IZ/MA  
Vial number: 6  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.13             | 0.677          | 8.63         | 3.90        | 433.643        | 6.71      |
| 2  | 9.04             | 0.238          | 26.50        | 11.96       | 427.294        | 6.61      |
| 3  | 10.84            | 0.287          | 45.79        | 20.67       | 919.454        | 14.22     |
| 4  | 13.94            | 0.313          | 30.61        | 13.82       | 627.183        | 9.70      |
| 5  | 15.82            | 0.383          | 35.76        | 16.14       | 908.244        | 14.05     |
| 6  | 22.53            | 0.666          | 30.45        | 13.75       | 1506.127       | 23.30     |
| 7  | 25.94            | 0.583          | 43.78        | 19.77       | 1643.268       | 25.42     |
| 7  | 29.00            | 0.449          | 221.52       | 100.00      | 6465.213       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 2/22/2023 11:54:09 AM  
Printed by: wet


Ident: STD6  
Analysis from: 2/21/2023 2:40:53 PM  
File: \_2023-02-21\_

Last save: 2/21/2023 4:00:42 PM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33366

Last save: 2/21/2023 2:38:

SAMPLE:  
: IZ/MA  
Vial number: 7  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.13             | 0.634          | 20.89        | 4.32        | 979.528        | 7.02      |
| 2  | 9.05             | 0.234          | 56.44        | 11.68       | 902.087        | 6.46      |
| 3  | 10.85            | 0.291          | 98.56        | 20.40       | 2013.979       | 14.43     |
| 4  | 13.98            | 0.309          | 65.85        | 13.63       | 1337.635       | 9.58      |
| 5  | 15.82            | 0.375          | 77.05        | 15.95       | 1923.216       | 13.78     |
| 6  | 22.53            | 0.639          | 69.97        | 14.48       | 3311.775       | 23.73     |
| 7  | 25.97            | 0.574          | 94.30        | 19.52       | 3488.521       | 25.00     |
| 7  | 29.00            | 0.436          | 483.06       | 100.00      | 13956.741      | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 2/22/2023 11:54:18 AM  
Printed by: wet


Ident: STD7  
Analysis from: 2/21/2023 3:12:31 PM  
File: \_2023-02-21\_

Last save: 2/21/2023 4:00:42 PM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33367

Last save: 2/21/2023 2:38:

SAMPLE:  
: IZ/MA  
Vial number: 8  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.20             | 0.671          | 27.03        | 4.26        | 1308.751       | 7.14      |
| 2  | 9.10             | 0.234          | 73.41        | 11.57       | 1175.426       | 6.41      |
| 3  | 10.90            | 0.292          | 127.97       | 20.18       | 2629.874       | 14.35     |
| 4  | 14.05            | 0.309          | 86.08        | 13.57       | 1753.981       | 9.57      |
| 5  | 15.89            | 0.373          | 101.56       | 16.01       | 2527.479       | 13.79     |
| 6  | 22.61            | 0.628          | 93.58        | 14.76       | 4343.957       | 23.70     |
| 7  | 26.08            | 0.571          | 124.56       | 19.64       | 4590.316       | 25.04     |
| 7  | 29.00            | 0.440          | 634.19       | 100.00      | 18329.785      | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 2/22/2023 11:54:53 AM  
Printed by: wet


Ident: ICV  
Analysis from: 2/21/2023 5:02:05 PM  
File: \_2023-02-21\_

Last save: 2/21/2023 5:30:05 PM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33372

Last save: 2/21/2023 4:07:

SAMPLE:  
:  
Vial number: 23  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.14             | 0.670          | 9.05         | 3.81        | 446.779        | 6.47      |
| 2  | 9.07             | 0.236          | 28.58        | 12.03       | 459.826        | 6.65      |
| 3  | 10.89            | 0.285          | 50.64        | 21.32       | 1014.118       | 14.68     |
| 4  | 14.01            | 0.314          | 32.79        | 13.80       | 675.678        | 9.78      |
| 5  | 15.90            | 0.383          | 38.30        | 16.13       | 971.732        | 14.06     |
| 6  | 22.57            | 0.664          | 31.83        | 13.40       | 1581.625       | 22.89     |
| 7  | 26.14            | 0.589          | 46.31        | 19.50       | 1760.194       | 25.47     |
| 7  | 28.00            | 0.449          | 237.51       | 100.00      | 6909.952       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 2/22/2023 11:55:14 AM  
Printed by: wet


Ident: ICB  
Analysis from: 2/21/2023 6:03:26 PM  
File: \_2023-02-21\_

Last save: 2/21/2023 6:31:26 PM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33374

Last save: 2/21/2023 4:07:

SAMPLE:  
: IZ/MA  
Vial number: 10  
Volume: 20.0  $\mu$ L  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net  
METROHM LTD

Report date: 2/23/2023 10:15:29 AM  
Printed by: wet


Ident: CCV  
Analysis from: 2/22/2023 10:28:29 AM  
File: \_2023-02-22\_

Last save: 2/22/2023 10:56:29 AM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33380

Last save: 2/21/2023 4:07:

SAMPLE:  
: IZ/MA  
Vial number: 29  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.21             | 0.607          | 10.36        | 4.43        | 485.092        | 6.78      |
| 2  | 9.34             | 0.249          | 28.54        | 12.21       | 485.529        | 6.78      |
| 3  | 11.27            | 0.306          | 48.25        | 20.63       | 1032.022       | 14.42     |
| 4  | 14.59            | 0.342          | 31.54        | 13.49       | 707.970        | 9.89      |
| 5  | 16.63            | 0.416          | 36.55        | 15.63       | 1009.254       | 14.10     |
| 6  | 22.89            | 0.679          | 31.69        | 13.55       | 1604.901       | 22.42     |
| 7  | 26.58            | 0.609          | 46.89        | 20.05       | 1832.815       | 25.61     |
| 7  | 28.00            | 0.459          | 233.81       | 100.00      | 7157.582       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 2/23/2023 10:15:38 AM  
Printed by: wet


Ident: CCB  
Analysis from: 2/22/2023 11:01:47 AM  
File: \_2023-02-22\_

Last save: 2/22/2023 11:29:47 AM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33381

Last save: 2/21/2023 4:07:

SAMPLE:  
: IZ/MA  
Vial number: 12  
Volume: 20.0  $\mu$ L  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net  
METROHM LTD

Report date: 2/23/2023 10:15:46 AM  
Printed by: wet


Ident: LB124209BLW  
Analysis from: 2/22/2023 11:32:24 AM  
File: \_2023-02-22\_

Last save: 2/23/2023 10:13:50 AM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33382

Last save: 2/21/2023 4:07:

SAMPLE:  
: IZ/MA  
Vial number: 13  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net  
METROHM LTD

Report date: 2/23/2023 10:15:53 AM  
Printed by: wet


Ident: LB124209BSW  
Analysis from: 2/22/2023 12:03:02 PM  
File: \_2023-02-22\_

Last save: 2/23/2023 10:13:50 AM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33383

Last save: 2/21/2023 4:07:

SAMPLE:  
: IZ/MA  
Vial number: 14  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.18             | 0.620          | 10.43        | 4.29        | 491.580        | 6.72      |
| 2  | 9.23             | 0.244          | 29.54        | 12.14       | 487.082        | 6.66      |
| 3  | 11.10            | 0.299          | 49.77        | 20.44       | 1042.345       | 14.25     |
| 4  | 14.33            | 0.330          | 32.98        | 13.55       | 715.648        | 9.78      |
| 5  | 16.31            | 0.402          | 38.56        | 15.84       | 1027.689       | 14.05     |
| 6  | 22.78            | 0.668          | 33.81        | 13.89       | 1682.987       | 23.01     |
| 7  | 26.45            | 0.601          | 48.34        | 19.86       | 1867.462       | 25.53     |
| 7  | 28.00            | 0.452          | 243.43       | 100.00      | 7314.793       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 2/23/2023 10:16:02 AM  
Printed by: wet


Ident: O1634-01  
Analysis from: 2/22/2023 12:34:04 PM  
File: \_2023-02-22\_

Last save: 2/22/2023 1:02:05 PM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33384

Last save: 2/21/2023 4:07:

SAMPLE:  
: IZ/MA  
Vial number: 15  
Volume: 20.0  $\mu$ L  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net  
METROHM LTD

Report date: 2/23/2023 10:16:09 AM  
Printed by: wet


Ident: O1634-01MS  
Analysis from: 2/22/2023 1:04:42 PM  
File: \_2023-02-22\_

Last save: 2/22/2023 1:32:43 PM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33385

Last save: 2/21/2023 4:07:

SAMPLE:  
: IZ/MA  
Vial number: 16  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.18             | 0.689          | 8.50         | 3.83        | 437.164        | 6.61      |
| 2  | 9.17             | 0.242          | 26.87        | 12.12       | 440.561        | 6.66      |
| 3  | 11.02            | 0.294          | 45.99        | 20.74       | 946.308        | 14.30     |
| 4  | 14.22            | 0.326          | 30.23        | 13.63       | 645.893        | 9.76      |
| 5  | 16.16            | 0.397          | 36.14        | 16.30       | 951.817        | 14.38     |
| 6  | 22.72            | 0.670          | 29.89        | 13.48       | 1497.955       | 22.64     |
| 7  | 26.37            | 0.598          | 44.12        | 19.90       | 1697.477       | 25.65     |
| 7  | 28.00            | 0.459          | 221.74       | 100.00      | 6617.175       | 100.00    |

This report has been created by IC Net  
METROHM LTD

This report has been created by IC Net  
METROHM LTD

Report date: 2/23/2023 10:16:23 AM  
Printed by: wet

Ident: O1646-01  
Analysis from: 2/22/2023 2:05:58 PM  
File: \_2023-02-22\_

Last save: 2/22/2023 2:33:58 PM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33387

Last save: 2/21/2023 4:07:

SAMPLE:  
: IZ/MA  
Vial number: 18  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 2  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 3  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 4  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 5  | 16.16            | 0.405          | 2.80         | 99.95       | 74.749         | 100.00    |
| 6  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 7  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 7  | 28.00            | 0.058          | 2.80         | 99.95       | 74.749         | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 2/23/2023 10:16:30 AM  
Printed by: wet


Ident: 01646-02  
Analysis from: 2/22/2023 2:36:36 PM  
File: \_2023-02-22\_

Last save: 2/22/2023 3:04:36 PM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33388

Last save: 2/21/2023 4:07:

SAMPLE:  
: IZ/MA  
Vial number: 19  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 2  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 3  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 4  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 5  | 16.23            | 0.406          | 2.09         | 99.93       | 55.619         | 100.00    |
| 6  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 7  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 7  | 28.00            | 0.058          | 2.09         | 99.93       | 55.619         | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 2/23/2023 10:16:37 AM  
Printed by: wet


Ident: 01646-03  
Analysis from: 2/22/2023 3:07:14 PM  
File: \_2023-02-22\_

Last save: 2/22/2023 3:35:14 PM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33389

Last save: 2/21/2023 4:07:

SAMPLE:  
: IZ/MA  
Vial number: 20  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 2  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 3  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 4  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 5  | 16.17            | 0.402          | 2.37         | 99.96       | 62.185         | 100.00    |
| 6  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 7  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 7  | 28.00            | 0.057          | 2.37         | 99.96       | 62.185         | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 2/23/2023 10:16:56 AM  
Printed by: wet


Ident: CCV  
Analysis from: 2/22/2023 4:10:35 PM  
File: \_2023-02-22\_

Last save: 2/22/2023 4:38:35 PM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33391

Last save: 2/21/2023 4:07:

SAMPLE:  
: IZ/MA  
Vial number: 14  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.17             | 0.678          | 9.58         | 3.83        | 482.617        | 6.55      |
| 2  | 9.13             | 0.239          | 30.30        | 12.11       | 490.663        | 6.66      |
| 3  | 10.96            | 0.290          | 52.05        | 20.80       | 1058.334       | 14.36     |
| 4  | 14.12            | 0.320          | 34.21        | 13.67       | 718.103        | 9.74      |
| 5  | 16.04            | 0.388          | 40.27        | 16.09       | 1038.680       | 14.09     |
| 6  | 22.70            | 0.660          | 34.67        | 13.85       | 1715.575       | 23.28     |
| 7  | 26.36            | 0.589          | 49.14        | 19.64       | 1865.383       | 25.31     |
| 7  | 28.00            | 0.452          | 250.21       | 100.00      | 7369.355       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 2/23/2023 10:17:03 AM  
Printed by: wet


Ident: CCB  
Analysis from: 2/22/2023 4:41:36 PM  
File: \_2023-02-22\_

Last save: 2/22/2023 5:09:36 PM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33392

Last save: 2/21/2023 4:07:

SAMPLE:  
: IZ/MA  
Vial number: 12  
Volume: 20.0  $\mu$ L  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net  
METROHM LTD

Report date: 2/23/2023 10:17:10 AM  
Printed by: wet


Ident: 01627-08  
Analysis from: 2/22/2023 5:12:17 PM  
File: \_2023-02-22\_

Last save: 2/22/2023 5:40:17 PM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33393

Last save: 2/21/2023 4:07:

SAMPLE:  
: IZ/MA  
Vial number: 21  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.20             | 0.727          | 1.32         | 3.17        | 73.562         | 5.90      |
| 2  | 9.18             | 0.242          | 5.09         | 12.21       | 83.044         | 6.66      |
| 3  | 11.04            | 0.290          | 8.12         | 19.47       | 163.942        | 13.14     |
| 4  | 14.19            | 0.327          | 5.59         | 13.41       | 118.441        | 9.49      |
| 5  | 16.17            | 0.400          | 6.73         | 16.13       | 177.530        | 14.23     |
| 6  | 22.84            | 0.698          | 3.66         | 8.77        | 192.213        | 15.41     |
| 7  | 26.52            | 0.609          | 11.18        | 26.82       | 438.751        | 35.17     |
| 7  | 28.00            | 0.471          | 41.69        | 99.99       | 1247.484       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 2/23/2023 10:17:17 AM  
 Printed by: wet


Ident: 01627-07  
 Analysis from: 2/22/2023 5:42:55 PM  
 File: \_2023-02-22\_

Last save: 2/22/2023 6:10:55 PM

Method: AnionsIC2-022123.mtw  
 Run operator: wet  
 Analysis number: 33394

Last save: 2/21/2023 4:07:

SAMPLE:  
 : IZ/MA  
 Vial number: 22  
 Volume: 20.0 µL  
 Dilution: 1.00  
 Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.15             | 0.690          | 0.59         | 3.00        | 28.532         | 4.91      |
| 2  | 9.13             | 0.238          | 2.42         | 12.27       | 37.509         | 6.46      |
| 3  | 10.99            | 0.290          | 3.86         | 19.58       | 77.739         | 13.38     |
| 4  | 14.14            | 0.324          | 2.60         | 13.18       | 54.062         | 9.31      |
| 5  | 16.12            | 0.400          | 3.22         | 16.34       | 84.477         | 14.54     |
| 6  | 22.72            | 0.699          | 1.57         | 7.97        | 82.490         | 14.20     |
| 7  | 26.39            | 0.614          | 5.44         | 27.62       | 216.184        | 37.21     |
| 7  | 28.00            | 0.465          | 19.70        | 99.96       | 580.994        | 100.00    |

This report has been created by IC Net  
 METROHM LTD

Report date: 2/23/2023 10:17:26 AM  
Printed by: wet


Ident: O1627-07RE 0.25ML  
Analysis from: 2/22/2023 6:13:33 PM  
File: \_2023-02-22\_

Last save: 2/23/2023 9:20:30 AM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33395

Last save: 2/21/2023 4:07:

SAMPLE:  
: IZ/MA  
Vial number: 23  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.15             | 0.690          | 0.30         | 2.83        | 13.712         | 4.24      |
| 2  | 9.16             | 0.236          | 1.16         | 11.12       | 17.604         | 5.44      |
| 3  | 11.03            | 0.294          | 1.75         | 16.67       | 37.400         | 11.57     |
| 4  | 14.20            | 0.319          | 1.19         | 11.39       | 23.979         | 7.42      |
| 5  | 16.21            | 0.401          | 1.50         | 14.32       | 39.033         | 12.07     |
| 6  | 22.79            | 0.718          | 0.61         | 5.85        | 32.949         | 10.19     |
| 7  | 26.45            | 0.621          | 3.95         | 37.73       | 158.710        | 49.08     |
| 7  | 28.00            | 0.468          | 10.46        | 99.92       | 323.387        | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 2/23/2023 10:17:33 AM  
Printed by: wet


Ident: CCV  
Analysis from: 2/22/2023 6:44:15 PM  
File: \_2023-02-22\_

Last save: 2/22/2023 7:12:15 PM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33396

Last save: 2/21/2023 4:07:

SAMPLE:  
: IZ/MA  
Vial number: 14  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.23             | 0.693          | 9.68         | 3.94        | 494.708        | 6.74      |
| 2  | 9.20             | 0.241          | 29.67        | 12.07       | 486.345        | 6.63      |
| 3  | 11.05            | 0.292          | 50.15        | 20.40       | 1030.199       | 14.04     |
| 4  | 14.24            | 0.322          | 33.34        | 13.56       | 704.952        | 9.61      |
| 5  | 16.17            | 0.390          | 39.34        | 16.00       | 1020.375       | 13.90     |
| 6  | 22.86            | 0.666          | 35.18        | 14.31       | 1758.330       | 23.96     |
| 7  | 26.55            | 0.591          | 48.47        | 19.72       | 1843.869       | 25.13     |
| 7  | 28.00            | 0.456          | 245.82       | 100.00      | 7338.777       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 2/23/2023 10:17:40 AM  
Printed by: wet


Ident: CCB  
Analysis from: 2/22/2023 7:14:54 PM  
File: \_2023-02-22\_

Last save: 2/22/2023 7:42:54 PM

Method: AnionsIC2-022123.mtw  
Run operator: wet  
Analysis number: 33397

Last save: 2/21/2023 4:07:

SAMPLE:  
: IZ/MA  
Vial number: 13  
Volume: 20.0  $\mu$ L  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net  
METROHM LTD