



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 3/3/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 03/02/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 03/03/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB124328

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O1742-01	RBR251669	1	1.14	8.57	9.71	9.66	99.4	
O1742-02	RBR251669-EPH-2	2	1.14	8.60	9.74	9.73	99.9	
O1743-01	34390	3	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O1743-02	34392	4	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O1743-03	34393	5	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O1744-01	72-11939	6	1.18	8.59	9.77	8.88	89.6	
O1744-02	72-11939-EPH-2	7	1.13	8.45	9.58	7.97	80.9	
O1747-01	HEF461Z-1-1	8	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1747-02	HEF461Z-1-2	9	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1747-03	HEF461Z-2-1	10	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1747-04	HEF461Z-2-2	11	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1747-05	BC226596-1-1	12	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1747-06	BC226596-1-2	13	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1747-07	EFB379J-1-1	14	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1747-08	EFB379J-1-2	15	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1747-09	EFB6600-1-1	16	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1747-10	EFB6600-1-2	17	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1748-01	347	18	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O1748-02	345	19	1.00	1.00	2.00	2.00	100.0	debris
O1749-01	CF-3	20	1.44	8.22	9.66	8.19	82.1	
O1749-02	CF-4	21	1.13	8.80	9.93	8.34	81.9	
O1749-03	CF-5	22	1.13	8.43	9.56	8.1	82.7	
O1749-04	CF-6	23	1.14	8.65	9.79	8.14	80.9	
O1749-05	CF-7	24	1.19	8.55	9.74	8.02	79.9	
O1749-06	CF-8	25	1.15	8.78	9.93	8.57	84.5	
O1749-07	CF-9	26	1.15	8.81	9.96	8.52	83.7	
O1749-08	CF-10	27	1.19	8.55	9.74	8.43	84.7	
O1749-09	CF-11	28	1.18	8.35	9.53	8.39	86.3	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 3/3/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 03/02/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 03/03/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB124328

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O1749-10	CF-12	29	1.19	8.47	9.66	8.27	83.6	
O1749-11	CF-13	30	1.18	8.61	9.79	8.58	85.9	
O1749-12	CF-14	31	1.17	8.76	9.93	8.69	85.8	
O1749-13	CF-15	32	1.19	8.79	9.98	8.66	85.0	
O1749-14	CF-16	33	1.16	8.47	9.63	8.43	85.8	
O1749-15	CF-17	34	1.14	8.83	9.97	8.8	86.7	
O1749-16	CF-18	35	1.14	8.44	9.58	8.39	85.9	
O1758-01	FMC-0081	36	1.00	1.00	2.00	2.00	100.0	DEBRIS
O1758-02	FMC-1601	37	1.00	1.00	2.00	2.00	100.0	DEBRIS
O1758-08	FMC-1221	55	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O1758-09	FMC-1222	56	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O1760-01	RBR200058	38	1.02	8.61	9.63	7.7	77.6	
O1760-02	RBR200058-EPH-2	39	1.04	8.75	9.79	8.7	87.5	
O1760-03	72-11998	40	1.03	8.90	9.93	9.16	91.3	
O1760-04	72-11998-EPH-2	41	1.03	8.88	9.91	9.00	89.8	
O1761-01	BC259215-1-1	42	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1761-02	BC259215-1-2	43	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1761-03	BC248150-1-1	44	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1761-04	BC248150-1-2	45	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1761-05	BC248150-2-1	46	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1761-06	BC248150-2-2	47	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1763-01	CM-A-CRUSHED-MASONRY-P ILE-1	48	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
O1764-01	CS-A-CRUSHED-STONE-PIL E-1	49	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
O1765-01	HEC935T-1-1	50	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1765-02	HEC935T-1-2	51	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1765-03	KMA652M-1-1	52	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1765-04	KMA652M-1-2	53	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 3/3/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 03/02/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 03/03/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB124328

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O1766-01	D2934	54	1.00	1.00	2.00	2.00	100.0	DEBRIS
O1767-01	CE-45-030123	57	1.13	8.46	9.59	8.48	86.9	
O1767-03	CE-45-030123	58	1.15	8.80	9.95	8.97	88.9	
O1767-05	CE-54-030123	59	1.19	8.65	9.84	9.29	93.6	
O1767-07	CE-54-030123	60	1.13	8.45	9.58	9.22	95.7	
O1767-09	CE-63-030123	61	1.11	8.78	9.89	9.27	92.9	
O1767-11	CE-63-030123	62	1.18	8.45	9.63	9.00	92.5	
O1767-13	CE-33-030123	63	1.16	8.54	9.7	9.27	95.0	
O1767-15	CE-33-030123	64	1.16	8.81	9.97	8.89	87.7	
O1770-01	VNJ-223	65	1.15	8.75	9.9	8.73	86.6	
O1776-01	TR-6-030223	66	1.18	8.46	9.64	9.13	94.0	
O1776-02	TR-6-030223-EPH-2	67	1.11	8.67	9.78	8.87	89.5	
O1777-01	TRE-1710	68	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O1777-02	TRE-1711	69	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$