



PERCENT SOLID

Supervisor: Iwona

Analyst: JIGNESH

Date: 3/16/2023

OVENTEMP IN Celsius(°C): 107

Time IN: 17:00

In Date: 03/15/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:15

Out Date: 03/16/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB124506

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O1924-01	RB21136	20	1.14	8.37	9.51	9.32	97.7	
O1924-02	RB21136-EPH-2	21	1.19	8.62	9.81	9.61	97.7	
O1924-03	CHRT26700	22	1.18	8.76	9.94	9.77	98.1	
O1924-04	CHRT26700-EPH-2	23	1.11	8.55	9.66	9.57	98.9	
O1924-05	CHRT27051	24	1.18	8.41	9.59	9.00	93.0	
O1924-06	CHRT27051-EPH-2	25	1.17	8.81	9.98	9.11	90.1	
O1934-01	HIH006X-1-1	26	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1934-02	HIH006X-1-2	27	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1934-03	HIH006X-2-1	28	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1934-04	HIH006X-2-2	29	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1934-05	BC216351-1-1	30	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1934-06	BC216351-1-2	31	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1934-07	BC137632-1-1	32	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1934-08	BC137632-1-2	33	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1934-09	BC226480-1-1	34	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1934-10	BC226480-1-2	35	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1934-11	GHA560N-1-1	36	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1934-12	GHA560N-1-2	37	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1934-13	GHA528P-1-1	38	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1934-14	GHA528P-1-2	39	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1935-01	2676	40	1.00	1.00	2.00	2.00	100.0	PLASTIC -gloves
O1936-01	OR-03-03152023	1	1.14	8.83	9.97	9.41	93.7	
O1936-02	OR-03-03152023-EPH-2	2	1.14	8.83	9.97	9.41	93.7	
O1938-01	KMH2352-1-1	4	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1938-02	KMH2352-1-2	5	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1938-03	JEH502N-1-1	6	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1938-04	JEH502N-1-2	7	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1939-01	D2945	8	1.00	1.00	2.00	2.00	100.0	DEBRIS



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 3/16/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 03/15/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 03/16/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB124506

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O1939-02	D2946	9	1.00	1.00	2.00	2.00	100.0	DEBRIS
O1939-03	D2943	10	1.12	8.46	9.58	8.33	85.2	
O1939-04	20071	11	1.18	8.35	9.53	8.81	91.4	
O1939-05	20071-EPH-2	12	1.13	8.69	9.82	7.84	77.2	
O1942-01	350	13	1.14	8.36	9.5	8.92	93.1	
O1943-01	BC204970-1-1	14	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1943-02	BC204970-1-2	15	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1943-03	BC107734-1-1	16	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1943-04	BC107734-1-2	17	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1943-05	HEE436L-1-1	18	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1943-06	HEE436L-1-2	19	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1944-01	MC-1	3	1.15	8.59	9.74	8.86	89.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$