



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 3/20/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 14:00
In Date: 03/17/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 03/18/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB124538

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O1955-01	DCAP9	1	1.00	1.00	2.00	2.00	100.0	P.E.SAMPLE
O1955-02	DCAQ0	2	1.00	1.00	2.00	2.00	100.0	P.E.SAMPLE
O1955-03	MDCAQ2	3	1.15	8.81	9.96	7.26	69.4	
O1955-04	MDCAQ3	4	1.12	8.51	9.63	7.67	77.0	
O1955-05	MDCAQ4	5	1.14	8.63	9.77	7.36	72.1	
O1955-06	MDCAQ6	6	1.13	8.64	9.77	7.35	72.0	
O1955-07	MDCAQ7	7	1.16	8.81	9.97	7.56	72.6	
O1955-08	MDCAQ8	8	1.15	8.35	9.5	7.24	72.9	
O1955-09	MDCAQ8D	9	1.15	8.35	9.5	7.24	72.9	
O1955-10	MDCAQ8S	10	1.15	8.35	9.5	7.24	72.9	
O1955-11	MDCAQ9	11	1.19	8.56	9.75	7.01	68.0	
O1955-12	MDCAR0	12	1.15	8.59	9.74	6.92	67.2	
O1955-13	MDCAR1	13	1.17	8.81	9.98	8.54	83.7	
O1955-14	MDCAR2	14	1.16	8.75	9.91	7.6	73.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$