

Clear table

| Instrument ID: IC-2 |                      |                       |                        |                        |                        |                         |                        |             |               | Analyst : IZ             |         | Method: 300.0 / 9056A |  |
|---------------------|----------------------|-----------------------|------------------------|------------------------|------------------------|-------------------------|------------------------|-------------|---------------|--------------------------|---------|-----------------------|--|
| ident               | concentra<br>tion F- | concentratio<br>n CL- | concentratio<br>on NO2 | concentratio<br>on BR- | concentratio<br>on NO3 | concentratio<br>on HPO4 | concentratio<br>on SO4 | file name   | date time     | Initial wt/<br>Final Vol | Analyst |                       |  |
| STD1                | 0.0000               | 0.0000                | 0.0000                 | 0.0000                 | 0.0000                 | 0.0000                  | 0.0000                 | 2023-03-22_ | 3/22/23 11:54 |                          | IZ/MA   |                       |  |
| STD2                | 0.4980               | 0.6810                | 0.6980                 | 2.3200                 | 0.5820                 | 0.9910                  | 3.5750                 | 2023-03-22_ | 3/22/23 12:26 |                          | IZ/MA   |                       |  |
| STD3                | 0.8100               | 1.2370                | 1.2280                 | 4.1060                 | 1.0300                 | 2.0120                  | 6.1800                 | 2023-03-22_ | 3/22/23 12:57 |                          | IZ/MA   |                       |  |
| STD4                | 0.9650               | 1.4990                | 1.5110                 | 5.0150                 | 1.2470                 | 2.4290                  | 7.4860                 | 2023-03-22_ | 3/22/23 13:29 |                          | IZ/MA   |                       |  |
| STD5                | 1.9180               | 2.8760                | 2.8560                 | 9.5370                 | 2.3830                 | 5.0910                  | 14.2000                | 2023-03-22_ | 3/22/23 14:01 |                          | IZ/MA   |                       |  |
| STD6                | 3.9240               | 5.8860                | 5.8680                 | 19.5610                | 4.8950                 | 10.0090                 | 29.3250                | 2023-03-22_ | 3/22/23 14:32 |                          | IZ/MA   |                       |  |
| STD7                | 5.0850               | 7.6220                | 7.6400                 | 25.4620                | 6.3630                 | 12.4680                 | 38.2350                | 2023-03-22_ | 3/22/23 15:04 |                          | IZ/MA   |                       |  |
| ICV                 | 2.0060               | 3.2520                | 3.2630                 | 10.8310                | 2.6890                 | 4.7380                  | 16.0140                | 2023-03-22_ | 3/22/23 16:43 |                          | IZ/MA   |                       |  |
| ICB                 | 0.0000               | 0.0000                | 0.0000                 | 0.0000                 | 0.0000                 | 0.0000                  | 0.0000                 | 2023-03-22_ | 3/22/23 17:46 |                          | IZ/MA   |                       |  |
| CCV                 | 2.0700               | 3.1330                | 3.0040                 | 10.3520                | 2.5850                 | 4.8810                  | 15.3100                | 2023-03-27_ | 3/27/23 19:14 |                          | IZ/MA   |                       |  |
| CCB                 | 0.0000               | 0.0000                | 0.0000                 | 0.0000                 | 0.0000                 | 0.0000                  | 0.0000                 | 2023-03-27_ | 3/27/23 19:45 |                          | IZ/MA   |                       |  |
| LB124712BLS         | 0.0000               | 0.0000                | 0.0000                 | 0.0000                 | 0.0000                 | 0.0000                  | 0.0000                 | 2023-03-27_ | 3/27/23 20:17 | 5.00g/100m               | IZ/MA   |                       |  |
| LB124712BSS         | 2.0970               | 3.1990                | 3.0240                 | 10.3730                | 2.5930                 | 4.8310                  | 15.3500                | 2023-03-27_ | 3/27/23 20:49 | 5.00g/100m               | IZ/MA   |                       |  |
| O2023-01            | 0.0000               | 11.1580               | 0.0000                 | 0.0000                 | 0.0000                 | 0.0000                  | 5.6100                 | 2023-03-27_ | 3/27/23 21:20 | 5.04g/100m               | IZ/MA   |                       |  |
| O2023-03            | 0.0000               | 1.1840                | 0.0000                 | 0.0000                 | 0.0000                 | 0.1690                  | 1.3280                 | 2023-03-27_ | 3/27/23 21:52 | 5.02g/100m               | IZ/MA   |                       |  |
| O2023-05            | 0.0000               | 5.1790                | 0.0000                 | 0.0000                 | 0.3180                 | 0.0000                  | 2.8570                 | 2023-03-27_ | 3/27/23 22:24 | 5.06g/100m               | IZ/MA   |                       |  |
| O2070-01            | 0.0000               | 5.2040                | 0.0000                 | 0.0000                 | 0.0000                 | 0.0000                  | 2.5580                 | 2023-03-27_ | 3/27/23 22:55 | 5.05g/100m               | IZ/MA   |                       |  |
| O2070-03            | 0.0000               | 2.3290                | 0.0000                 | 0.0000                 | 0.0000                 | 0.0000                  | 2.0940                 | 2023-03-27_ | 3/27/23 23:27 | 5.07g/100m               | IZ/MA   |                       |  |
| O2070-03MS          | 1.8470               | 5.5560                | 2.9340                 | 9.9410                 | 2.5040                 | 4.4260                  | 16.0000                | 2023-03-27_ | 3/27/23 23:59 | 5.08g/100m               | IZ/MA   |                       |  |
| O2070-03MSD         | 1.7720               | 5.5600                | 2.9320                 | 9.9700                 | 2.5120                 | 4.1980                  | 15.8890                | 2023-03-28_ | 3/28/23 0:30  | 5.05g/100m               | IZ/MA   |                       |  |
| CCV                 | 2.0690               | 3.1420                | 3.0680                 | 10.3910                | 2.5970                 | 5.0300                  | 15.3540                | 2023-03-28_ | 3/28/23 1:02  |                          | IZ/MA   |                       |  |
| CCB                 | 0.0000               | 0.0000                | 0.0000                 | 0.0000                 | 0.0000                 | 0.0000                  | 0.0000                 | 2023-03-28_ | 3/28/23 1:34  |                          | IZ/MA   |                       |  |


Clear table

Instrument ID: IC-2      Analyst: IZ      Method: 300.0 / 9056A

| ident | concentra<br>tion F- | concentratio<br>n CL- | concentratio<br>on NO2 | concentratio<br>on BR- | concentratio<br>on NO3 | concentratio<br>on HPO4 | concentratio<br>on SO4 | file name   | date time     | Initial wt/<br>Final Vol | Analyst |
|-------|----------------------|-----------------------|------------------------|------------------------|------------------------|-------------------------|------------------------|-------------|---------------|--------------------------|---------|
| STD1  | 0.0000               | 0.0000                | 0.0000                 | 0.0000                 | 0.0000                 | 0.0000                  | 0.0000                 | 2023-03-22_ | 3/22/23 11:54 |                          | IZ/MA   |
| STD2  | 0.4980               | 0.6810                | 0.6980                 | 2.3200                 | 0.5820                 | 0.9910                  | 3.5750                 | 2023-03-22_ | 3/22/23 12:26 |                          | IZ/MA   |
| STD3  | 0.8100               | 1.2370                | 1.2280                 | 4.1060                 | 1.0300                 | 2.0120                  | 6.1800                 | 2023-03-22_ | 3/22/23 12:57 |                          | IZ/MA   |
| STD4  | 0.9650               | 1.4990                | 1.5110                 | 5.0150                 | 1.2470                 | 2.4290                  | 7.4860                 | 2023-03-22_ | 3/22/23 13:29 |                          | IZ/MA   |
| STD5  | 1.9180               | 2.8760                | 2.8560                 | 9.5370                 | 2.3830                 | 5.0910                  | 14.2000                | 2023-03-22_ | 3/22/23 14:01 |                          | IZ/MA   |
| STD6  | 3.9240               | 5.8860                | 5.8680                 | 19.5610                | 4.8950                 | 10.0090                 | 29.3250                | 2023-03-22_ | 3/22/23 14:32 |                          | IZ/MA   |
| STD7  | 5.0850               | 7.6220                | 7.6400                 | 25.4620                | 6.3630                 | 12.4680                 | 38.2350                | 2023-03-22_ | 3/22/23 15:04 |                          | IZ/MA   |
| ident | True<br>Value F-     | True Value<br>CL-     | True Value<br>NO2      | True Value<br>BR-      | True Value<br>NO3      | True Value<br>HPO4      | True Value<br>SO4      |             |               |                          |         |
| STD1  | 0.0000               | 0.0000                | 0.0000                 | 0.0000                 | 0.0000                 | 0.0000                  | 0.0000                 |             |               |                          |         |
| STD2  | 0.4000               | 0.6000                | 0.6000                 | 2.0000                 | 0.5000                 | 1.0000                  | 3.0000                 |             |               |                          |         |
| STD3  | 0.8000               | 1.2000                | 1.2000                 | 4.0000                 | 1.0000                 | 2.0000                  | 6.0000                 |             |               |                          |         |
| STD4  | 1.0000               | 1.5000                | 1.5000                 | 5.0000                 | 1.2500                 | 2.5000                  | 7.5000                 |             |               |                          |         |
| STD5  | 2.0000               | 3.0000                | 3.0000                 | 10.0000                | 2.5000                 | 5.0000                  | 15.0000                |             |               |                          |         |
| STD6  | 4.0000               | 6.0000                | 6.0000                 | 20.0000                | 5.0000                 | 10.0000                 | 30.0000                |             |               |                          |         |
| STD7  | 5.0000               | 7.5000                | 7.5000                 | 25.0000                | 6.2500                 | 12.5000                 | 37.0000                |             |               |                          |         |
| ident | Relative<br>Error F- | Relative<br>Error CL- | Relative<br>Error NO2  | Relative<br>Error BR-  | Relative<br>Error NO3  | Relative<br>Error HPO4  | Relative<br>Error SO4  |             |               |                          |         |
| STD1  |                      |                       |                        |                        |                        |                         |                        |             |               |                          |         |
| STD2  | 24.5000              | 13.5000               | 16.3333                | 16.0000                | 16.4000                | -0.9000                 | 19.1667                |             |               |                          |         |
| STD3  | 1.2500               | 3.0833                | 2.3333                 | 2.6500                 | 3.0000                 | 0.6000                  | 3.0000                 |             |               |                          |         |
| STD4  | -3.5000              | -0.0667               | 0.7333                 | 0.3000                 | -0.2400                | -2.8400                 | -0.1867                |             |               |                          |         |
| STD5  | -4.1000              | -4.1333               | -4.8000                | -4.6300                | -4.6800                | 1.8200                  | -5.3333                |             |               |                          |         |
| STD6  | -1.9000              | -1.9000               | -2.2000                | -2.1950                | -2.1000                | 0.0900                  | -2.2500                |             |               |                          |         |
| STD7  | 1.7000               | 1.6267                | 1.8667                 | 1.8480                 | 1.8080                 | -0.2560                 | 3.3378                 |             |               |                          |         |

# CALIBRATION OF COMPONENT F-

Method: AnionsIC2-032223.mtw  
 Equation:  $Q = 0.079352 \cdot A + 3.7981$   
 RSD: 3.977 %  
 Correlation coefficient: 0.999137




K3 = 0      K2 = 0      K1 = 0.079352      K0 = 3.7981  
 Base: Area  
 Ref.channel: ch1  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height | Area  | Conc. | Vol/Dil | Retention | Used | File |
|-------|--------|-------|-------|---------|-----------|------|------|
| 1     | 0      | 0     | 0     | 0       | 0         |      | 0    |
| 2     | 1.24   | 77.73 | 0.4   | 20      | 6.113     |      |      |
| 3     | 2.655  | 156.2 | 0.8   | 20      | 6.113     |      |      |
| 4     | 3.558  | 195.3 | 1     | 20      | 6.113     |      |      |
| 5     | 7.666  | 435.6 | 2     | 20      | 6.113     |      |      |
| 6     | 17.01  | 941.2 | 4     | 20      | 6.113     |      |      |
| 7     | 23.02  | 1234  | 5     | 20      | 6.113     |      |      |

CALIBRATION OF COMPONENT CL-

Method: AnionsIC2-032223.mtw  
 Equation:  $Q = 0.128342 \cdot A + 3.42866$   
 RSD: 3.434 %  
 Correlation coefficient: 0.999357




K3 = 0      K2 = 0      K1 = 0.128342      K0 = 3.42866  
 Base: Area  
 Ref.channel: chl  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height | Area  | Conc. | Vol/Dil | Retention | Used | File |
|-------|--------|-------|-------|---------|-----------|------|------|
| 1     | 0      | 0     | 0     | 0       | 0         |      | 0    |
| 2     | 4.787  | 79.42 | 0.6   | 20      | 9.089     |      |      |
| 3     | 10.02  | 166   | 1.2   | 20      | 9.089     |      |      |
| 4     | 12.57  | 206.8 | 1.5   | 20      | 9.089     |      |      |
| 5     | 25.26  | 421.4 | 3     | 20      | 9.089     |      |      |
| 6     | 54.43  | 890.5 | 6     | 20      | 9.089     |      |      |
| 7     | 71.27  | 1161  | 7.5   | 20      | 9.089     |      |      |

# CALIBRATION OF COMPONENT NO2

Method: AnionsIC2-032223.mtw  
 Equation:  $Q = 0.0559228 \cdot A + 4.05126$   
 RSD: 3.955 %  
 Correlation coefficient: 0.999147




K3 = 0      K2 = 0      K1 = 0.0559228      K0 = 4.05126  
 Base: Area  
 Ref.channel: chl  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height | Area  | Conc. | Vol/Dil | Retention | Used File |
|-------|--------|-------|-------|---------|-----------|-----------|
| 1     | 0      | 0     | 0     | 0       | 0         | 0         |
| 2     | 8.7    | 177   | 0.6   | 20      | 10.92     |           |
| 3     | 18.14  | 366.8 | 1.2   | 20      | 10.92     |           |
| 4     | 22.79  | 467.9 | 1.5   | 20      | 10.92     |           |
| 5     | 45.96  | 948.8 | 3     | 20      | 10.92     |           |
| 6     | 97.18  | 2026  | 6     | 20      | 10.92     |           |
| 7     | 125    | 2660  | 7.5   | 20      | 10.92     |           |

# CALIBRATION OF COMPONENT BR-

Method: AnionsIC2-032223.mtw  
 Equation:  $Q = 0.283585 \cdot A + 15.1401$   
 RSD: 3.898 %  
 Correlation coefficient: 0.999172




K3 = 0      K2 = 0      K1 = 0.283585      K0 = 15.1401  
 Base: Area  
 Ref.channel: chl  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height | Area  | Conc. | Vol/Dil | Retention | Used File |
|-------|--------|-------|-------|---------|-----------|-----------|
| 1     | 0      | 0     | 0     | 0       | 0         | 0         |
| 2     | 5.175  | 110.2 | 2     | 20      | 14.05     |           |
| 3     | 11.07  | 236.2 | 4     | 20      | 14.05     |           |
| 4     | 13.94  | 300.3 | 5     | 20      | 14.05     |           |
| 5     | 28.56  | 619.2 | 10    | 20      | 14.05     |           |
| 6     | 62.13  | 1326  | 20    | 20      | 14.05     |           |
| 7     | 81.26  | 1742  | 25    | 20      | 14.05     |           |

# CALIBRATION OF COMPONENT NO3

Method: AnionsIC2-032223.mtw  
 Equation:  $Q = 0.0497148 \cdot A + 3.34859$   
 RSD: 3.867 %  
 Correlation coefficient: 0.999185




K3 = 0      K2 = 0      K1 = 0.0497148      K0 = 3.34859  
 Base: Area  
 Ref.channel: ch1  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height | Area  | Conc. | Vol/Dil | Retention | Used | File |
|-------|--------|-------|-------|---------|-----------|------|------|
| 1     | 0      | 0     | 0     | 0       | 0         |      | 0    |
| 2     | 6.277  | 166.9 | 0.5   | 20      | 15.97     |      |      |
| 3     | 13.12  | 347.1 | 1     | 20      | 15.97     |      |      |
| 4     | 16.42  | 434.2 | 1.25  | 20      | 15.97     |      |      |
| 5     | 33.37  | 891.3 | 2.5   | 20      | 15.97     |      |      |
| 6     | 72.25  | 1902  | 5     | 20      | 15.97     |      |      |
| 7     | 94.55  | 2492  | 6.25  | 20      | 15.97     |      |      |

# CALIBRATION OF COMPONENT HPO4

Method: AnionsIC2-032223.mtw  
 Equation:  $Q = 0.0591873 \cdot A + 0.736942$   
 RSD: 1.100 %  
 Correlation coefficient: 0.999934




K3 = 0      K2 = 0      K1 = 0.0591873      K0 = 0.736942  
 Base: Area  
 Ref.channel: ch1  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height | Area  | Conc. | Vol/Dil | Retention | Used | File |
|-------|--------|-------|-------|---------|-----------|------|------|
| 1     | 0      | 0     | 0     | 0       | 0         |      | 0    |
| 2     | 5.663  | 322.4 | 1     | 20      | 22.61     |      |      |
| 3     | 12.23  | 667.4 | 2     | 20      | 22.61     |      |      |
| 4     | 15.05  | 808.3 | 2.5   | 20      | 22.61     |      |      |
| 5     | 33.07  | 1708  | 5     | 20      | 22.61     |      |      |
| 6     | 69.62  | 3370  | 10    | 20      | 22.61     |      |      |
| 7     | 89.21  | 4201  | 12.5  | 20      | 22.61     |      |      |

# CALIBRATION OF COMPONENT SO4

Method: AnionsIC2-032223.mtw  
 Equation:  $Q = 0.162926 \cdot A + 20.1134$   
 RSD: 4.283 %  
 Correlation coefficient: 0.999000



K3 = 0      K2 = 0      K1 = 0.162926      K0 = 20.1134  
 Base: Area  
 Ref.channel: ch1  
 ISTD:  
 Formula: Linear  
 Weight: 1

| Level | Height | Area  | Conc. | Vol/Dil | Retention | Used | File |
|-------|--------|-------|-------|---------|-----------|------|------|
| 1     | 0      | 0     | 0     | 0       | 0         |      | 0    |
| 2     | 8.205  | 315.4 | 3     | 20      | 26.01     |      |      |
| 3     | 16.88  | 635.1 | 6     | 20      | 26.01     |      |      |
| 4     | 21.18  | 795.4 | 7.5   | 20      | 26.01     |      |      |
| 5     | 43.49  | 1620  | 15    | 20      | 26.01     |      |      |
| 6     | 95.29  | 3476  | 30    | 20      | 26.01     |      |      |
| 7     | 125.9  | 4570  | 37.5  | 20      | 26.01     |      |      |

Report date: 3/23/2023 10:30:18 AM  
Printed by: wet


Ident: STD1  
Analysis from: 3/22/2023 11:54:40 AM  
File: \_2023-03-22\_

Last save: 3/22/2023 3:42:45 PM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33651

Last save: 3/22/2023 3:01:

SAMPLE:  
: IZ/MA  
Vial number: 2  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net  
METROHM LTD

Report date: 3/23/2023 10:30:34 AM  
Printed by: wet


Ident: STD2  
Analysis from: 3/22/2023 12:26:18 PM  
File: \_2023-03-22\_

Last save: 3/22/2023 3:42:45 PM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33652

Last save: 3/22/2023 3:01:

SAMPLE:  
: IZ/MA  
Vial number: 3  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.30             | 0.846          | 1.24         | 3.09        | 77.731         | 6.22      |
| 2  | 9.30             | 0.245          | 4.79         | 11.95       | 79.419         | 6.36      |
| 3  | 11.19            | 0.294          | 8.70         | 21.73       | 177.026        | 14.17     |
| 4  | 14.39            | 0.326          | 5.17         | 12.92       | 110.244        | 8.83      |
| 5  | 16.43            | 0.397          | 6.28         | 15.67       | 166.892        | 13.36     |
| 6  | 22.84            | 0.752          | 5.66         | 14.14       | 322.370        | 25.81     |
| 7  | 26.30            | 0.582          | 8.20         | 20.49       | 315.375        | 25.25     |
| 7  | 29.00            | 0.492          | 40.04        | 99.99       | 1249.057       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 3/23/2023 10:30:43 AM  
Printed by: wet


Ident: STD3  
Analysis from: 3/22/2023 12:57:56 PM  
File: \_2023-03-22\_

Last save: 3/22/2023 3:42:45 PM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33653

Last save: 3/22/2023 3:01:

SAMPLE:  
: IZ/MA  
Vial number: 4  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.25             | 0.802          | 2.65         | 3.16        | 156.247        | 6.07      |
| 2  | 9.25             | 0.243          | 10.02        | 11.91       | 165.985        | 6.45      |
| 3  | 11.14            | 0.293          | 18.13        | 21.56       | 366.820        | 14.25     |
| 4  | 14.34            | 0.324          | 11.07        | 13.16       | 236.179        | 9.17      |
| 5  | 16.36            | 0.394          | 13.12        | 15.60       | 347.138        | 13.48     |
| 6  | 22.70            | 0.721          | 12.23        | 14.54       | 667.368        | 25.92     |
| 7  | 26.15            | 0.571          | 16.88        | 20.07       | 635.145        | 24.67     |
| 7  | 29.00            | 0.478          | 84.11        | 99.99       | 2574.882       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 3/23/2023 10:30:50 AM  
 Printed by: wet


Ident: STD4  
 Analysis from: 3/22/2023 1:29:34 PM  
 File: \_2023-03-22\_

Last save: 3/22/2023 3:42:45 PM

Method: AnionsIC2-032223.mtw  
 Run operator: wet  
 Analysis number: 33654

Last save: 3/22/2023 3:01:

SAMPLE:  
 : IZ/MA  
 Vial number: 5  
 Volume: 20.0 µL  
 Dilution: 1.00  
 Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.25             | 0.750          | 3.56         | 3.37        | 195.251        | 6.09      |
| 2  | 9.29             | 0.244          | 12.57        | 11.91       | 206.846        | 6.45      |
| 3  | 11.19            | 0.297          | 22.79        | 21.60       | 467.873        | 14.58     |
| 4  | 14.40            | 0.325          | 13.94        | 13.21       | 300.287        | 9.36      |
| 5  | 16.43            | 0.395          | 16.42        | 15.56       | 434.155        | 13.53     |
| 6  | 22.75            | 0.709          | 15.05        | 14.26       | 808.344        | 25.20     |
| 7  | 26.21            | 0.569          | 21.17        | 20.07       | 795.436        | 24.79     |
| 7  | 29.00            | 0.470          | 105.49       | 99.99       | 3208.192       | 100.00    |

This report has been created by IC Net  
 METROHM LTD

Report date: 3/23/2023 10:30:55 AM  
Printed by: wet


Ident: STD5  
Analysis from: 3/22/2023 2:01:12 PM  
File: \_2023-03-22\_

Last save: 3/22/2023 3:42:45 PM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33655

Last save: 3/22/2023 3:01:

SAMPLE:  
: IZ/MA  
Vial number: 6  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.33             | 0.783          | 7.66         | 3.53        | 435.585        | 6.56      |
| 2  | 9.36             | 0.245          | 25.26        | 11.62       | 421.395        | 6.34      |
| 3  | 11.27            | 0.299          | 45.96        | 21.14       | 948.790        | 14.28     |
| 4  | 14.53            | 0.325          | 28.56        | 13.14       | 619.181        | 9.32      |
| 5  | 16.56            | 0.397          | 33.37        | 15.35       | 891.319        | 13.42     |
| 6  | 22.87            | 0.681          | 33.07        | 15.21       | 1707.881       | 25.71     |
| 7  | 26.37            | 0.564          | 43.49        | 20.01       | 1619.669       | 24.38     |
| 7  | 29.00            | 0.471          | 217.37       | 100.00      | 6643.821       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 3/23/2023 10:31:02 AM  
Printed by: wet


Ident: STD6  
Analysis from: 3/22/2023 2:32:50 PM  
File: \_2023-03-22\_

Last save: 3/22/2023 3:42:45 PM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33656

Last save: 3/22/2023 3:01:

SAMPLE:  
: IZ/MA  
Vial number: 7  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.31             | 0.769          | 17.01        | 3.63        | 941.242        | 6.76      |
| 2  | 9.31             | 0.239          | 54.43        | 11.63       | 890.495        | 6.39      |
| 3  | 11.21            | 0.302          | 97.18        | 20.77       | 2026.270       | 14.54     |
| 4  | 14.47            | 0.319          | 62.13        | 13.28       | 1326.133       | 9.52      |
| 5  | 16.46            | 0.390          | 72.25        | 15.44       | 1901.780       | 13.65     |
| 6  | 22.70            | 0.641          | 69.62        | 14.88       | 3369.622       | 24.19     |
| 7  | 26.21            | 0.551          | 95.30        | 20.37       | 3476.388       | 24.95     |
| 7  | 29.00            | 0.459          | 467.90       | 100.00      | 13931.929      | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 3/23/2023 10:31:09 AM  
Printed by: wet


Ident: STD7  
Analysis from: 3/22/2023 3:04:28 PM  
File: \_2023-03-22\_

Last save: 3/22/2023 3:42:45 PM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33657

Last save: 3/22/2023 3:01:

SAMPLE:  
: IZ/MA  
Vial number: 8  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.30             | 0.737          | 23.02        | 3.77        | 1233.709       | 6.83      |
| 2  | 9.35             | 0.238          | 71.27        | 11.68       | 1161.075       | 6.43      |
| 3  | 11.26            | 0.310          | 125.02       | 20.49       | 2659.744       | 14.73     |
| 4  | 14.56            | 0.319          | 81.26        | 13.32       | 1742.331       | 9.65      |
| 5  | 16.55            | 0.390          | 94.55        | 15.49       | 2492.438       | 13.80     |
| 6  | 22.74            | 0.625          | 89.21        | 14.62       | 4200.744       | 23.26     |
| 7  | 26.27            | 0.548          | 125.87       | 20.63       | 4570.049       | 25.30     |
| 7  | 29.00            | 0.452          | 610.20       | 100.00      | 18060.090      | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 3/23/2023 10:31:32 AM  
Printed by: wet


Ident: ICV  
Analysis from: 3/22/2023 4:43:00 PM  
File: \_2023-03-22\_

Last save: 3/22/2023 5:12:01 PM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33660

Last save: 3/22/2023 3:33:

SAMPLE:  
: IZ/MA  
Vial number: 9  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.29             | 0.739          | 8.37         | 3.46        | 457.682        | 6.37      |
| 2  | 9.38             | 0.243          | 28.98        | 11.98       | 479.992        | 6.68      |
| 3  | 11.32            | 0.302          | 52.68        | 21.78       | 1094.392       | 15.23     |
| 4  | 14.63            | 0.328          | 32.36        | 13.38       | 710.503        | 9.88      |
| 5  | 16.68            | 0.402          | 37.53        | 15.52       | 1014.298       | 14.11     |
| 6  | 22.79            | 0.651          | 32.27        | 13.34       | 1588.596       | 22.10     |
| 7  | 26.31            | 0.560          | 49.70        | 20.55       | 1842.401       | 25.63     |
| 7  | 29.00            | 0.461          | 241.90       | 100.00      | 7187.864       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 3/23/2023 10:31:39 AM  
Printed by: wet


Ident: ICB  
Analysis from: 3/22/2023 5:46:27 PM  
File: \_2023-03-22\_

Last save: 3/22/2023 6:15:27 PM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33662

Last save: 3/22/2023 3:33:

SAMPLE:  
: IZ/MA  
Vial number: 10  
Volume: 20.0  $\mu$ L  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net  
METROHM LTD

Report date: 3/28/2023 11:53:20 AM  
Printed by: wet


Ident: CCV  
Analysis from: 3/27/2023 7:14:15 PM  
File: \_2023-03-27\_

Last save: 3/27/2023 7:43:15 PM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33724

Last save: 3/22/2023 3:33:

SAMPLE:  
: IZ/MA  
Vial number: 11  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.38             | 0.839          | 7.89         | 3.51        | 473.906        | 6.79      |
| 2  | 9.43             | 0.251          | 27.04        | 12.01       | 461.566        | 6.61      |
| 3  | 11.36            | 0.306          | 47.57        | 21.13       | 1001.777       | 14.35     |
| 4  | 14.68            | 0.335          | 30.41        | 13.51       | 676.680        | 9.70      |
| 5  | 16.73            | 0.409          | 35.45        | 15.75       | 972.477        | 13.93     |
| 6  | 23.10            | 0.697          | 31.41        | 13.95       | 1636.723       | 23.45     |
| 7  | 26.57            | 0.592          | 45.31        | 20.13       | 1755.905       | 25.16     |
| 7  | 29.00            | 0.490          | 225.08       | 100.00      | 6979.034       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 3/28/2023 11:53:31 AM  
Printed by: wet


Ident: CCB  
Analysis from: 3/27/2023 7:45:52 PM  
File: \_2023-03-27\_

Last save: 3/27/2023 8:14:52 PM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33725

Last save: 3/22/2023 3:33:

SAMPLE:  
: IZ/MA  
Vial number: 12  
Volume: 20.0  $\mu$ L  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net  
METROHM LTD

Report date: 3/29/2023 9:32:14 AM  
Printed by: wet


Ident: LB124712BLS  
Analysis from: 3/27/2023 8:17:29 PM  
File: \_2023-03-27\_

Last save: 3/29/2023 9:32:08 AM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33726

Last save: 3/22/2023 3:33:

SAMPLE: 5.00g/100mL  
: IZ/MA  
Vial number: 13  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net  
METROHM LTD

Report date: 3/28/2023 11:54:00 AM  
Printed by: wet


Ident: LB124712BSS  
Analysis from: 3/27/2023 8:49:07 PM  
File: \_2023-03-27\_

Last save: 3/28/2023 11:53:54 AM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33727

Last save: 3/22/2023 3:33:

SAMPLE: 5.00g/100mL  
IZ/MA  
Vial number: 14  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.47             | 0.842          | 7.96         | 3.50        | 480.634        | 6.87      |
| 2  | 9.52             | 0.255          | 27.19        | 11.96       | 471.794        | 6.74      |
| 3  | 11.47            | 0.304          | 48.36        | 21.27       | 1009.115       | 14.42     |
| 4  | 14.77            | 0.331          | 30.96        | 13.61       | 678.189        | 9.69      |
| 5  | 16.79            | 0.403          | 36.07        | 15.86       | 975.847        | 13.95     |
| 6  | 23.23            | 0.691          | 31.30        | 13.77       | 1620.093       | 23.16     |
| 7  | 26.70            | 0.590          | 45.53        | 20.03       | 1760.833       | 25.17     |
| 7  | 29.00            | 0.488          | 227.37       | 100.00      | 6996.505       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 3/28/2023 11:54:15 AM  
Printed by: wet


Ident: O2023-01  
Analysis from: 3/27/2023 9:20:50 PM  
File: \_2023-03-27\_

Last save: 3/28/2023 11:50:49 AM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33728

Last save: 3/22/2023 3:33:

SAMPLE: 5.04g/100mL  
IZ/MA  
Vial number: 29  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 2  | 9.32             | 0.237          | 105.49       | 88.16       | 1712.004       | 75.18     |
| 3  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 4  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 5  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 6  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 7  | 26.35            | 0.607          | 14.17        | 11.84       | 565.253        | 24.82     |
| 7  | 29.00            | 0.120          | 119.66       | 100.00      | 2277.257       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 3/28/2023 11:54:23 AM  
Printed by: wet


Ident: O2023-03  
Analysis from: 3/27/2023 9:52:28 PM  
File: \_2023-03-27\_

Last save: 3/28/2023 11:50:49 AM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33729

Last save: 3/22/2023 3:33:

SAMPLE: 5.02g/100mL  
: IZ/MA  
Vial number: 30  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 2  | 9.33             | 0.252          | 9.28         | 85.17       | 157.854        | 65.15     |
| 3  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 4  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 5  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 6  | 23.05            | 0.771          | 0.74         | 6.78        | 44.816         | 18.50     |
| 7  | 26.56            | 0.686          | 0.87         | 8.01        | 39.618         | 16.35     |
| 7  | 29.00            | 0.244          | 10.89        | 99.96       | 242.288        | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 3/28/2023 11:54:38 AM  
Printed by: wet

Ident: O2023-05  
Analysis from: 3/27/2023 10:24:06 PM  
File: \_2023-03-27\_

Last save: 3/28/2023 11:50:49 AM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33730

Last save: 3/22/2023 3:33:

SAMPLE: 5.06g/100mL  
IZ/MA  
Vial number: 31  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 2  | 9.36             | 0.246          | 46.34        | 85.75       | 780.411        | 73.06     |
| 3  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 4  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 5  | 16.64            | 0.412          | 2.20         | 4.08        | 60.620         | 5.67      |
| 6  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 7  | 26.61            | 0.631          | 5.50         | 10.17       | 227.217        | 21.27     |
| 7  | 29.00            | 0.184          | 54.04        | 100.00      | 1068.248       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 3/28/2023 11:54:46 AM  
Printed by: wet


Ident: O2070-01  
Analysis from: 3/27/2023 10:55:44 PM  
File: \_2023-03-27\_

Last save: 3/28/2023 11:50:49 AM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33731

Last save: 3/22/2023 3:33:

SAMPLE: 5.05g/100mL  
IZ/MA  
Vial number: 32  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 2  | 9.27             | 0.243          | 47.25        | 91.14       | 784.212        | 80.45     |
| 3  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 4  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 5  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 6  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 7  | 26.36            | 0.630          | 4.59         | 8.86        | 190.586        | 19.55     |
| 7  | 29.00            | 0.125          | 51.84        | 100.00      | 974.797        | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 3/28/2023 11:54:54 AM  
Printed by: wet


Ident: O2070-03  
Analysis from: 3/27/2023 11:27:22 PM  
File: \_2023-03-27\_

Last save: 3/28/2023 11:50:49 AM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33732

Last save: 3/22/2023 3:33:

SAMPLE: 5.07g/100mL  
IZ/MA  
Vial number: 33  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 2  | 9.27             | 0.247          | 20.08        | 86.35       | 336.157        | 71.56     |
| 3  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 4  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 5  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 6  | 0.00             | 0.000          | 0.00         | 0.00        | 0.000          | 0.00      |
| 7  | 26.39            | 0.637          | 3.17         | 13.62       | 133.580        | 28.44     |
| 7  | 29.00            | 0.126          | 23.24        | 99.98       | 469.738        | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 3/28/2023 11:55:03 AM  
Printed by: wet


Ident: O2070-03MS  
Analysis from: 3/27/2023 11:59:00 PM  
File: \_2023-03-27\_

Last save: 3/28/2023 11:50:49 AM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33733

Last save: 3/22/2023 3:33:

SAMPLE: 5.08g/100mL  
IZ/MA  
Vial number: 34  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.35             | 0.818          | 7.12         | 2.92        | 417.670        | 5.85      |
| 2  | 9.31             | 0.245          | 50.20        | 20.58       | 839.107        | 11.74     |
| 3  | 11.18            | 0.302          | 46.86        | 19.21       | 976.860        | 13.67     |
| 4  | 14.40            | 0.325          | 30.12        | 12.35       | 647.699        | 9.07      |
| 5  | 16.35            | 0.396          | 35.46        | 14.54       | 939.817        | 13.15     |
| 6  | 23.04            | 0.718          | 27.69        | 11.35       | 1482.976       | 20.76     |
| 7  | 26.50            | 0.605          | 46.47        | 19.05       | 1840.618       | 25.76     |
| 7  | 29.00            | 0.487          | 243.92       | 100.00      | 7144.747       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 3/28/2023 11:55:11 AM  
Printed by: wet


Ident: 02070-03MSD  
Analysis from: 3/28/2023 12:30:38 AM  
File: \_2023-03-28\_

Last save: 3/28/2023 11:50:49 AM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33734

Last save: 3/22/2023 3:33:

SAMPLE: 5.05g/100mL  
IZ/MA  
Vial number: 35  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.28             | 0.783          | 7.05         | 2.88        | 398.853        | 5.66      |
| 2  | 9.23             | 0.243          | 50.77        | 20.78       | 839.672        | 11.93     |
| 3  | 11.09            | 0.300          | 47.18        | 19.31       | 976.176        | 13.86     |
| 4  | 14.27            | 0.322          | 30.49        | 12.48       | 649.766        | 9.23      |
| 5  | 16.20            | 0.393          | 35.89        | 14.69       | 943.320        | 13.40     |
| 6  | 22.88            | 0.714          | 26.44        | 10.82       | 1406.150       | 19.97     |
| 7  | 26.32            | 0.601          | 46.46        | 19.02       | 1826.951       | 25.95     |
| 7  | 29.00            | 0.479          | 244.28       | 100.00      | 7040.889       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 3/28/2023 11:55:18 AM  
Printed by: wet


Ident: CCV  
Analysis from: 3/28/2023 1:02:24 AM  
File: \_2023-03-28\_

Last save: 3/28/2023 1:31:24 AM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33735

Last save: 3/22/2023 3:33:

SAMPLE:  
: IZ/MA  
Vial number: 11  
Volume: 20.0 µL  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

| No | Retention<br>min | Width/2<br>min | Height<br>mV | Height<br>% | Area<br>mV*sec | Area<br>% |
|----|------------------|----------------|--------------|-------------|----------------|-----------|
| 1  | 6.25             | 0.719          | 8.94         | 3.83        | 473.591        | 6.70      |
| 2  | 9.23             | 0.245          | 27.88        | 11.94       | 462.912        | 6.55      |
| 3  | 11.08            | 0.301          | 49.19        | 21.06       | 1024.673       | 14.50     |
| 4  | 14.26            | 0.318          | 32.23        | 13.80       | 679.473        | 9.62      |
| 5  | 16.18            | 0.388          | 37.62        | 16.11       | 977.206        | 13.83     |
| 6  | 22.88            | 0.698          | 32.39        | 13.87       | 1687.346       | 23.88     |
| 7  | 26.35            | 0.594          | 45.27        | 19.39       | 1761.366       | 24.93     |
| 7  | 29.00            | 0.466          | 233.53       | 100.00      | 7066.568       | 100.00    |

This report has been created by IC Net  
METROHM LTD

Report date: 3/28/2023 11:55:26 AM  
Printed by: wet


Ident: CCB  
Analysis from: 3/28/2023 1:34:02 AM  
File: \_2023-03-28\_

Last save: 3/28/2023 2:03:01 AM

Method: AnionsIC2-032223.mtw  
Run operator: wet  
Analysis number: 33736

Last save: 3/22/2023 3:33:

SAMPLE:  
: IZ/MA  
Vial number: 12  
Volume: 20.0  $\mu$ L  
Dilution: 1.00  
Amount: 1.0000



Quantitation method: Custom

No peaks

This report has been created by IC Net  
METROHM LTD