



PERCENT SOLID

Supervisor: Iwona

Analyst: JIGNESH

Date: 4/24/2023

OVENTEMP IN Celsius(°C): 107

Time IN: 14:00

In Date: 04/21/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:00

Out Date: 04/22/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB125081

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O2418-01	EW8Y7	1	1.18	8.61	9.79	8.05	79.8	
O2418-02	EW8Y8	2	1.13	8.60	9.73	8.5	85.7	
O2418-03	EW8Y9	3	1.13	8.40	9.53	6.93	69.0	
O2418-04	EW8Z0	4	1.13	8.67	9.8	8.31	82.8	
O2418-05	EW8Z1	5	1.18	8.39	9.57	8.18	83.4	
O2418-06	EW8Z2	6	1.18	8.76	9.94	7.71	74.5	
O2418-07	EW8Z3	7	1.18	8.53	9.71	8.51	85.9	
O2418-08	EW8Z4	8	1.14	8.81	9.95	7.01	66.6	
O2418-09	EW8Z5	9	1.11	8.69	9.8	7.16	69.6	
O2418-10	EW8Z5MS	10	1.11	8.69	9.8	7.16	69.6	
O2418-11	EW8Z5MSD	11	1.11	8.69	9.8	7.16	69.6	
O2418-12	EW8Z6	12	1.19	8.62	9.81	7.7	75.5	
O2418-13	EW8Z7	13	1.16	8.35	9.51	7.74	78.8	
O2418-14	EW8Z8	14	1.13	8.46	9.59	7.45	74.7	
O2418-15	EW8Z9	15	1.11	8.49	9.6	8.1	82.3	
O2418-16	EW900	16	1.15	8.81	9.96	7.94	77.1	
O2418-17	EW901	17	1.14	8.85	9.99	7.84	75.7	
O2418-18	VHBLK001	18	1.00	1.00	2.00	2.00	100.0	VHBLK
O2418-19	EW8Z5	19	1.11	8.69	9.8	7.16	69.6	
O2418-20	EW8Z5MS	20	1.11	8.69	9.8	7.16	69.6	
O2418-21	EW8Z5MSD	21	1.11	8.69	9.8	7.16	69.6	
O2418-22	EW8Z6	22	1.19	8.62	9.81	7.7	75.5	
O2418-23	EW8Z7	23	1.16	8.35	9.51	7.74	78.8	
O2418-24	EW8Z8	24	1.13	8.46	9.59	7.45	74.7	
O2418-25	EW8Z9	25	1.11	8.49	9.6	8.1	82.3	
O2418-26	EW900	26	1.15	8.81	9.96	7.94	77.1	
O2418-27	EW901	27	1.14	8.85	9.99	7.84	75.7	
O2418-28	EW902	28	1.18	8.52	9.7	8.28	83.3	



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 4/24/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 14:00
In Date: 04/21/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 04/22/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB125081

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O2418-29	EW903	29	1.1	8.48	9.58	8.18	83.5	
O2418-30	EW904	30	1.18	8.79	9.97	8.54	83.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$