



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 5/25/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 18:00
In Date: 05/24/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:37
Out Date: 05/25/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB125623

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O2908-01	SB10A	1	1.14	8.38	9.52	8.12	83.3	
O2908-02	SB10B	2	1.15	8.81	9.96	8.52	83.7	
O2908-03	WC10	3	1.18	8.42	9.6	8.14	82.7	
O2908-04	SB09	4	1.18	8.72	9.9	8.63	85.4	
O2908-05	WC09	5	1.18	8.79	9.97	9.2	91.2	
O2908-06	DUP01	6	1.15	8.40	9.55	8.26	84.6	
O2908-07	SB13	7	1.16	8.63	9.79	8.77	88.2	
O2908-08	WC13	8	1.11	8.66	9.77	8.75	88.2	
O2913-01	TP-70A	9	1.16	8.62	9.78	8.76	88.2	
O2913-03	TP-70A-EPH-2	10	1.12	8.40	9.52	8.17	83.9	
O2915-01	01-A-01-B-01-C	11	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O2920-01	280608	29	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
O2920-02	280608-EPH-2	30	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
O2920-03	RBR-400027	31	1.15	8.42	9.57	8.47	86.9	
O2920-04	RBR-400027-EPH-2	32	1.19	8.78	9.97	8.55	83.8	
O2921-01	PCB-051223-103-08	12	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O2921-02	PCB-051223-103-09	13	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O2921-03	PCB-051223-103-10	14	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O2921-04	PCB-051223-103-11	15	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O2921-05	PCB-051223-103-12	16	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O2921-06	PCB-051223-103-013	17	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O2921-07	PCB-051223-103-14	18	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O2925-01	CONCRETE-052423	19	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
O2931-01	HEE103S-1-1	37	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2931-02	HEE103S-1-2	38	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2931-03	HEE103S-2-1	39	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2931-04	HEE103S-2-2	40	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2931-05	BC274768-1-1	41	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 5/25/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 18:00
In Date: 05/24/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:37
Out Date: 05/25/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB125623

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O2931-06	BC274768-1-2	42	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2931-07	BC270896-1-1	43	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2931-08	BC270896-1-2	44	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2931-09	BC107490-1-1	45	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2931-10	BC107490-1-2	46	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2931-11	BC107060-1-1	47	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2931-12	BC107060-1-2	48	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2931-13	KQA098A-1-1	49	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2931-14	KQA098A-1-2	50	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2933-01	17A	51	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2933-02	17B	52	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2934-01	BC247601-1-1	53	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2934-02	BC247601-1-2	54	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2934-03	KMA788P-1-1	55	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2934-04	KMA788P-1-2	56	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2934-05	KMA864L-1-1	57	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2934-06	KMA864L-1-2	58	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2934-07	BC245029-1-1	59	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2934-08	BC245029-1-2	60	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2934-09	BC122669-1-1	61	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2934-10	BC122669-1-2	62	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2935-01	ETGI-346	33	1.14	8.82	9.96	8.75	86.3	
O2935-02	ETGI-346-EPH-2	34	1.15	8.40	9.55	9.17	95.5	
O2938-01	SU-04-052423	23	1.12	8.65	9.77	9.1	92.3	
O2938-02	SU-04-052423-EPH-2	24	1.14	8.84	9.98	9.11	90.2	
O2939-01	HEE595N-1-1	25	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2939-02	HEE595N-1-2	26	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2939-03	HEE595N-2-1	27	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 5/25/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 18:00
In Date: 05/24/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:37
Out Date: 05/25/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB125623

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O2939-04	HEE595N-2-2	28	1.00	1.00	2.00	2.00	100.0	PILC REEL
O2940-01	RBR-251671	35	1.17	8.73	9.9	9.38	94.0	
O2940-02	RBR-251671-EPH-2	36	1.16	8.48	9.64	9.29	95.9	
O2941-01	EME-TOPSOIL-001	20	1.14	8.65	9.79	8.24	82.1	
O2941-02	TILCON-CLEANFILL-001	21	1.13	8.57	9.7	9.13	93.3	
O2941-03	GHOST-CLEANFILL-001	22	1.18	8.64	9.82	9.42	95.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$