



PERCENT SOLID

Supervisor: Sohil
Analyst: jignesh
Date: 7/5/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 06/22/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:31
Out Date: 06/23/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB126127

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3375-01	SB-5 (0-2)	1	1.14	8.74	9.88	8.55	84.8	
O3375-02	SB-5 (8-10)	2	1.15	8.83	9.98	9.31	92.4	
O3375-03	SB-2 (0-2)	3	1.16	8.81	9.97	8.37	81.8	
O3375-04	SB-2 (8-10)	4	1.17	8.74	9.91	8.75	86.7	
O3375-05	O3375-04MS	5	1.16	8.75	9.91	8.75	86.7	
O3375-06	O3375-04MSD	6	1.16	8.75	9.91	8.75	86.7	
O3375-07	SB-3 (0-2)	7	1.18	8.76	9.94	8.91	88.2	
O3375-08	SB-3 (8-10)	8	1.17	8.71	9.88	8.97	89.6	
O3380-01	HEH6520-1-1	21	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3380-02	HEH6520-1-2	22	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3380-03	HEH6520-2-1	23	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3380-04	HEH6520-2-2	24	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3381-01	COMP1-0622	25	1.17	8.49	9.66	8.71	88.8	
O3381-02	COMP2-0622	26	1.14	8.63	9.77	9.00	91.1	
O3381-03	8352	27	1.00	1.00	2.00	2.00	100.0	debris
O3388-01	TR-05-062223	28	1.15	8.81	9.96	9.14	90.7	
O3388-02	TR-05-062223-EPH-2	29	1.13	8.59	9.72	9.02	91.9	
O3392-01	RY-1	9	1.13	8.69	9.82	9.17	92.5	
O3392-02	RY-1	10	1.1	8.42	9.52	9.03	94.2	
O3392-03	SRY-1	11	1.1	8.61	9.71	9.00	91.8	
O3392-04	SRY-1	12	1.13	8.79	9.92	9.27	92.6	
O3392-05	SRY-3	13	1.14	8.55	9.69	9.22	94.5	
O3392-06	SRY-3	14	1.14	8.51	9.65	9.12	93.8	
O3392-07	SB-64	15	1.12	8.73	9.85	9.18	92.3	
O3392-08	SB-64	16	1.15	8.84	9.99	9.26	91.7	
O3393-01	SF-9	17	1.12	8.66	9.78	9.00	91.0	
O3393-02	SF-9	18	1.12	8.48	9.6	8.75	90.0	
O3393-03	SF-10	19	1.14	8.68	9.82	9.38	94.9	



PERCENT SOLID

Supervisor: Sohil
Analyst: jignesh
Date: 7/5/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 06/22/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:31
Out Date: 06/23/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB126127

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3393-04	SF-10	20	1.11	8.69	9.8	9.34	94.7	
O3394-01	ORA-1651	30	1.15	8.49	9.64	8.58	87.5	
O3394-03	ORA-1651-EPH-2	31	1.13	8.39	9.52	8.66	89.7	
O3398-01	WC-1	32	1.14	8.85	9.99	9.05	89.4	
O3398-02	WC-1-EPH-2	33	1.14	8.84	9.98	9.00	88.9	
O3399-01	CONC-WC-1	34	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 7/3/2023

OVENTEMP IN Celsius(°C): 107

Time IN: 17:00

In Date: 06/22/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:31

Out Date: 06/23/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB126127

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3375-01	SB-5(0-2)	1	1.14	8.74	9.88	8.55	84.8	
O3375-02	SB-5(8-10)	2	1.15	8.83	9.98	9.31	92.4	
O3375-03	SB-2(0-2)	3	1.16	8.81	9.97	8.37	81.8	
O3375-04	SB-2(8-10)	4	1.17	8.74	9.91	8.75	86.7	
O3375-05	O3375-04MS	5	1.16	8.75	9.91	8.75	86.7	
O3375-06	O3375-04MSD	6	1.16	8.75	9.91	8.75	86.7	
O3375-07	SB-3(0-2)	7	1.18	8.76	9.94	8.91	88.2	
O3375-08	SB-3(8-10)	8	1.17	8.71	9.88	8.97	89.6	
O3380-01	HEH6520-1-1	21	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3380-02	HEH6520-1-2	22	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3380-03	HEH6520-2-1	23	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3380-04	HEH6520-2-2	24	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3381-01	COMP1-0622	25	1.17	8.49	9.66	8.71	88.8	
O3381-02	COMP2-0622	26	1.14	8.63	9.77	9.00	91.1	
O3381-03	8352	27	1.00	1.00	2.00	2.00	100.0	debris
O3388-01	TR-05-062223	28	1.15	8.81	9.96	9.14	90.7	
O3388-02	TR-05-062223-EPH-2	29	1.13	8.59	9.72	9.02	91.9	
O3392-01	RY-1	9	1.13	8.69	9.82	9.17	92.5	
O3392-02	RY-1	10	1.1	8.42	9.52	9.03	94.2	
O3392-03	SRY-1	11	1.1	8.61	9.71	9.00	91.8	
O3392-04	SRY-1	12	1.13	8.79	9.92	9.27	92.6	
O3392-05	SRY-3	13	1.14	8.55	9.69	9.22	94.5	
O3392-06	SRY-3	14	1.14	8.51	9.65	9.12	93.8	
O3392-07	SB-64	15	1.12	8.73	9.85	9.18	92.3	
O3392-08	SB-64	16	1.15	8.84	9.99	9.26	91.7	
O3393-01	SF-9	17	1.12	8.66	9.78	9.00	91.0	
O3393-02	SF-9	18	1.12	8.48	9.6	8.75	90.0	
O3393-03	SF-10	19	1.14	8.68	9.82	9.38	94.9	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 7/3/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 06/22/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:31
Out Date: 06/23/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB126127

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3393-04	SF-10	20	1.11	8.69	9.8	9.34	94.7	
O3394-01	ORA-1651	30	1.15	8.49	9.64	8.58	87.5	
O3394-03	ORA-1651-EPH-2	31	1.13	8.39	9.52	8.66	89.7	
O3398-01	WC-1	32	1.14	8.85	9.99	9.05	89.4	
O3398-02	WC-1-EPH-2	33	1.14	8.84	9.98	9.00	88.9	
O3399-01	CONC-WC-1	34	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$