



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 6/26/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 13:25
In Date: 06/24/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 06/25/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-OVEN

QC:LB126171

Table with 9 columns: Lab ID, Client SampleID, Dish #, Dish Wt (g) (A), Sample Wt (g), Dish + Sample Wt (g) (B), Dish+Dry Sample Wt (g) (C), % Solid, Comments. It contains 19 rows of data for samples YBTX8 through YBW21.

% Solid = (C-A) * 100 / (B-A)