



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 6/29/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 06/28/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 06/29/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB126228

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3453-01	P-062623RB01	1	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O3453-02	P-062623RB02	2	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O3453-03	P-062623RB03	3	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O3453-04	P-062623RB04	4	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O3454-01	AT029P-MW01-SB02-00G	5	1.12	8.59	9.71	8.93	90.9	
O3454-02	AT029P-MW02S-SB02-00G	6	1.19	8.47	9.66	8.85	90.4	
O3454-03	AT029P-MW02D-SB02-00G	7	1.15	8.80	9.95	9.07	90.0	
O3455-01	SS092P-MW01SDSB02-00G	8	1.13	8.85	9.98	9.45	94.0	
O3455-02	SS092P-MW01I-SB02-00G	9	1.14	8.47	9.61	8.99	92.7	
O3455-03	SS092P-MW01S-SB02-00G	10	1.17	8.65	9.82	9.12	91.9	
O3462-01	SU-02-62823	11	1.16	8.40	9.56	8.81	91.1	
O3462-02	SU-02-62823-EFH-2	12	1.16	8.40	9.56	8.81	91.1	
O3463-01	BC274271-1-1	13	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3463-02	BC274271-1-2	14	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3463-03	BC260544-1-1	15	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3463-04	BC260544-1-2	16	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3464-01	364	17	1.14	8.59	9.73	9.22	94.1	
O3464-02	365	18	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O3465-01	3102	19	1.18	8.56	9.74	8.62	86.9	
O3465-02	3107-3109	20	1.12	8.87	9.99	6.95	65.7	
O3466-01	BC283069-1-1	21	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3466-02	BC283069-1-2	22	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3466-03	BC247799-1-1	23	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3466-04	BC247799-1-2	24	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3466-05	BC279304-1-1	25	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3466-06	BC279304-1-2	26	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3469-01	2849-51-52-53	27	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O3469-02	2854	28	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 6/29/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 06/28/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 06/29/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB126228

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3470-01	BC216351-1-1	29	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3470-02	BC216351-1-2	30	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3470-03	BC216351-2-1	31	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3470-04	BC216351-2-2	32	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3471-01	201	33	1.16	8.51	9.67	9.46	97.5	
O3471-02	203	34	1.15	8.39	9.54	8.98	93.3	
O3471-03	VNJ-232	35	1.17	8.50	9.67	9.02	92.4	
O3471-04	VNJ-232-EPH-2	36	1.11	8.77	9.88	9.26	92.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$