



PERCENT SOLID

Supervisor: mohammad

Analyst: JIGNESH

Date: 7/3/2023

OVENTEMP IN Celsius(°C): 107

Time IN: 17:05

In Date: 06/30/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:30

Out Date: 07/01/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB126287

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3472-01	A46K8	1	1.00	1.00	2.00	2.00	100.0	F.B.
O3472-02	A46P1	2	1.14	8.82	9.96	8.34	81.6	
O3472-03	A46P2	3	1.15	8.40	9.55	7.58	76.5	
O3472-04	A46P6	4	1.00	1.00	2.00	2.00	100.0	F.B.
O3472-05	VHBLK001	5	1.00	1.00	2.00	2.00	100.0	VHBLK
O3472-06	A46Q2	6	1.17	8.80	9.97	8.5	83.3	
O3472-07	A46Q2MS	7	1.17	8.80	9.97	8.5	83.3	
O3472-08	A46Q2MSD	8	1.17	8.80	9.97	8.5	83.3	
O3472-09	A46Q3	9	1.17	8.81	9.98	6.5	60.5	
O3472-10	A46Q4	10	1.15	8.45	9.6	8.55	87.6	
O3472-11	A46Q5	11	1.16	8.49	9.65	7.93	79.7	
O3472-12	A46S4	12	1.16	8.80	9.96	4.41	36.9	
O3472-13	A46S5	13	1.11	8.47	9.58	7.55	76.0	
O3472-14	A46S7	14	1.12	8.67	9.79	8.00	79.4	
O3472-15	A46S9	15	1.13	8.84	9.97	5.8	52.8	
O3472-16	A46T0	16	1.19	8.71	9.9	8.68	86.0	
O3472-17	A4656	17	1.00	1.00	2.00	2.00	100.0	P.E. SAMPLE
O3472-18	A46Q9	18	1.15	8.82	9.97	3.47	26.3	
O3472-19	A46R0	19	1.16	8.83	9.99	6.3	58.2	
O3472-20	A46R1	20	1.12	8.56	9.68	8.92	91.1	
O3472-21	A46R2	21	1.15	8.80	9.95	7.22	69.0	
O3472-22	A46R3	22	1.16	8.61	9.77	4.45	38.2	
O3472-23	A46R4	23	1.11	8.63	9.74	6.11	57.9	
O3472-24	A46R5	24	1.14	8.82	9.96	8.28	81.0	
O3472-25	A46P1	25	1.14	8.82	9.96	8.34	81.6	
O3472-26	A46P2	26	1.15	8.40	9.55	7.58	76.5	
O3472-27	A4664	27	1.00	1.00	2.00	2.00	100.0	P.E. SAMPLE
O3472-28	A4683	28	1.00	1.00	2.00	2.00	100.0	P.E. SAMPLE



PERCENT SOLID

Supervisor: mohammad

Analyst: JIGNESH

Date: 7/3/2023

OVENTEMP IN Celsius(°C): 107

Time IN: 17:05

In Date: 06/30/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:30

Out Date: 07/01/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB126287

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3472-29	A46Q2	29	1.17	8.80	9.97	8.5	83.3	
O3472-30	A46Q2MS	30	1.17	8.80	9.97	8.5	83.3	
O3472-31	A46Q2MSD	31	1.17	8.80	9.97	8.5	83.3	
O3472-32	A46S7	32	1.12	8.67	9.79	8.00	79.4	
O3472-33	A46T0	33	1.19	8.71	9.9	8.68	86.0	
O3472-34	A46S4	34	1.16	8.80	9.96	4.41	36.9	
O3472-35	A46S5	35	1.11	8.47	9.58	7.55	76.0	
O3472-36	A46S9	36	1.13	8.84	9.97	5.8	52.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$