



PERCENT SOLID

Supervisor: Iwona

Analyst: JIGNESH

Date: 7/18/2023

OVENTEMP IN Celsius(°C): 107

Time IN: 17:15

In Date: 07/17/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:11

Out Date: 07/18/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB126475

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3638-01	TP1	1	1.15	8.81	9.96	9.14	90.7	
O3638-03	TP1-A	2	1.15	8.60	9.75	8.52	85.7	
O3638-04	TP1-B	3	1.13	8.60	9.73	8.88	90.1	
O3638-05	TP1-C	4	1.13	8.58	9.71	8.74	88.7	
O3638-06	TP1-D	5	1.12	8.51	9.63	8.58	87.7	
O3638-07	TP1-E	6	1.13	8.72	9.85	9.05	90.8	
O3638-08	TP1-F	7	1.16	8.47	9.63	8.47	86.3	
O3638-09	TP1-G	8	1.14	8.79	9.93	9.16	91.2	
O3638-10	TP1-H	9	1.12	8.44	9.56	8.55	88.0	
O3638-11	TP2	10	1.12	8.65	9.77	8.52	85.5	
O3638-13	TP2-A	11	1.15	8.81	9.96	8.99	89.0	
O3638-14	TP2-B	12	1.13	8.45	9.58	8.38	85.8	
O3638-15	TP2-C	13	1.12	8.66	9.78	8.55	85.8	
O3638-16	TP2-D	14	1.18	8.50	9.68	8.68	88.2	
O3638-17	TP2-E	15	1.18	8.54	9.72	8.52	85.9	
O3638-18	TP2-F	16	1.19	8.66	9.85	8.61	85.7	
O3638-19	TP2-G	17	1.15	8.80	9.95	8.69	85.7	
O3638-20	TP2-H	18	1.12	8.41	9.53	8.24	84.7	
O3638-21	TP6	19	1.14	8.82	9.96	9.31	92.6	
O3638-23	TP6-A	20	1.16	8.80	9.96	9.29	92.4	
O3638-24	TP6-B	21	1.11	8.55	9.66	8.72	89.0	
O3638-25	TP6-C	22	1.19	8.63	9.82	9.14	92.1	
O3638-26	TP6-D	23	1.12	8.71	9.83	8.84	88.6	
O3638-27	TP6-E	24	1.14	8.42	9.56	8.86	91.7	
O3638-28	TP6-F	25	1.11	8.86	9.97	9.24	91.8	
O3638-29	TP6-G	26	1.12	8.58	9.7	8.67	88.0	
O3638-30	TP6-H	27	1.14	8.84	9.98	9.22	91.4	
O3638-31	TP7	28	1.18	8.72	9.9	8.96	89.2	



PERCENT SOLID

Supervisor: Iwona

Analyst: JIGNESH

Date: 7/18/2023

OVENTEMP IN Celsius(°C): 107

Time IN: 17:15

In Date: 07/17/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:11

Out Date: 07/18/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB126475

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3638-33	TP7-A	29	1.12	8.61	9.73	8.68	87.8	
O3638-34	TP7-B	30	1.11	8.52	9.63	8.6	87.9	
O3638-35	TP7-C	31	1.18	8.81	9.99	9.37	93.0	
O3638-36	TP7-D	32	1.11	8.86	9.97	8.76	86.3	
O3638-37	TP7-E	33	1.19	8.57	9.76	8.73	88.0	
O3638-38	TP7-F	34	1.19	8.52	9.71	8.72	88.4	
O3638-39	TP7-G	35	1.13	8.59	9.72	8.71	88.2	
O3638-40	TP7-H	36	1.15	8.83	9.98	9.14	90.5	
O3639-01	TP3	37	1.12	8.70	9.82	8.57	85.6	
O3639-03	TP3-A	38	1.12	8.68	9.8	8.48	84.8	
O3639-04	TP3-B	39	1.13	8.85	9.98	9.06	89.6	
O3639-05	TP3-C	40	1.18	8.72	9.9	8.53	84.3	
O3639-06	TP3-D	41	1.13	8.55	9.68	8.66	88.1	
O3639-07	TP3-E	42	1.19	8.41	9.6	8.48	86.7	
O3639-08	TP3-F	43	1.19	8.63	9.82	8.7	87.0	
O3639-09	TP3-G	44	1.18	8.68	9.86	8.74	87.1	
O3639-10	TP3-H	45	1.11	8.50	9.61	8.65	88.7	
O3639-11	TP4	46	1.19	8.53	9.72	8.67	87.7	
O3639-13	TP4-A	47	1.15	8.51	9.66	8.36	84.7	
O3639-14	TP4-B	48	1.14	8.55	9.69	8.6	87.3	
O3639-15	TP4-C	49	1.12	8.68	9.8	8.6	86.2	
O3639-16	TP4-D	50	1.19	8.72	9.91	8.81	87.4	
O3639-17	TP4-E	51	1.15	8.54	9.69	8.64	87.7	
O3639-18	TP4-F	52	1.12	8.65	9.77	8.72	87.9	
O3639-19	TP4-G	53	1.14	8.83	9.97	8.86	87.4	
O3639-20	TP4-H	54	1.12	8.85	9.97	8.65	85.1	
O3639-21	TP5	55	1.12	8.46	9.58	8.77	90.4	
O3639-23	TP5-A	56	1.14	8.81	9.95	9.28	92.4	



PERCENT SOLID

Supervisor: Iwona

Analyst: JIGNESH

Date: 7/18/2023

OVENTEMP IN Celsius(°C): 107

Time IN: 17:15

In Date: 07/17/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:11

Out Date: 07/18/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB126475

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3639-24	TP5-B	57	1.17	8.60	9.77	8.97	90.7	
O3639-25	TP5-C	58	1.15	8.81	9.96	9.21	91.5	
O3639-26	TP5-D	59	1.15	8.38	9.53	8.79	91.2	
O3639-27	TP5-E	60	1.19	8.43	9.62	8.91	91.6	
O3639-28	TP5-F	61	1.19	8.55	9.74	9.00	91.3	
O3639-29	TP5-G	62	1.18	8.60	9.78	8.91	89.9	
O3639-30	TP5-H	63	1.12	8.47	9.59	8.7	89.5	
O3639-31	TP10	64	1.12	8.77	9.89	8.64	85.7	
O3639-33	TP10-A	65	1.15	8.60	9.75	8.63	87.0	
O3639-34	TP10-B	66	1.15	8.59	9.74	8.65	87.3	
O3639-35	TP10-C	67	1.17	8.60	9.77	8.63	86.7	
O3639-36	TP10-D	68	1.16	8.83	9.99	9.47	94.1	
O3639-37	TP10-E	69	1.19	8.57	9.76	8.63	86.8	
O3639-38	TP10-F	70	1.16	8.82	9.98	9.5	94.6	
O3639-39	TP10-G	71	1.16	8.77	9.93	8.78	86.9	
O3639-40	TP10-H	72	1.15	8.81	9.96	9.14	90.7	
O3640-01	TP9	73	1.12	8.47	9.59	8.72	89.7	
O3640-03	TP9-A	74	1.17	8.53	9.7	8.77	89.1	
O3640-04	TP9-B	75	1.15	8.82	9.97	9.34	92.9	
O3640-05	TP9-C	76	1.12	8.56	9.68	8.83	90.1	
O3640-06	TP9-D	77	1.11	8.69	9.8	8.88	89.4	
O3640-07	TP9-E	78	1.13	8.69	9.82	8.92	89.6	
O3640-08	TP9-F	79	1.16	8.83	9.99	9.41	93.4	
O3640-09	TP9-G	80	1.19	8.59	9.78	8.94	90.2	
O3640-10	TP9-H	81	1.19	8.60	9.79	8.93	90.0	
O3640-11	TP8	82	1.12	8.80	9.92	8.91	88.5	
O3640-13	TP8-A	83	1.11	8.88	9.99	9.55	95.0	
O3640-14	TP8-B	84	1.19	8.72	9.91	8.81	87.4	



PERCENT SOLID

Supervisor: Iwona

Analyst: JIGNESH

Date: 7/18/2023

OVENTEMP IN Celsius(°C): 107

Time IN: 17:15

In Date: 07/17/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:11

Out Date: 07/18/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB126475

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3640-15	TP8-C	85	1.14	8.84	9.98	9.56	95.2	
O3640-16	TP8-D	86	1.19	8.70	9.89	8.79	87.4	
O3640-17	TP8-E	87	1.15	8.84	9.99	9.03	89.1	
O3640-18	TP8-F	88	1.12	8.45	9.57	8.58	88.3	
O3640-19	TP8-G	89	1.19	8.63	9.82	8.89	89.2	
O3640-20	TP8-H	90	1.14	8.53	9.67	8.59	87.3	
O3640-21	TP12	91	1.11	8.87	9.98	8.81	86.8	
O3640-23	TP12-A	92	1.16	8.60	9.76	8.56	86.0	
O3640-24	TP12-B	93	1.19	8.63	9.82	8.84	88.6	
O3640-25	TP12-C	94	1.17	8.52	9.69	8.66	87.9	
O3640-26	TP12-D	95	1.14	8.66	9.8	8.76	88.0	
O3640-27	TP12-E	96	1.16	8.48	9.64	8.65	88.3	
O3640-28	TP12-F	97	1.16	8.77	9.93	8.82	87.3	
O3640-29	TP12-G	98	1.12	8.45	9.57	8.41	86.3	
O3640-30	TP12-H	99	1.18	8.57	9.75	8.72	88.0	
O3640-31	TP11	100	1.12	8.62	9.74	8.7	87.9	
O3640-33	TP11-A	101	1.12	8.87	9.99	9.54	94.9	
O3640-34	TP11-B	102	1.14	8.83	9.97	8.7	85.6	
O3640-35	TP11-C	103	1.12	8.61	9.73	8.48	85.5	
O3640-36	TP11-D	104	1.15	8.79	9.94	9.1	90.4	
O3640-37	TP11-E	105	1.15	8.40	9.55	8.7	89.9	
O3640-38	TP11-F	106	1.15	8.82	9.97	8.96	88.5	
O3640-39	TP11-G	107	1.15	8.84	9.99	9.31	92.3	
O3640-40	TP11-H	108	1.17	8.59	9.76	8.82	89.1	
O3641-01	TP13	109	1.14	8.83	9.97	8.4	82.2	
O3641-03	TP13-A	110	1.19	8.47	9.66	8.34	84.4	
O3641-04	TP13-B	111	1.1	8.68	9.78	8.33	83.3	
O3641-05	TP13-C	112	1.13	8.37	9.5	8.15	83.9	



PERCENT SOLID

Supervisor: Iwona

Analyst: JIGNESH

Date: 7/18/2023

OVENTEMP IN Celsius(°C): 107

Time IN: 17:15

In Date: 07/17/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:11

Out Date: 07/18/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB126475

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3641-06	TP13-D	113	1.17	8.58	9.75	8.28	82.9	
O3641-07	TP13-E	114	1.13	8.69	9.82	8.42	83.9	
O3641-08	TP13-F	115	1.19	8.47	9.66	8.02	80.6	
O3641-09	TP13-G	116	1.14	8.37	9.51	8.16	83.9	
O3641-10	TP13-H	117	1.12	8.63	9.75	8.4	84.4	
O3641-11	TP14	118	1.18	8.52	9.7	8.4	84.7	
O3641-13	TP14-A	119	1.12	8.72	9.84	8.67	86.6	
O3641-14	TP14-B	120	1.14	8.80	9.94	8.62	85.0	
O3641-15	TP14-C	121	1.16	8.41	9.57	8.31	85.0	
O3641-16	TP14-D	122	1.17	8.58	9.75	8.5	85.4	
O3641-17	TP14-E	123	1.17	8.45	9.62	8.33	84.7	
O3641-18	TP14-F	124	1.15	8.82	9.97	9.09	90.0	
O3641-19	TP14-G	125	1.18	8.54	9.72	8.54	86.2	
O3641-20	TP14-H	126	1.18	8.66	9.84	8.24	81.5	
O3641-21	TP16	127	1.18	8.54	9.72	8.6	86.9	
O3641-23	TP16-A	128	1.16	8.74	9.9	8.78	87.2	
O3641-24	TP16-B	129	1.19	8.55	9.74	7.9	78.5	
O3641-25	TP16-C	130	1.17	8.77	9.94	8.75	86.4	
O3641-26	TP16-D	131	1.19	8.53	9.72	8.53	86.0	
O3641-27	TP16-E	132	1.15	8.80	9.95	8.28	81.0	
O3641-28	TP16-F	133	1.19	8.50	9.69	7.71	76.7	
O3641-29	TP16-G	134	1.19	8.47	9.66	7.4	73.3	
O3641-30	TP16-H	135	1.15	8.37	9.52	7.47	75.5	
O3641-31	TP15	136	1.12	8.71	9.83	8.37	83.2	
O3641-33	TP15-A	137	1.12	8.46	9.58	7.86	79.7	
O3641-34	TP15-B	138	1.15	8.38	9.53	7.93	80.9	
O3641-35	TP15-C	139	1.17	8.50	9.67	8.06	81.1	
O3641-36	TP15-D	140	1.16	8.68	9.84	8.8	88.0	



PERCENT SOLID

Supervisor: Iwona

Analyst: JIGNESH

Date: 7/18/2023

OVENTEMP IN Celsius(°C): 107

Time IN: 17:15

In Date: 07/17/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:11

Out Date: 07/18/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB126475

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3641-37	TP15-E	141	1.16	8.78	9.94	8.57	84.4	
O3641-38	TP15-F	142	1.15	8.84	9.99	8.4	82.0	
O3641-39	TP15-G	143	1.18	8.57	9.75	7.99	79.5	
O3641-40	TP15-H	144	1.15	8.80	9.95	8.32	81.5	
O3643-01	01A-01B-01C	145	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O3644-01	01A-01B-01C	146	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O3645-01	SB-02- (3-5)	147	1.15	8.81	9.96	8.75	86.3	
O3645-02	SB-04- (1-5)	148	1.15	8.36	9.51	8.78	91.3	
O3645-03	SB-07- (1-3)	149	1.13	8.75	9.88	8.09	79.5	
O3645-04	SB-08- (10.5-2.0)	150	1.13	8.40	9.53	7.4	74.6	
O3645-05	SB-09- (2.0-4.0)	151	1.19	8.71	9.9	8.11	79.4	
O3645-06	SB-10- (0.5-2.0)	152	1.19	8.43	9.62	8.21	83.3	
O3645-07	DUP	153	1.18	8.81	9.99	8.3	80.8	
O3645-09	O3645-02MS	154	1.15	8.36	9.51	8.78	91.3	
O3645-10	O3645-02MSD	155	1.15	8.36	9.51	8.78	91.3	
O3646-01	GATE-1	164	1.14	8.85	9.99	8.52	83.4	
O3646-03	GATE-1-EPH-2	165	1.12	8.85	9.97	8.53	83.7	
O3646-04	YARD-1	166	1.18	8.45	9.63	8.8	90.2	
O3646-06	YARD-1-EPH-2	167	1.1	8.79	9.89	8.36	82.6	
O3647-01	SAMPLE-1	156	1.16	8.50	9.66	8.89	90.9	
O3648-01	OK-01-071723	168	1.13	8.45	9.58	7.93	80.5	
O3648-02	OK-01-071723-EPH-2	169	1.13	8.45	9.58	7.93	80.5	
O3649-01	X182-S5	157	1.15	8.50	9.65	8.14	82.2	
O3649-02	X182-B4	158	1.15	8.81	9.96	8.58	84.3	
O3649-03	X182-S3	159	1.12	8.72	9.84	8.9	89.2	
O3649-04	X182-B3	160	1.13	8.53	9.66	8.51	86.5	
O3649-05	X182-B3-DUP	161	1.19	8.71	9.9	8.85	87.9	
O3649-06	X182-S4	162	1.16	8.83	9.99	9.39	93.2	



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 7/18/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 07/17/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 07/18/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB126475

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3650-01	NWB-1933	163	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3654-01	72-12016	170	1.19	8.57	9.76	8.96	90.7	
O3654-03	72-12016-EPH-2	171	1.18	8.54	9.72	9.15	93.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$