

**PERCENT SOLID**

Supervisor: Iwona

Analyst: JIGNESH

Date: 7/20/2023

OVENTEMP IN Celsius(°C): 107

Time IN: 17:15

In Date: 07/19/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:15

Out Date: 07/20/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB126510

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3566-02	61R-1B	1	1.12	8.60	9.72	8.99	91.5	
O3566-03	61R-1C	2	1.15	8.40	9.55	8.91	92.4	
O3566-05	61R-2B	3	1.16	8.50	9.66	8.63	87.9	
O3566-06	61R-2C	4	1.12	8.71	9.83	8.56	85.4	
O3566-08	61R-3B	5	1.15	8.60	9.75	9.3	94.8	
O3566-09	61R-3C	6	1.12	8.40	9.52	8.24	84.8	
O3566-11	61R-4B	7	1.19	8.51	9.7	9.04	92.2	
O3566-12	61R-4C	8	1.16	8.65	9.81	9.06	91.3	
O3566-14	61R-5B	9	1.18	8.60	9.78	9.05	91.5	
O3566-15	61R-5C	10	1.16	8.61	9.77	8.73	87.9	
O3566-17	61R-6B	11	1.15	8.39	9.54	8.37	86.1	
O3566-18	61R-6C	12	1.11	8.75	9.86	8.07	79.5	
O3566-20	61R-7B	13	1.17	8.50	9.67	8.92	91.2	
O3566-21	61R-7C	14	1.19	8.43	9.62	8.26	83.9	
O3566-23	61R-8B	15	1.17	8.60	9.77	9.14	92.7	
O3566-24	61R-8C	16	1.1	8.48	9.58	8.98	92.9	
O3566-26	61R-9B	17	1.15	8.75	9.9	9.32	93.4	
O3566-27	61R-9C	18	1.18	8.71	9.89	8.98	89.6	
O3566-29	61R-10B	19	1.15	8.76	9.91	8.42	83.0	
O3567-01	61R-10C	20	1.16	8.54	9.7	8.27	83.3	
O3567-03	61R-11B	21	1.15	8.81	9.96	8.4	82.3	
O3567-04	61R-11C	22	1.13	8.67	9.8	8.17	81.2	
O3567-06	61R-12B	23	1.19	8.55	9.74	8.84	89.5	
O3567-07	61R-12C	24	1.11	8.80	9.91	8.62	85.3	
O3567-09	61R-13B	25	1.19	8.47	9.66	8.22	83.0	
O3567-10	61R-13C	26	1.17	8.43	9.6	8.7	89.3	
O3567-12	61R-14B	27	1.18	8.62	9.8	8.26	82.1	
O3567-13	61R-14C	28	1.11	8.78	9.89	8.45	83.6	



PERCENT SOLID

Supervisor: Iwona

Analyst: JIGNESH

Date: 7/20/2023

OVENTEMP IN Celsius(°C): 107

Time IN: 17:15

In Date: 07/19/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:15

Out Date: 07/20/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB126510

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3567-15	61R-15B	29	1.19	8.61	9.8	8.27	82.2	
O3567-16	61R-15C	30	1.16	8.50	9.66	8.41	85.3	
O3567-18	61R-16B	31	1.13	8.74	9.87	8.5	84.3	
O3567-19	61R-16C	32	1.19	8.62	9.81	8.37	83.3	
O3567-21	61R-17B	33	1.15	8.82	9.97	9.27	92.1	
O3567-22	61R-17C	34	1.17	8.64	9.81	8.68	86.9	
O3568-21	61R-DUP4	35	1.13	8.48	9.61	8.94	92.1	
O3568-22	61R-DUP5	36	1.19	8.53	9.72	9.22	94.1	
O3568-24	61R-DUP7	37	1.16	8.70	9.86	8.63	85.9	
O3607-05	BG-1	80	1.13	8.49	9.62	8.02	81.2	
O3662-01	BUR-1276	38	1.15	8.84	9.99	9.24	91.5	
O3663-01	1	39	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3663-02	2	40	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3663-03	3	41	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3663-04	4	42	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3663-05	5	43	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3663-06	6	44	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3663-07	7	45	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3663-08	8	46	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3663-09	9	47	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3663-10	10	48	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3663-11	11	49	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3663-12	12	50	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3663-13	13	51	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3663-14	14	52	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3663-15	15	53	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3663-16	16	54	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3663-17	17	55	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 7/20/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 07/19/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 07/20/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB126510

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3663-18	18	56	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3663-19	19	57	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3663-20	20	58	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O3664-01	AUD-1021-1020-COMP	59	1.14	8.44	9.58	8.98	92.9	
O3664-02	AUD-1023	60	1.19	8.55	9.74	7.55	74.4	
O3665-01	40265	61	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3665-02	40266	62	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3665-03	40267	63	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3666-01	P001-04A	64	1.13	8.12	9.25	8.93	96.1	
O3666-02	P001-04B	65	1.15	5.85	7.00	4.64	59.7	
O3666-03	P001-04C	66	1.12	8.28	9.4	9.28	98.6	
O3666-04	P001-04D	67	1.19	8.42	9.61	9.05	93.3	
O3666-05	P001-04E	68	1.12	8.75	9.87	8.67	86.3	
O3666-06	P001-04F	69	1.16	4.88	6.04	5.28	84.4	
O3666-07	P001-04G	70	1.12	5.75	6.87	5.6	77.9	
O3666-08	P001-04H	71	1.15	8.40	9.55	5.63	53.3	
O3666-09	P001-04I	72	1.12	8.39	9.51	8.12	83.4	
O3666-10	P001-04J	73	1.15	8.62	9.77	8.14	81.1	
O3666-11	P001-04K	74	1.19	8.71	9.9	9.00	89.7	
O3666-12	P001-04L	75	1.18	7.04	8.22	4.86	52.3	
O3666-13	P001-WD04-COMP1-01	76	1.18	8.58	9.76	6.52	62.2	
O3666-14	P001-WD04-COMP2-01	77	1.12	8.81	9.93	6.84	64.9	
O3666-15	P001-WD04-COMP3-01	78	1.16	8.64	9.8	5.92	55.1	
O3666-16	P001-WD04-COMP4-01	79	1.15	7.35	8.5	6.27	69.7	
O3668-01	SLUDGE-DRUM	94	1.14	8.71	9.85	6.24	58.6	
O3668-03	332	89	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
O3668-04	332-EPH-2	90	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
O3669-01	DEL-AVE-1	83	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 7/20/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 07/19/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 07/20/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB126510

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3669-02	DEL-AVE-2	84	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
O3680-01	TP-3	81	1.14	8.85	9.99	8.84	87.0	
O3680-03	TP-3-EPH-2	82	1.12	8.84	9.96	8.94	88.5	
O3682-01	ELZ-23-00066-00070	85	1.15	8.38	9.53	7.55	76.4	
O3682-02	ELZ-23-00061-00062-00069	86	1.00	1.00	2.00	2.00	100.0	debris
O3682-03	ELZ-23-00056-00057-00058	87	1.11	8.58	9.69	9.48	97.6	
O3684-01	SF6-WASTE	91	1.00	1.00	2.00	2.00	100.0	plastic sample
O3685-01	CHRT-24152	92	1.17	8.45	9.62	9.49	98.5	
O3685-02	CHRT-24152-EPH-2	93	1.19	8.55	9.74	9.51	97.3	
O3687-01	CO-38882	88	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$