



PERCENT SOLID

Supervisor: Iwona

Analyst: JIGNESH

Date: 7/21/2023

OVENTEMP IN Celsius(°C): 107

Time IN: 17:00

In Date: 07/20/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:22

Out Date: 07/21/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB126526

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3674-01	WO-7	1	1.14	8.76	9.9	8.51	84.1	
O3674-02	WO-8	2	1.16	8.64	9.8	8.57	85.8	
O3674-03	WO-9	3	1.15	8.81	9.96	9.06	89.8	
O3677-01	SWC-1	4	1.16	8.83	9.99	7.4	70.7	
O3678-01	LWC-1	5	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O3678-02	LWC-2	6	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O3678-03	LWC-3	7	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O3681-01	RPTP-1-CS-1	8	1.12	8.55	9.67	8.73	89.0	
O3681-02	RPTP-1-GRAB-3	9	1.12	5.06	6.18	5.61	88.7	
O3681-03	RPTP-1-GRAB-4	10	1.15	5.11	6.26	5.49	84.9	
O3681-04	RPTP-2-CS-2	11	1.17	8.61	9.78	8.78	88.4	
O3681-05	RPTP-2-GRAB-1	12	1.13	5.63	6.76	6.06	87.6	
O3681-06	RPTP-2-GRAB-2	13	1.14	4.36	5.5	4.92	86.7	
O3681-07	RPTP-3-CS-3	14	1.19	8.74	9.93	9.00	89.4	
O3681-08	RPTP-3-GRAB-5	15	1.15	4.90	6.05	5.7	92.9	
O3681-09	RPTP-3-GRAB-6	16	1.16	5.46	6.62	5.87	86.3	
O3681-10	RPTP-4-CS-4	17	1.18	8.74	9.92	9.4	94.1	
O3681-11	RPTP-4-GRAB-7	18	1.12	4.86	5.98	5.7	94.2	
O3681-12	RPTP-4-GRAB-8	19	1.18	4.92	6.1	5.82	94.3	
O3681-13	RPTP-5-CS-5	20	1.19	8.40	9.59	9.08	93.9	
O3681-14	RPTP-5-GRAB-9	21	1.18	4.23	5.41	5.09	92.4	
O3681-15	RPTP-5-GRAB-10	22	1.19	5.00	6.19	5.78	91.8	
O3681-16	RPTP-6-CS-6	23	1.18	8.44	9.62	9.12	94.1	
O3681-17	RPTP-6-GRAB-11	24	1.19	4.55	5.74	5.44	93.4	
O3681-18	RPTP-6-GRAB-12	25	1.15	4.38	5.53	5.27	94.1	
O3681-19	RPTP-7-CS-7	26	1.19	8.47	9.66	9.1	93.4	
O3681-20	RPTP-7-GRAB-13	27	1.13	4.47	5.6	5.29	93.1	
O3681-21	RPTP-7-GRAB-14	28	1.18	4.49	5.67	5.29	91.5	



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 7/21/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 07/20/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 07/21/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB126526

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3681-22	RPTP-8-CS-8	29	1.18	8.59	9.77	9.25	93.9	
O3681-23	RPTP-8-GRAB-15	30	1.11	4.82	5.93	5.54	91.9	
O3681-24	RPTP-8-GRAB-16	31	1.14	4.48	5.62	5.62	100.0	
O3686-01	SB-21-0-2-20230719	32	1.15	8.81	9.96	9.38	93.4	
O3686-02	SB-21-16-18-20230719	33	1.16	8.83	9.99	9.7	96.7	
O3686-03	SB-26-0-2-20230719	34	1.12	8.68	9.8	9.21	93.2	
O3686-04	SB-26-15-17-20230719	35	1.12	8.46	9.58	9.22	95.7	
O3694-01	WC-1	36	1.16	8.50	9.66	7.94	79.8	
O3694-02	WC-1-EFH-A	37	1.18	8.63	9.81	8.03	79.4	
O3706-01	B-P9A	38	1.14	8.83	9.97	8.74	86.1	
O3706-02	B-P9B	39	1.12	8.84	9.96	9.1	90.3	
O3706-03	B-P10A	40	1.18	8.75	9.93	8.67	85.6	
O3706-04	B-P10B	41	1.11	8.46	9.57	8.28	84.8	
O3706-05	B-P20A	42	1.13	8.84	9.97	9.08	89.9	
O3706-06	B-P20B	43	1.13	8.81	9.94	8.96	88.9	
O3706-07	B-P21A	44	1.17	8.80	9.97	8.19	79.8	
O3706-08	B-P21B	45	1.16	8.79	9.95	8.66	85.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$