



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 7/27/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 07/26/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:27
Out Date: 07/27/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB126607

Table with 9 columns: Lab ID, Client SampleID, Dish #, Dish Wt(g) (A), Sample Wt(g), Dish + Sample Wt(g) (B), Dish+Dry Sample Wt(g) (C), % Solid, Comments. It contains 23 rows of data, with the last row (03785-23) marked as 'P.E. SAMPLE'.

% Solid = (C-A) * 100 / (B-A)