



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 8/11/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 08/10/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 08/11/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB126811

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3956-01	B-PP22	1	1.19	8.71	9.9	8.65	85.6	
O3956-02	B-PP23A	2	1.15	8.38	9.53	7.88	80.3	
O3956-03	B-PP23B	3	1.17	8.73	9.9	8.49	83.8	
O3958-01	NEVINS-SB02-Z106-107	4	1.17	8.51	9.68	9.04	92.5	
O3958-02	NEVINS-SB02-Z93-94	5	1.11	8.57	9.68	8.72	88.8	
O3958-03	NEVINS-SB02-Z64-65	6	1.18	8.45	9.63	8.57	87.5	
O3958-04	NEVINS-SB02-Z46-47	7	1.18	8.81	9.99	9.04	89.2	
O3958-05	NEVINS-SB02-Z37-39	8	1.19	8.78	9.97	8.82	86.9	
O3958-06	NEVINS-SB02-Z21-23	9	1.12	8.81	9.93	8.6	84.9	
O3958-07	NEVINS-SB02-Z10-12	10	1.15	8.40	9.55	8.44	86.8	
O3961-01	ETGI-339	21	1.16	8.50	9.66	8.79	89.8	
O3961-02	ETGI-339-EPH-2	22	1.18	8.53	9.71	8.26	83.0	
O3961-03	20071	23	1.17	8.38	9.55	8.55	88.1	
O3961-04	20071-EPH-2	24	1.15	8.39	9.54	8.3	85.2	
O3966-01	KMA772-A-END-1-1	35	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3966-02	KMA772-A-END-1-2	36	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3966-03	HEE695R-END-1-1	37	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3966-04	HEE695R-END-1-2	38	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3966-05	KMA576Q-END-1-1	39	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3966-06	KMA576Q-END-1-2	40	1.00	1.00	2.00	2.00	100.0	PILC REEL
O3969-01	COMP-1	32	1.00	1.00	2.00	2.00	100.0	debris
O3969-02	COMP-2	33	1.00	1.00	2.00	2.00	100.0	debris
O3969-03	SOM-23-00072	34	1.00	1.00	2.00	2.00	100.0	debris
O3971-01	B-PP02A	11	1.15	8.42	9.57	8.92	92.3	
O3971-02	B-PP01A	12	1.16	8.47	9.63	8.81	90.3	
O3971-03	WC-PP01	13	1.19	8.59	9.78	8.36	83.5	
O3971-04	B-PP02B	14	1.19	8.78	9.97	9.22	91.5	
O3971-05	B-PP041A	15	1.19	8.64	9.83	9.1	91.6	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 8/11/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 08/10/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 08/11/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB126811

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O3971-06	B-PP04B	16	1.19	8.61	9.8	8.6	86.1	
O3971-07	DUP-05	17	1.19	8.80	9.99	9.44	93.7	
O3971-08	B-PP03A	18	1.18	8.57	9.75	8.76	88.4	
O3971-09	B-PP08A	19	1.16	8.51	9.67	8.76	89.3	
O3971-10	B-PP08B	20	1.15	8.84	9.99	8.84	87.0	
O3972-01	367	25	1.17	8.60	9.77	7.45	73.0	
O3972-02	367-EPH-2	26	1.15	8.82	9.97	7.55	72.6	
O3972-03	356	27	1.00	1.00	2.00	2.00	100.0	debris
O3972-04	356-EPH-2	28	1.00	1.00	2.00	2.00	100.0	debris
O3972-05	2876-2879	29	1.14	8.81	9.95	9.26	92.2	
O3974-01	EO-01-081023	30	1.14	8.66	9.8	9.11	92.0	
O3974-02	EO-01-081023-EPH-2	31	1.14	8.66	9.8	9.11	92.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$