



PERCENT SOLID

Supervisor: Sohil
Analyst: jignesh
Date: 8/18/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 08/17/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:20
Out Date: 08/18/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB126897

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O4046-01	SB-02 (0-5)	1	1.18	8.39	9.57	8.92	92.3	
O4046-02	SB-06 (0-5)	2	1.15	8.84	9.99	9.24	91.5	
O4046-03	SB-01 (0-5)	3	1.15	8.47	9.62	9.31	96.3	
O4046-04	SB-07 (0-5)	4	1.15	8.78	9.93	9.44	94.4	
O4047-01	NEVINS-SB05-Z26-27	5	1.16	8.81	9.97	8.32	81.3	
O4047-02	NEVINS-SB05-Z60-61	6	1.19	8.53	9.72	8.17	81.8	
O4047-03	NEVINS-SB05-Z60-61-Q	7	1.12	8.78	9.9	8.36	82.5	
O4047-04	NEVINS-SB05-Z83-84	8	1.18	8.43	9.61	8.42	85.9	
O4054-01	WASTE	9	1.19	8.77	9.96	8.08	78.6	
O4054-02	VOC	10	1.15	8.80	9.95	8.26	80.8	
O4054-03	1	11	1.13	8.45	9.58	7.88	79.9	
O4054-04	2	12	1.17	8.57	9.74	8.04	80.2	
O4054-05	3	13	1.13	8.76	9.89	8.22	80.9	
O4054-06	4	14	1.16	8.56	9.72	8.00	79.9	
O4054-07	5	15	1.11	8.80	9.91	8.34	82.2	
O4055-01	SS051P-MW03-SB01-00G	39	1.18	8.55	9.73	8.71	88.1	
O4055-02	SS051P-MW01D-SB01-00G	40	1.14	8.64	9.78	8.69	87.4	
O4056-01	P001-PCBULK-001-01	16	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O4056-02	P001-PCBULK-001-02	17	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O4056-03	P001-PCBULK-002-01	18	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O4056-04	P001-SS004-0002-01	19	1.12	8.45	9.57	6.82	67.5	
O4056-05	P001-SS009-0002-01	20	1.14	8.70	9.84	7.08	68.3	
O4056-06	P001-SS010-0002-01	21	1.15	8.80	9.95	8.07	78.6	
O4056-07	P001-SS011-0002-01	22	1.18	8.46	9.64	6.99	68.7	
O4056-08	P001-SS012-0002-01	23	1.12	8.45	9.57	7.05	70.2	
O4056-09	P001-SS013-0002-01	24	1.17	8.53	9.7	7.29	71.7	
O4056-10	P001-SS014-0002-01	25	1.19	8.72	9.91	8.48	83.6	
O4056-11	P001-SS015-0002-01	26	1.15	8.36	9.51	6.98	69.7	



PERCENT SOLID

Supervisor: Sohil
Analyst: jignesh
Date: 8/18/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 08/17/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:20
Out Date: 08/18/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB126897

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O4056-12	P001-SS016-0002-01	27	1.19	8.55	9.74	8.12	81.1	
O4056-13	P001-SS017-0002-01	28	1.19	8.46	9.65	8.12	81.9	
O4056-14	P001-SS019-0001-01	29	1.17	8.58	9.75	6.76	65.2	
O4056-15	O4056-14MS	30	1.17	8.58	9.75	6.76	65.2	
O4056-16	O4056-14MSD	31	1.17	8.58	9.75	6.76	65.2	
O4056-17	P001-SS019-0001-02	32	1.15	8.62	9.77	6.81	65.7	
O4056-18	P001-SS022-0002-01	33	1.15	8.79	9.94	6.49	60.8	
O4056-19	P001-SS023-0002-01	34	1.15	8.66	9.81	8.34	83.0	
O4056-20	P001-SS024-0002-01	35	1.15	8.40	9.55	7.75	78.6	
O4056-21	P001-SS025-0002-01	36	1.19	8.53	9.72	7.68	76.1	
O4056-22	P001-SS026-0002-01	37	1.19	8.45	9.64	8.19	82.8	
O4056-23	P001-SS027-0002-01	38	1.19	8.65	9.84	8.38	83.1	
O4057-01	P001-SS028-0002-01	41	1.18	8.46	9.64	8.57	87.4	
O4057-02	P001-SS031-0002-01	42	1.15	8.75	9.9	7.25	69.7	
O4057-03	P001-SS033-0002-01	43	1.15	8.57	9.72	8.17	81.9	
O4057-04	P001-SS034-0002-01	44	1.14	8.48	9.62	7.7	77.4	
O4057-05	P001-SS035-0002-01	45	1.14	8.80	9.94	8.08	78.9	
O4057-06	P001-SS040-0002-01	46	1.14	8.70	9.84	6.9	66.2	
O4057-07	P001-SS042-0002-01	47	1.19	8.78	9.97	7.36	70.3	
O4057-08	P001-SS047-0002-01	48	1.17	8.74	9.91	5.22	46.3	
O4057-09	P001-SS049-0002-01	49	1.16	8.49	9.65	5.7	53.5	
O4057-10	P001-SS053-0002-01	50	1.15	8.83	9.98	6.66	62.4	
O4057-11	P001-SS056-0002-01	51	1.14	8.57	9.71	7.89	78.8	
O4057-12	P001-SS058-0002-01	52	1.15	8.78	9.93	7.28	69.8	
O4057-13	P001-SS061-0002-01	53	1.17	8.60	9.77	7.65	75.3	
O4057-14	P001-SS062-0002-01	54	1.14	8.84	9.98	7.24	69.0	
O4057-15	P001-SS064-0002-01	55	1.17	8.67	9.84	7.66	74.9	
O4057-16	P001-SS065-0002-01	56	1.15	8.71	9.86	8.04	79.1	



PERCENT SOLID

Supervisor: Sohil
Analyst: jignesh
Date: 8/18/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 08/17/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:20
Out Date: 08/18/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB126897

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O4057-17	P001-SS066-0002-01	57	1.12	8.64	9.76	7.67	75.8	
O4057-18	O4057-17MS	58	1.12	8.64	9.76	7.67	75.8	
O4057-19	O4057-17MSD	59	1.12	8.64	9.76	7.67	75.8	
O4057-20	P001-SS066-0002-02	60	1.14	8.80	9.94	7.83	76.0	
O4057-21	P001-SS067-0002-01	61	1.11	8.63	9.74	7.38	72.7	
O4057-22	P001-SS068-0002-01	62	1.19	8.43	9.62	8.09	81.9	
O4057-23	P001-SS071-0002-01	63	1.13	8.47	9.6	8.31	84.8	
O4058-01	SU-04-081723	67	1.15	8.37	9.52	8.86	92.1	
O4058-02	SU-04-081723-EPH-2	68	1.19	8.53	9.72	9.4	96.2	
O4060-01	AUD-1002	69	1.17	8.59	9.76	9.38	95.6	
O4069-01	CONC-WC-1	64	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
O4070-01	WC-1	65	1.15	8.84	9.99	8.39	81.9	
O4070-02	WC-1-EPH-2	66	1.15	8.84	9.99	8.39	81.9	
O4071-01	SB-1	70	1.17	8.62	9.79	9.74	99.4	
O4071-02	SB-2	71	1.15	8.39	9.54	9.5	99.5	
O4071-03	SB-3	72	1.16	8.55	9.71	9.66	99.4	
O4071-04	SB-4	73	1.16	8.76	9.92	9.78	98.4	
O4071-05	SB-5	74	1.13	8.67	9.8	9.75	99.4	
O4071-06	SB-6	75	1.18	8.49	9.67	9.6	99.2	
O4071-07	SB-7	76	1.17	8.68	9.85	9.83	99.8	
O4071-08	SB-8	77	1.19	8.62	9.81	9.75	99.3	
O4071-09	SB-9	78	1.18	8.50	9.68	9.62	99.3	
O4071-10	SB-10	79	1.15	8.40	9.55	9.28	96.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$