



PERCENT SOLID

Supervisor: Sohil
Analyst: jignesh
Date: 8/22/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 08/21/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:14
Out Date: 08/22/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB126941

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O4090-01	TP-I	1	1.14	8.78	9.92	8.27	81.2	
O4090-03	TP-I-EPH-2	2	1.12	8.82	9.94	8.19	80.2	
O4090-04	TP-K	3	1.19	8.51	9.7	8.47	85.5	
O4090-06	TP-K-EPH-2	4	1.16	8.82	9.98	9.12	90.2	
O4090-07	TP-L	5	1.12	8.87	9.99	8.65	84.9	
O4090-09	TP-L-EPH-2	6	1.19	8.79	9.98	8.68	85.2	
O4090-10	TP-N	7	1.12	8.82	9.94	8.83	87.4	
O4090-12	TP-N-EPH-2	8	1.14	8.81	9.95	8.82	87.2	
O4091-01	TP-A	9	1.18	8.50	9.68	8.35	84.4	
O4091-03	TP-A-EPH-2	10	1.19	8.47	9.66	8.52	86.5	
O4091-04	TP-B	11	1.19	8.57	9.76	8.31	83.1	
O4091-06	TP-B-EPH-2	12	1.11	8.81	9.92	8.53	84.2	
O4091-07	TP-C	13	1.16	8.81	9.97	8.35	81.6	
O4091-09	TP-C-EPH-2	14	1.17	8.59	9.76	8.89	89.9	
O4091-10	TP-E	15	1.13	8.53	9.66	8.71	88.9	
O4091-12	TP-E-EPH-2	16	1.15	8.81	9.96	8.67	85.4	
O4091-13	TP-H	17	1.16	8.81	9.97	9.24	91.7	
O4091-15	TP-H-EPH-2	18	1.12	8.65	9.77	8.62	86.7	
O4091-16	TP-G	19	1.19	8.79	9.98	9.24	91.6	
O4091-18	TP-G-EPH-2	20	1.15	8.41	9.56	8.69	89.7	
O4091-19	TP-F	21	1.12	8.80	9.92	8.78	87.0	
O4091-21	TP-F-EPH-2	22	1.1	8.76	9.86	8.84	88.4	
O4093-01	NEVINS-SB04-Z45-46	23	1.18	8.48	9.66	8.14	82.1	
O4093-02	NEVINS-SB04-Z52-53	24	1.13	8.69	9.82	8.71	87.2	
O4093-03	NEVINS-SB04-Z83-84	25	1.15	8.82	9.97	8.19	79.8	
O4093-04	NEVINS-SB04-Z94-95	26	1.19	8.78	9.97	9.6	95.8	
O4098-01	P080923RB01	27	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O4098-02	P080923RB02	28	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE



PERCENT SOLID

Supervisor: Sohil
Analyst: jignesh
Date: 8/22/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 08/21/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:14
Out Date: 08/22/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB126941

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O4099-01	CONC-WC-1	29	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
O4100-01	WC-1	30	1.14	8.49	9.63	8.43	85.9	
O4100-02	WC-1-EPH-2	31	1.14	8.49	9.63	8.43	85.9	
O4103-01	OR-2-082123	32	1.14	8.47	9.61	9.14	94.5	
O4103-02	OR-2-082123-EPH-2	33	1.15	8.53	9.68	9.22	94.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$