



PERCENT SOLID

Supervisor: Sohil
Analyst: jignesh
Date: 8/29/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 08/28/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 08/29/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB127042

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O4189-01	SB-02 (22-23)	1	1.12	8.85	9.97	5.13	45.3	
O4189-02	SB-06 (16-18)	2	1.12	8.39	9.51	7.96	81.5	
O4189-03	SB-07 (14-15)	3	1.13	8.70	9.83	8.67	86.7	
O4189-04	SB-01 (48-50)	4	1.14	8.76	9.9	8.42	83.1	
O4196-01	OR-34-COMP-01	5	1.15	8.71	9.86	9.08	91.0	
O4196-03	OR-345-VOC-01	6	1.15	8.71	9.86	9.08	91.0	
O4196-05	OR-345-01	7	1.15	8.65	9.8	9.34	94.7	
O4196-06	OR-345-02	8	1.15	8.84	9.99	9.27	91.9	
O4196-07	OR-345-03	9	1.15	8.80	9.95	9.01	89.3	
O4196-08	OR-345-COMP-02	10	1.14	8.62	9.76	9.02	91.4	
O4196-10	OR-345-VOC-02	11	1.14	8.62	9.76	9.02	91.4	
O4196-12	OR-345-04	12	1.13	8.84	9.97	9.34	92.9	
O4196-13	OR-345-05	13	1.14	8.49	9.63	8.95	92.0	
O4196-14	OR-345-06	14	1.16	8.71	9.87	9.14	91.6	
O4196-15	OR-345-COMP-03	15	1.15	8.68	9.83	9.19	92.6	
O4196-17	OR-345-VOC-03	16	1.15	8.68	9.83	9.19	92.6	
O4196-19	OR-345-07	17	1.14	8.82	9.96	9.14	90.7	
O4196-20	OR-345-08	18	1.15	8.47	9.62	8.46	86.3	
O4196-21	OR-345-09	19	1.14	8.69	9.83	8.73	87.3	
O4196-22	OR-345-COMP-04	20	1.14	8.50	9.64	9.22	95.1	
O4196-24	OR-345-VOC-04	21	1.15	8.49	9.64	9.22	95.1	
O4196-26	OR-345-10	22	1.15	8.49	9.64	9.39	97.1	
O4196-27	OR-345-11	23	1.14	8.83	9.97	9.58	95.6	
O4196-28	OR-345-12	24	1.15	8.49	9.64	9.21	94.9	
O4196-29	OR-345-COMP-05	25	1.15	8.49	9.64	9.21	94.9	
O4196-31	OR-345-VOC-05	26	1.15	8.49	9.64	9.21	94.9	
O4196-33	OR-345-13	27	1.14	8.74	9.88	9.11	91.2	
O4196-34	OR-345-14	28	1.11	8.85	9.96	9.06	89.8	



PERCENT SOLID

Supervisor: Sohil
Analyst: jignesh
Date: 8/29/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 08/28/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 08/29/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB127042

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O4196-35	OR-345-15	29	1.12	8.65	9.77	8.72	87.9	
O4197-01	OR-309-OK-370-COMP-20	30	1.1	8.72	9.82	9.5	96.3	
O4197-03	OR-309-OK-370-VOC-20	31	1.1	8.72	9.82	9.5	96.3	
O4197-05	OR-309-OK-370-39	32	1.12	8.79	9.91	9.63	96.8	
O4197-06	OR-309-OK-370-40	33	1.13	8.74	9.87	9.45	95.2	
O4197-07	OR-309-OK-370-41	34	1.1	8.43	9.53	9.14	95.4	
O4197-08	OR-309-OK-370-COMP-21	35	1.11	8.86	9.97	9.7	97.0	
O4197-10	OR-309-OK-370-VOC-21	36	1.11	8.86	9.97	9.7	97.0	
O4197-12	OR-309-OK-370-42	37	1.12	8.83	9.95	9.59	95.9	
O4197-13	OR-309-OK-370-43	38	1.1	8.68	9.78	9.25	93.9	
O4197-14	OR-309-OK-370-44	39	1.13	8.86	9.99	8.93	88.0	
O4197-15	OR-309-OK-370-COMP-22	40	1.13	8.83	9.96	8.65	85.2	
O4197-17	OR-309-OK-370-VOC-22	41	1.13	8.83	9.96	8.65	85.2	
O4197-19	OR-309-OK-370-45	42	1.12	8.63	9.75	8.56	86.2	
O4197-20	OR-309-OK-370-46	43	1.13	8.55	9.68	8.77	89.4	
O4197-21	OR-309-OK-370-47	44	1.15	8.59	9.74	8.65	87.3	
O4197-22	OR-309-OK-370-COMP-23	45	1.14	8.44	9.58	9.25	96.1	
O4197-24	OR-309-OK-370-VOC-23	46	1.14	8.44	9.58	9.25	96.1	
O4197-26	OR-309-OK-370-48	47	1.14	8.84	9.98	9.62	95.9	
O4197-27	OR-309-OK-370-49	48	1.13	8.55	9.68	9.08	93.0	
O4197-28	OR-309-OK-370-50	49	1.14	8.74	9.88	9.61	96.9	
O4198-01	OR-309-OK-370-COMP-24	50	1.13	8.67	9.8	9.36	94.9	
O4198-03	OR-309-OK-370-VOC-24	51	1.13	8.67	9.8	9.36	94.9	
O4198-05	OR-309-OK-370-51	52	1.14	8.73	9.87	9.02	90.3	
O4198-06	OR-309-OK-370-52	53	1.15	8.69	9.84	9.6	97.2	
O4198-07	OR-309-OK-370-COMP-25	54	1.15	8.75	9.9	9.5	95.4	
O4198-09	OR-309-OK-370-VOC-25	55	1.15	8.75	9.9	9.5	95.4	
O4198-11	OR-309-OK-370-53	56	1.14	8.68	9.82	9.18	92.6	



PERCENT SOLID

Supervisor: Sohil
Analyst: jignesh
Date: 8/29/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 08/28/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 08/29/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB127042

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O4198-12	OR-309-OK-370-54	57	1.15	8.65	9.8	9.43	95.7	
O4198-13	OR-309-OK-370-55	58	1.15	8.65	9.8	9.43	95.7	
O4198-14	OR-309-OK-370-COMP-26	59	1.15	8.73	9.88	9.17	91.9	
O4198-16	OR-309-OK-370-VOC-26	60	1.14	8.74	9.88	9.17	91.9	
O4198-18	OR-309-OK-370-56	61	1.15	8.50	9.65	9.31	96.0	
O4198-19	OR-309-OK-370-57	62	1.12	8.84	9.96	9.25	92.0	
O4199-01	KMA664D-1-1	63	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4199-02	KMA664D-1-2	64	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4200-01	COMP-1	65	1.14	8.40	9.54	8.11	83.0	
O4200-02	COMP-2	66	1.14	8.75	9.89	8.99	89.7	
O4200-03	COMP-3	67	1.16	8.44	9.6	9.15	94.7	
O4200-04	8372	68	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O4202-01	C-1	69	1.14	8.58	9.72	8.04	80.4	
O4202-02	C-2	70	1.15	8.83	9.98	8.71	85.6	
O4202-03	C-3	71	1.13	8.69	9.82	8.53	85.2	
O4203-01	SB-1 (12-12.5)	72	1.13	8.84	9.97	8.38	82.0	
O4203-02	SB-2 (12-12.5)	73	1.12	8.73	9.85	8.1	80.0	
O4203-03	SB-3 (12-12.5)	74	1.14	8.54	9.68	8.47	85.8	
O4203-04	SB-4 (12-12.5)	75	1.15	8.45	9.6	8.49	86.9	
O4203-05	SB-5 (12-12.5)	76	1.14	8.83	9.97	8.62	84.7	
O4203-06	SB-6 (12-12.5)	77	1.16	8.77	9.93	8.54	84.2	
O4203-07	SB-7 (12-12.5)	78	1.12	8.45	9.57	8.25	84.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$