



PERCENT SOLID

Supervisor: Sohil
Analyst: JIGNESH
Date: 8/31/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:35
In Date: 08/30/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 08/31/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB127081

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O4231-01	0706	1	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O4232-01	SBG-1	2	1.15	8.64	9.79	9.5	96.6	
O4232-02	SBG-2	3	1.19	8.70	9.89	9.77	98.6	
O4232-03	SBG-3	4	1.19	8.57	9.76	9.68	99.1	
O4232-04	SBG-4	5	1.15	8.38	9.53	9.43	98.8	
O4232-05	SBG-5	6	1.19	8.46	9.65	9.5	98.2	
O4232-06	SBG-6	7	1.16	8.40	9.56	9.44	98.6	
O4232-07	SBG-7	8	1.12	8.56	9.68	9.6	99.1	
O4232-08	SBG-8	9	1.19	8.38	9.57	9.24	96.1	
O4232-09	SB-11	10	1.19	8.43	9.62	9.56	99.3	
O4232-10	SB-12	11	1.1	8.69	9.79	9.66	98.5	
O4242-01	DEBRIS-CONCRETE-COMP	13	1.00	1.00	2.00	2.00	100.0	DEBRIS-CONCRETE
O4244-02	BIN-BLOCK	12	1.00	1.00	2.00	2.00	100.0	concrete
O4246-01	BC226599LOZ-END-1-1	14	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-02	BC226599LOZ-END-1-2	15	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-03	BC280325LLZ-END-1-1	16	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-04	BC280325LLZ-END-1-2	17	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-05	BC280325LLZ-END-2-1	18	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-06	BC280325LLZ-END-2-2	19	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-07	BC247799LOK-END-1-1	20	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-08	BC247799LOK-END-1-2	21	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-09	BC274271LOK-END-1-1	22	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-10	BC274271LOK-END-1-2	23	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-11	BC274271LOK-END-2-1	24	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-12	BC274271LOK-END-2-2	25	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-13	BC704970LOK-END-1-1	26	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-14	BC704970LOK-END-1-2	27	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-15	BC704970LOK-END-2-1BC704970LOK-END-2-1	28	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: Sohil
Analyst: JIGNESH
Date: 8/31/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:35
In Date: 08/30/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 08/31/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB127081

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O4246-16	BC704970LOK-END-2-2	29	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-17	BC122903LOK-END-1-1	30	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-18	BC122903LOK-END-1-2	31	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-19	BC136902LOK-END-1-1	32	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-20	BC136902LOK-END-1-2	33	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-21	BC136902LOK-END-2-1	34	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-22	BC136902LOK-END-2-2	35	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-23	DBA684K-END-1-1	36	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4246-24	DBA684K-END-1-2	37	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4249-01	SU-02-083023	38	1.13	8.69	9.82	9.14	92.2	
O4249-02	SU-02-083023-EPH-2	39	1.13	8.69	9.82	9.14	92.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$