



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 9/26/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 09/25/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:16
Out Date: 09/26/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB127463

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O4596-01	0817	9	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O4597-01	340	10	1.15	8.82	9.97	8.77	86.4	
O4598-01	CONCRETE-BLOCK	11	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
O4599-01	KMA809L-1-1	12	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4599-02	KMA809L-1-2	13	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4599-03	KMA809L-2-1	14	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4599-04	KMA809L-2-2	15	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4600-01	COMP-1	16	1.00	1.00	2.00	2.00	100.0	debris
O4600-02	COMP-2	17	1.00	1.00	2.00	2.00	100.0	debris
O4600-03	COMP-3	18	1.00	1.00	2.00	2.00	100.0	debris
O4601-01	TP-1	1	1.18	8.58	9.76	8.7	87.6	
O4601-03	TP-1-EPH	2	1.18	8.58	9.76	8.7	87.6	
O4601-04	TP-2	3	1.16	8.48	9.64	8.45	86.0	
O4601-06	TP-2-EPH	4	1.16	8.48	9.64	8.45	86.0	
O4602-01	001	5	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O4602-02	002	6	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O4602-03	003	7	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O4603-01	010	8	1.15	8.56	9.71	8.24	82.8	
O4604-01	TR-04-092523	19	1.14	8.63	9.77	8.62	86.7	
O4604-02	TR-04-092523-E2	20	1.11	8.87	9.98	8.46	82.9	
O4605-01	KMA872E-END-1-1	21	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4605-02	KMA872E-END-1-2	22	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4605-03	KMA872E-END-2-1	23	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4605-04	KMA872E-END-2-2	24	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4605-05	JEC429U-END-1-1	25	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4605-06	JEC429U-END-1-2	26	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4605-07	BC259215L0K-END-1-1	27	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4605-08	BC259215L0K-END-1-2	28	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 9/26/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 09/25/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:16
Out Date: 09/26/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB127463

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O4605-09	BC136902-END-1-1	29	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4605-10	BC136902-END-1-2	30	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4605-11	BC258874L02-END-1-1	31	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4605-12	BC258874L02-END-1-2	32	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4605-13	BC204970LOK-END-1-1	33	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4605-14	BC204970LOK-END-1-2	34	1.00	1.00	2.00	2.00	100.0	PILC REEL
O4606-01	ETGI-327	35	1.19	8.43	9.62	8.9	91.5	
O4606-02	ETGI-327	36	1.16	8.79	9.95	8.34	81.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$