



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 9/27/2023

OVENTEMP IN Celsius(°C): 106
Time IN: 17:25
In Date: 09/26/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 09/27/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB127476

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O4612-01	01-A-01-B-01-C	1	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O4612-02	02-A-02-B-02-C	2	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O4612-03	03-A-03-B-03-C	3	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O4612-04	04-A-04-B-04-C	4	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O4612-05	05-A-05-B-05-C	5	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O4612-06	06-A-06-B-06-C	6	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O4612-07	07-A-07-B-07-C	7	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O4612-08	08-A-08-B-08-C	8	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O4613-01	1029	9	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O4613-02	1030	10	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O4613-03	1031	11	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O4613-04	1018	12	1.13	8.44	9.57	7.4	74.3	
O4613-05	1034	13	1.00	1.00	2.00	2.00	100.0	silica gel
O4614-01	AUD-1053	14	1.15	8.43	9.58	8.94	92.4	
O4614-02	AUD-1055	15	1.13	8.54	9.67	9.65	99.8	
O4614-03	AUD-1048	16	1.19	8.45	9.64	9.09	93.5	
O4614-04	AUD-1049	17	1.16	8.66	9.82	9.33	94.3	
O4614-05	AUD-1056	18	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O4615-01	AUD-40279	19	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O4615-02	AUD-40280	20	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O4615-03	AUD-40281	21	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O4615-04	AUD-40282	22	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O4615-05	AUD-40283	23	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O4615-06	AUD-40284	24	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O4615-07	AUD-40285	25	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O4615-08	AUD-40286	26	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O4619-01	TRE-100	27	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O4619-02	TRE-101	28	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 9/27/2023

OVENTEMP IN Celsius(°C): 106
Time IN: 17:25
In Date: 09/26/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 09/27/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB127476

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O4619-03	TRE-102	29	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O4620-01	CON-PAD	30	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
O4620-03	CON-WALL	31	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
O4620-05	CON-SUPPORT-A	32	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
O4620-07	CON-SUPPORT-B	33	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
O4620-11	STONE-IN-E2	34	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
O4620-14	STONE-OUT-E2	35	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
O4621-01	ETGI-287	36	1.17	8.59	9.76	8.61	86.6	
O4621-02	ETGI-287	37	1.11	8.66	9.77	8.64	87.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$