



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/12/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 10/11/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 10/12/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB127732

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O4725-02	K242-IB	24	1.12	8.73	9.85	7.3	70.8	
O4725-15	K242-6C	25	1.18	8.42	9.6	8.55	87.5	
O4725-16	K242-6B	26	1.12	8.76	9.88	8.56	84.9	
O4725-17	K242-7C	27	1.12	8.70	9.82	8.35	83.1	
O4725-18	K242-7B	28	1.16	8.47	9.63	8.6	87.8	
O4725-19	K242-8C	29	1.14	8.81	9.95	8.45	83.0	
O4726-01	K242-9C	30	1.14	8.78	9.92	8.12	79.5	
O4726-03	K242-10C	31	1.16	8.63	9.79	8.53	85.4	
O4726-05	K242-11C	32	1.12	8.61	9.73	8.43	84.9	
O4726-08	K242-12C	33	1.15	8.50	9.65	8.62	87.9	
O4726-13	K242-14C	34	1.17	8.60	9.77	8.36	83.6	
O4726-16	K242-15C	35	1.16	8.80	9.96	9.04	89.5	
O4726-19	K242-16C	36	1.16	8.50	9.66	8.48	86.1	
O4727-02	K242-17C	37	1.14	8.41	9.55	8.27	84.8	
O4727-05	K242-18C	38	1.13	8.73	9.86	8.79	87.7	
O4727-08	K242-19C	39	1.15	8.80	9.95	8.62	84.9	
O4727-09	K242-19B	40	1.13	8.70	9.83	9.44	95.5	
O4727-10	K242-20C	41	1.13	8.73	9.86	8.76	87.4	
O4727-11	K242-20B	42	1.13	8.64	9.77	8.86	89.5	
O4727-12	K242-21C	43	1.16	8.81	9.97	8.99	88.9	
O4727-13	K242-21B	44	1.16	8.80	9.96	8.77	86.5	
O4727-15	K242-22C	45	1.14	8.80	9.94	8.67	85.6	
O4727-16	K242-22B	46	1.11	8.78	9.89	8.76	87.1	
O4727-18	K242-23C	47	1.15	8.81	9.96	8.79	86.7	
O4727-19	K242-23B	48	1.15	8.67	9.82	8.6	85.9	
O4728-01	K242-24C	49	1.19	8.36	9.55	7.87	79.9	
O4728-02	K242-24B	50	1.16	8.74	9.9	8.68	86.0	
O4728-03	K242-25C	51	1.13	8.46	9.59	8.06	81.9	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/12/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 10/11/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 10/12/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB127732

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O4728-04	K242-25B	52	1.17	8.57	9.74	8.52	85.8	
O4729-04	K242-33C	53	1.15	8.40	9.55	8.03	81.9	
O4729-13	K242-36C	54	1.11	8.66	9.77	8.66	87.2	
O4729-20	K242-34C	55	1.13	8.55	9.68	8.4	85.0	
O4730-07	K242-38C	56	1.16	8.50	9.66	8.63	87.9	
O4730-23	DUP-9	57	1.18	8.61	9.79	8.44	84.3	
O4841-01	G-1-(5-10)	1	1.18	8.53	9.71	8.44	85.1	
O4841-02	M-3(0-5)	2	1.15	8.37	9.52	8.92	92.8	
O4842-01	G-1(5-10)	3	1.00	1.00	2.00	2.00	100.0	
O4842-02	M-3(0-5)	4	1.00	1.00	2.00	2.00	100.0	
O4844-01	OR-603-COMP-05	5	1.19	8.59	9.78	9.28	94.2	
O4844-03	OR-603-11	6	1.18	8.53	9.71	9.09	92.7	
O4844-04	OR-603-12	7	1.16	8.79	9.95	9.00	89.2	
O4844-05	OR-603-COMP-06	8	1.12	8.75	9.87	9.06	90.7	
O4844-07	OR-603-13	9	1.16	8.47	9.63	8.9	91.4	
O4844-08	OR-603-14	10	1.19	8.62	9.81	9.12	92.0	
O4844-09	OR-603-COMP-07	11	1.19	8.62	9.81	9.04	91.1	
O4844-11	OR-603-15	12	1.18	8.77	9.95	9.28	92.4	
O4844-12	OR-603-16	13	1.19	8.48	9.67	8.8	89.7	
O4844-13	OR-603-COMP-08	14	1.19	8.49	9.68	8.83	90.0	
O4844-15	OR-603-17	15	1.19	8.42	9.61	8.86	91.1	
O4844-16	OR-603-18	16	1.12	8.47	9.59	8.77	90.3	
O4846-01	WASTE	17	1.19	8.66	9.85	8.77	87.5	
O4854-01	TR-05-101123	58	1.1	8.77	9.87	9.16	91.9	
O4854-02	TR-05-101123-E2	59	1.17	8.75	9.92	9.17	91.4	
O4858-01	WC-1	18	1.13	8.73	9.86	8.66	86.3	
O4858-02	WC-1-EPH	19	1.13	8.73	9.86	8.66	86.3	
O4858-04	WC-2	20	1.17	8.60	9.77	8.77	88.4	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/12/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 10/11/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 10/12/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB127732

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O4858-05	WC-2-EPH	21	1.17	8.60	9.77	8.77	88.4	
O4858-07	WC-3	22	1.16	8.81	9.97	8.84	87.2	
O4858-08	WC-3-EPH	23	1.16	8.81	9.97	8.84	87.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$