



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 11/9/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:35
In Date: 11/08/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:28
Out Date: 11/09/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB128226

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O5279-01	TP-1	1	1.13	8.45	9.58	8.81	90.9	
O5279-02	TP-1-EPH	2	1.16	8.80	9.96	9.14	90.7	
O5279-03	TP-1-VOC	3	1.19	8.57	9.76	8.93	90.3	
O5279-05	MH-1	4	1.18	8.44	9.62	8.15	82.6	
O5279-06	MH-1-EPH	5	1.19	8.73	9.92	7.5	72.3	
O5279-07	MH-1-VOC	6	1.16	8.65	9.81	8.39	83.6	
O5285-01	001	7	1.15	8.82	9.97	4.98	43.4	
O5286-01	001	8	1.15	8.82	9.97	9.41	93.7	
O5286-02	002	9	1.15	8.82	9.97	9.41	93.7	
O5293-01	SB-1	10	1.15	8.78	9.93	9.21	91.8	
O5293-02	SB-2	11	1.16	8.48	9.64	8.29	84.1	
O5293-03	SB-3	12	1.18	8.45	9.63	8.52	86.9	
O5293-04	SB-4	13	1.18	8.70	9.88	8.82	87.8	
O5293-05	SB-5	14	1.15	8.82	9.97	9.1	90.1	
O5293-06	SB-6	15	1.15	8.77	9.92	9.06	90.2	
O5293-07	SB-7	16	1.19	8.78	9.97	9.34	92.8	
O5293-08	SB-8	17	1.16	8.48	9.64	7.98	80.4	
O5294-01	WT-1	18	1.18	8.56	9.74	8.8	89.0	
O5294-02	WT-2	19	1.19	8.79	9.98	7.11	67.3	
O5294-03	WT-3	20	1.15	8.46	9.61	7.32	72.9	
O5294-04	WT-4	21	1.19	8.63	9.82	9.06	91.2	
O5294-05	WT-5	22	1.13	8.72	9.85	9.28	93.5	
O5294-06	WT-6	23	1.18	8.48	9.66	8.68	88.4	
O5295-01	TP-2	24	1.16	8.80	9.96	8.27	80.8	
O5295-02	TP-2-EPH	25	1.19	8.56	9.75	8.65	87.1	
O5295-03	TP-2-VOC	26	1.18	8.81	9.99	9.12	90.1	
O5300-01	01-A-01-B-01-C	27	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O5300-02	02-A-02-B-02-C	28	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE



PERCENT SOLID

Supervisor: Iwona

Analyst: JIGNESH

Date: 11/9/2023

OVENTEMP IN Celsius(°C): 107

Time IN: 17:35

In Date: 11/08/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:28

Out Date: 11/09/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB128226

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O5311-01	TP-1	29	1.15	8.81	9.96	9.2	91.4	
O5311-02	TP-1-EPH	30	1.18	8.58	9.76	8.92	90.2	
O5311-04	TP-2	31	1.15	8.35	9.5	8.57	88.9	
O5311-05	TP-2-EPH	32	1.17	8.78	9.95	9.3	92.6	
O5311-07	TP-3	33	1.16	8.59	9.75	9.02	91.5	
O5311-08	TP-3-EPH	34	1.15	8.62	9.77	8.99	91.0	
O5311-10	TP-4	35	1.16	8.72	9.88	9.09	90.9	
O5311-11	TP-4-EPH	36	1.12	8.74	9.86	9.14	91.8	
O5311-13	TP-5	37	1.18	8.79	9.97	9.19	91.1	
O5311-14	TP-5-EPH	38	1.18	8.63	9.81	9.08	91.5	
O5311-16	TP-6	39	1.18	8.81	9.99	9.44	93.8	
O5311-17	TP-6-EPH	40	1.15	8.40	9.55	8.29	85.0	
O5314-01	01	41	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O5314-02	02	42	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O5314-03	03	43	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O5314-04	04	44	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O5314-05	05	45	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O5314-06	06	46	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O5314-07	07	47	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O5314-08	08	48	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O5314-09	09	49	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O5317-05	D3308	50	1.00	1.00	2.00	2.00	100.0	debris
O5317-06	D3317	51	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O5320-01	MW-03	52	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O5321-01	P001-RC01-01	53	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
O5321-02	O5321-01MS	54	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
O5321-03	O5321-01MSD	55	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
O5321-04	P001-RC01-02	56	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS



PERCENT SOLID

Supervisor: Iwona

Analyst: JIGNESH

Date: 11/9/2023

OVENTEMP IN Celsius(°C): 107

Time IN: 17:35

In Date: 11/08/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:28

Out Date: 11/09/2023

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB128226

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O5321-05	P001-RC02-01	57	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
O5321-06	P001-RC03-01	58	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
O5321-07	P001-RC04-01	59	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
O5321-08	P001-RC05-01	60	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
O5321-09	P001-RC06-01	61	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
O5321-10	P001-RC07-01	62	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
O5321-11	P001-RC08-01	63	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
O5321-12	P001-RC09-01	64	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
O5322-01	KMH370T-1-1	65	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5322-02	KMH370T-1-2	66	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5322-03	KMH777E-1-1	67	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5322-04	KMH777E-1-2	68	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5322-05	KMH389A-1-1	69	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5322-06	KMH389A-1-2	70	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5322-07	HEF218I-1-1	71	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5322-08	HEF218I-1-2	72	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5322-09	BC247601-1-1	73	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5322-10	BC247601-1-2	74	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5322-11	KMA985F-1-1	75	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5322-12	KMA985F-1-2	76	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5322-13	ECA281F-1-1	77	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5322-14	ECA281F-1-2	78	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5322-15	KMA9375-1-1	79	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5322-16	KMA9375-1-1	80	1.00	1.00	2.00	2.00	100.0	PILC REEL

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$