



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 11/13/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 11/10/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 11/11/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB128275

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O5359-01	48-WC-1	1	1.15	8.47	9.62	7.72	77.6	
O5360-01	MASTIC	2	1.00	1.00	2.00	2.00	100.0	MASTIC
O5360-02	TILE	3	1.00	1.00	2.00	2.00	100.0	TILE
O5361-01	OK-01-111023	4	1.16	8.49	9.65	8.98	92.1	
O5361-02	OK-01-111023-E2	5	1.11	8.67	9.78	9.11	92.3	
O5374-01	RB21197	34	1.14	8.53	9.67	8.6	87.5	
O5374-02	RB21197	35	1.13	8.78	9.91	8.7	86.2	
O5376-01	204	36	1.18	8.63	9.81	8.88	89.2	
O5376-03	210	37	1.19	8.73	9.92	9.04	89.9	
O5376-05	223	38	1.14	8.82	9.96	9.52	95.0	
O5377-01	360	39	1.18	8.45	9.63	8.52	86.9	
O5377-02	360	40	1.16	8.82	9.98	8.78	86.4	
O5377-03	289	41	1.17	8.59	9.76	8.74	88.1	
O5377-04	289	42	1.13	8.50	9.63	8.52	86.9	
O5378-01	MH-3	16	1.15	8.80	9.95	9.28	92.4	
O5378-02	MH-3-EPH	17	1.15	8.39	9.54	8.75	90.6	
O5378-03	MH-3-VOC	18	1.15	8.39	9.54	8.75	90.6	
O5379-01	WC-9	19	1.18	8.75	9.93	9.16	91.2	
O5379-02	WC-9-EPH	20	1.12	8.85	9.97	8.87	87.6	
O5379-03	WC-9-VOC	21	1.12	8.85	9.97	8.87	87.6	
O5379-05	WC-4	22	1.15	8.81	9.96	8.66	85.2	
O5379-06	WC-4-EPH	23	1.18	8.79	9.97	8.6	84.4	
O5379-07	WC-4-VOC	24	1.18	8.79	9.97	8.6	84.4	
O5379-09	WC-5	25	1.16	8.47	9.63	8.23	83.5	
O5379-10	WC-5-EPH	26	1.12	8.74	9.86	8.42	83.5	
O5379-11	WC-5-VOC	27	1.12	8.74	9.86	8.42	83.5	
O5380-01	WC-7	28	1.15	8.80	9.95	8.4	82.4	
O5380-02	WC-7-EPH	29	1.17	8.77	9.94	8.51	83.7	



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 11/13/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 11/10/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 11/11/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB128275

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O5380-03	WC-7-VOC	30	1.17	8.77	9.94	8.51	83.7	
O5380-05	WC-8	31	1.15	8.82	9.97	8.85	87.3	
O5380-06	WC-8-EPH	32	1.12	8.73	9.85	8.61	85.8	
O5380-07	WC-8-VOC	33	1.12	8.73	9.85	8.61	85.8	
O5381-05	SVOC-GPC-BLANK	6	1.00	1.00	2.00	2.00	100.0	
O5381-06	PEST-GPC-BLANK	7	1.00	1.00	2.00	2.00	100.0	
O5381-07	PEST-GPC-BLANK-SPIKE	8	1.00	1.00	2.00	2.00	100.0	
O5381-08	PCB-GPC-BLANK	9	1.00	1.00	2.00	2.00	100.0	
O5381-09	PCB-GPC-BLANK-SPIKE	10	1.00	1.00	2.00	2.00	100.0	
O5381-10	SVOC-GPC2-BLANK	11	1.00	1.00	2.00	2.00	100.0	
O5381-11	PEST-GPC2-BLANK	12	1.00	1.00	2.00	2.00	100.0	
O5381-12	PEST-GPC2-BLANK-SPIKE	13	1.00	1.00	2.00	2.00	100.0	
O5381-13	PCB-GPC-BLANK-SPIKE	14	1.00	1.00	2.00	2.00	100.0	
O5381-14	PCB-GCP2-BLANK-SPIKE	15	1.00	1.00	2.00	2.00	100.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$