



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 12/1/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 11/30/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:33
Out Date: 12/01/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB128506

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O5570-01	68X-15A-0-2	22	1.15	8.84	9.99	8.68	85.2	
O5570-02	68X-15A-2-4	23	1.12	8.45	9.57	8.48	87.1	
O5572-01	BC136902-LOK-END-1-1	1	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5572-02	BC136902-LOK-END-1-2	2	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5572-03	39811-END-1-1	3	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5572-04	39811-END-1-2	4	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5572-05	39811-END-2-1	5	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5572-06	39811-END-2-2	6	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5573-01	BLUE-VARIOUS-OIL-DRUM	7	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O5573-03	2943	8	1.12	8.84	9.96	7.00	66.5	
O5573-05	2945	9	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O5573-07	2946	10	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O5573-09	2947	11	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O5578-01	BC283069-LOK-END-1-1	12	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5578-02	BC283069-LOK-END-1-2	13	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5578-03	BC271332-LOZ-END-1-1	14	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5578-04	BC271332-LOZ-END-1-2	15	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5578-05	BC271332-LOZ-END-2-1	16	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5578-06	BC271332-LOZ-END-2-2	17	1.00	1.00	2.00	2.00	100.0	PILC REEL
O5579-01	3344	18	1.00	1.00	2.00	2.00	100.0	DEBRIS
O5579-05	3330	19	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O5580-01	218	20	1.00	1.00	2.00	2.00	100.0	debris
O5580-03	222	21	1.00	1.00	2.00	2.00	100.0	debris
O5586-01	VAULT-1	34	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
O5586-02	VAULT-1-EPH2	35	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
O5587-01	WE-MR-CON	36	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
O5588-01	TUDOR	24	1.18	8.61	9.79	8.53	85.4	
O5589-01	DRUM-A-OH-OILY-STONE	25	1.00	1.00	2.00	2.00	100.0	oily-stone



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 12/1/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 11/30/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:33
Out Date: 12/01/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB128506

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O5589-02	DRUM-B-OH-OIL	26	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O5590-01	SH-SOIL-A	27	1.17	8.60	9.77	8.28	82.7	
O5590-02	SH-SOIL-A-EPH2	28	1.19	8.59	9.78	8.53	85.4	
O5590-03	SH-SOIL-B	29	1.16	8.83	9.99	8.89	87.5	
O5590-04	SH-SOIL-B-EPH2	30	1.15	8.82	9.97	8.45	82.8	
O5591-01	33586-33587	31	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O5595-01	OR-03-113023	32	1.19	8.40	9.59	8.88	91.5	
O5595-02	OR-03-113023-EPH2	33	1.16	8.61	9.77	9.41	95.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$