



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 1/5/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:40
In Date: 01/04/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 01/05/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB129020

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P1015-01	MR-301P-0-2	1	1.15	8.81	9.96	8.46	83.0	
P1015-03	MR-301P-25-27	2	1.19	8.53	9.72	8.26	82.9	
P1018-01	001	3	1.12	8.46	9.58	8.47	86.9	
P1018-02	002	4	1.12	8.46	9.58	8.47	86.9	
P1027-01	BC247799-1-1	5	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1027-02	BC247799-1-2	6	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1027-03	BC24857-1-1	7	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1027-04	BC24857-1-2	8	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1029-01	225	9	1.00	1.00	2.00	2.00	100.0	debris
P1030-01	3366-3370	10	1.00	1.00	2.00	2.00	100.0	debris
P1030-02	3374	11	1.00	1.00	2.00	2.00	100.0	debris
P1031-01	KMH364M-1-1	12	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1031-02	KMH364M-1-2	13	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1031-03	BC270896-1-1	14	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1031-04	BC270896-1-2	15	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1031-05	KMA5830-1-1	16	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1031-06	KMA5830-1-2	17	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1031-07	BC260544-1-1	18	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1031-08	BC260544-1-2	19	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1032-01	378	20	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
P1032-02	379	21	1.00	1.00	2.00	2.00	100.0	debris
P1032-03	380	22	1.00	1.00	2.00	2.00	100.0	debris
P1033-01	HIB153Q-A	23	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1033-02	HIB153Q-B	24	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1033-03	W029022-A	25	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1033-04	W029022-B	26	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1033-05	BC226505-A	27	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1033-06	BC226505-B	28	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 1/5/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:40
In Date: 01/04/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 01/05/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB129020

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P1037-04	END-1	29	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1037-05	END-2	30	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1040-01	COMP-1	31	1.15	8.40	9.55	8.06	82.3	
P1040-02	COMP-2	32	1.13	8.68	9.81	7.71	75.8	
P1040-03	COMP-3	33	1.19	8.57	9.76	7.45	73.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$