



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 1/8/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 01/05/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 01/06/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB129048

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P1039-01	TREMONT-DRUM-27-POLY	1	1.00	1.00	2.00	2.00	100.0	plastic sample
P1039-02	TREMONT-DRUM-28-PPE	2	1.00	1.00	2.00	2.00	100.0	plastic sample
P1041-05	SVOC-GPC-BLANK	3	1.00	1.00	2.00	2.00	100.0	
P1041-06	PEST-GPC-BLANK	4	1.00	1.00	2.00	2.00	100.0	
P1041-07	PEST-GPC-BLANK-SPIKE	5	1.00	1.00	2.00	2.00	100.0	
P1041-08	PCB-GPC-BLANK	6	1.00	1.00	2.00	2.00	100.0	
P1041-09	PCB-GPC-BLANK-SPIKE	7	1.00	1.00	2.00	2.00	100.0	
P1041-10	SVOC-GPC2-BLANK	8	1.00	1.00	2.00	2.00	100.0	
P1041-11	PEST-GPC2-BLANK	9	1.00	1.00	2.00	2.00	100.0	
P1041-12	PEST-GPC2-BLANK-SPIKE	10	1.00	1.00	2.00	2.00	100.0	
P1041-13	PCB-GPC-BLANK-SPIKE	11	1.00	1.00	2.00	2.00	100.0	
P1041-14	PCB-GCP2-BLANK-SPIKE	12	1.00	1.00	2.00	2.00	100.0	
P1050-01	CHRT24029	13	1.13	8.76	9.89	9.26	92.8	
P1050-02	CHRT24029-E2	14	1.18	8.43	9.61	8.49	86.7	
P1051-01	340	15	1.19	8.65	9.84	8.74	87.3	
P1051-02	340-E2	16	1.19	8.62	9.81	8.76	87.8	
P1052-01	LEAD-10F6	17	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1052-02	LEAD-20F6	18	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1052-03	LEAD-30F6	19	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1052-04	LEAD-40F6	20	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1052-05	LEAD-50F6	21	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1052-06	LEAD-60F6	22	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1053-01	OILY-DEBRIS	23	1.00	1.00	2.00	2.00	100.0	OILY DEBRIS
P1053-04	OILY-DEBRIS	24	1.00	1.00	2.00	2.00	100.0	OILY DEBRIS
P1054-01	TR-06-01052024	25	1.15	8.38	9.53	8.29	85.2	
P1054-02	TR-06-01052024	26	1.14	8.62	9.76	8.45	84.8	
P1055-02	TOTE-1221	27	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
P1055-04	TOTE-1222	28	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 1/8/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 01/05/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 01/06/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB129048

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P1055-06	TOTE-1223	29	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
P1055-07	SPILL-DEBRIS-COMP	30	1.14	8.72	9.86	9.4	94.7	
P1056-01	MR-311P-0-2	31	1.12	8.76	9.88	8.7	86.5	
P1056-03	MR-311P-30-32	32	1.18	8.80	9.98	9.55	95.1	
P1057-01	CLI-2024-1A	33	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1057-02	CLI-2024-1B	34	1.00	1.00	2.00	2.00	100.0	PILC REEL

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$