



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 1/17/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 01/16/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 01/17/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB129159

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O5831-07	CHRT-20017	1	1.19	8.58	9.77	8.13	80.9	
P1150-01	MR-306-0-2	2	1.14	8.77	9.91	8.63	85.4	
P1150-03	MR-306-25-27	3	1.16	8.50	9.66	8.06	81.2	
P1152-03	275	4	1.15	8.40	9.55	8.34	85.6	
P1152-04	275-E2	5	1.16	8.81	9.97	8.68	85.4	
P1153-01	PAD-A	6	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
P1153-02	PAD-B	7	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
P1154-01	60039	8	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
P1155-01	BB-TESTPITS	9	1.16	8.80	9.96	8.55	84.0	
P1155-02	BB-TESTPITS-E2	10	1.18	8.81	9.99	8.49	83.0	
P1156-01	SB-09-9.0-9.5	11	1.14	8.82	9.96	9.22	91.6	
P1156-02	SB-09-19.5-20.0	12	1.13	8.82	9.95	8.3	81.3	
P1156-05	SB-8-9.0-9.5	13	1.13	8.70	9.83	8.73	87.4	
P1156-07	SB-8-19.5-20.0	14	1.18	8.76	9.94	8.83	87.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$