



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 1/22/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 16:40
In Date: 01/19/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 01/20/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB129183

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P1183-01	MR-305-0-2	1	1.19	8.63	9.82	7.86	77.3	
P1183-03	MR-305-25-27	2	1.15	8.54	9.69	7.91	79.2	
P1185-01	SB-7-19.5-20.0	3	1.15	8.83	9.98	8.35	81.5	
P1185-02	SB-7-9.5-10.0	4	1.16	8.50	9.66	9.14	93.9	
P1185-03	SB-6-10.5-11.0	5	1.15	8.82	9.97	9.36	93.1	
P1185-04	SB-6-19.5-20.0	6	1.13	8.54	9.67	7.74	77.4	
P1187-01	LOD-MDL-SOIL-01-QT1-2024	7	1.00	1.00	2.00	2.00	100.0	
P1187-02	LOQ-SOIL-02-QT1-2024	8	1.00	1.00	2.00	2.00	100.0	
P1187-03	MDL-SOIL-03-QT1-2024	9	1.00	1.00	2.00	2.00	100.0	
P1187-04	LOD-MDL-SOIL-04-QT1-2024	10	1.00	1.00	2.00	2.00	100.0	
P1187-05	LOQ-SOIL-05-QT1-2024	11	1.00	1.00	2.00	2.00	100.0	
P1187-06	MDL-SOIL-06-QT1-2024	12	1.00	1.00	2.00	2.00	100.0	
P1187-16	LLOQ-SOIL-QT1-2024	13	1.00	1.00	2.00	2.00	100.0	
P1188-01	SOUTH-JERSEY-PORT-SAMPLE-SET-2-A	14	1.15	8.80	9.95	8.77	86.6	
P1188-02	SOUTH-JERSEY-PORT-SAMPLE-SET-2-B	15	1.19	8.51	9.7	8.36	84.3	
P1197-05	SVOC-GPC-BLANK	16	1.00	1.00	2.00	2.00	100.0	
P1197-06	PEST-GPC-BLANK	17	1.00	1.00	2.00	2.00	100.0	
P1197-07	PEST-GPC-BLANK-SPIKE	18	1.00	1.00	2.00	2.00	100.0	
P1197-08	PCB-GPC-BLANK	19	1.00	1.00	2.00	2.00	100.0	
P1197-09	PCB-GPC-BLANK-SPIKE	20	1.00	1.00	2.00	2.00	100.0	
P1197-10	SVOC-GPC2-BLANK	21	1.00	1.00	2.00	2.00	100.0	
P1197-11	PEST-GPC2-BLANK	22	1.00	1.00	2.00	2.00	100.0	
P1197-12	PEST-GPC2-BLANK-SPIKE	23	1.00	1.00	2.00	2.00	100.0	
P1197-13	PCB-GPC2-BLANK	24	1.00	1.00	2.00	2.00	100.0	
P1197-14	PCB-GPC2-BLANK-SPIKE	25	1.00	1.00	2.00	2.00	100.0	
P1200-01	CHRT-24376	26	1.13	8.64	9.77	2.17	12.0	GEL MATRIX
P1200-02	CHRT-24376-E2	27	1.16	8.80	9.96	2.85	19.2	GEL MATRIX



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 1/22/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 16:40
In Date: 01/19/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 01/20/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB129183

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P1201-01	3393	28	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
P1201-02	3394	29	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
P1201-03	3395	30	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
P1202-01	KMA695Y-1-1	31	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1202-02	KMA695Y-1-2	32	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1202-03	KMA695Y-2-1	33	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1202-04	KMA695Y-2-2	34	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1202-05	KZZ375Z-1-1	35	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1202-06	KZZ375Z-1-2	36	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1202-07	JEH420H-1-1	37	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1202-08	JEH420H-1-2	38	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1202-09	3385-1	39	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1202-10	3385-2	40	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1202-11	3385-3	41	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1202-12	3385-4	42	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1202-13	3385-5	43	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1202-14	3385-6	44	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1202-15	3386-1	45	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1202-16	3386-2	46	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1202-17	3386-3	47	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1202-18	3386-4	48	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1202-19	3386-5	49	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1202-20	3386-6	50	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1203-01	72-11977	51	1.16	8.81	9.97	9.15	90.7	
P1203-02	72-11977-E2	52	1.16	8.83	9.99	9.13	90.3	
P1204-01	KMA895G-1-1	53	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1204-02	KMA895G-1-2	54	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1206-01	Y23080069-1-1	55	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 1/22/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 16:40
In Date: 01/19/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 01/20/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB129183

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P1206-02	Y23080069-1-2	56	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1206-03	Y23080069-2-1	57	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1206-04	Y23080069-2-2	58	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1206-05	Y23080077-1-1	59	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1206-06	Y23080077-1-2	60	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1206-07	ECA392U-1-1	61	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1206-08	ECA392U-1-2	62	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1206-09	ECA262B-1-1	63	1.00	1.00	2.00	2.00	100.0	PILC REEL
P1206-10	ECA262B-1-2	64	1.00	1.00	2.00	2.00	100.0	PILC REEL

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$