



PERCENT SOLID

Supervisor: Iwona
Analyst: Rubina
Date: 2/1/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 18:20
In Date: 01/31/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:45
Out Date: 02/01/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB129322

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P1304-01	BRIGHTON-WEST-8TH-STREET	8	1.12	8.64	9.76	8.97	90.9	
P1304-02	BRIGHTON-WEST-8TH-STREET	1	1.13	8.86	9.99	9.32	92.4	
P1306-01	SS-01	2	1.15	8.83	9.98	9.77	97.6	
P1307-01	TR-614-COMP-01	3	1.16	8.73	9.89	8.83	87.9	
P1307-04	TR-614-VOC-01	4	1.16	8.73	9.89	8.83	87.9	
P1307-05	TR-614-01	5	1.19	8.77	9.96	8.93	88.3	
P1307-06	TR-614-02	6	1.19	8.78	9.97	9.14	90.5	
P1307-07	TR-614-03	7	1.16	8.63	9.79	8.73	87.7	
P1307-08	TR-614-COMP-02	9	1.15	8.81	9.96	9.37	93.3	
P1307-11	TR-614-VOC-02	10	1.15	8.81	9.96	9.37	93.3	
P1307-12	TR-614-04	11	1.19	8.55	9.74	8.84	89.5	
P1307-13	TR-614-05	12	1.19	8.80	9.99	9.44	93.7	
P1307-14	TR-614-06	13	1.15	8.40	9.55	8.4	86.3	
P1307-15	TR-614-COMP-03	14	1.16	8.82	9.98	8.86	87.3	
P1307-18	TR-614-VOC-03	15	1.16	8.82	9.98	8.86	87.3	
P1307-19	TR-614-07	16	1.17	8.60	9.77	8.26	82.4	
P1307-20	TR-614-08	17	1.18	8.79	9.97	8.85	87.3	
P1307-21	TR-614-09	18	1.18	8.44	9.62	8.53	87.1	
P1311-01	TP-1	23	1.14	8.58	9.72	8.36	84.1	
P1311-02	TP-1	24	1.14	8.58	9.72	8.36	84.1	
P1311-03	TP-2	25	1.15	8.81	9.96	8.8	86.8	
P1311-04	TP-2	26	1.15	8.81	9.96	8.8	86.8	
P1311-05	TP-3	27	1.15	8.81	9.96	8.58	84.3	
P1311-06	TP-3	28	1.15	8.81	9.96	8.58	84.3	
P1311-07	TP-4	29	1.11	8.39	9.5	7.99	82.0	
P1311-08	TP-4	30	1.11	8.39	9.5	7.99	82.0	
P1312-03	19690	37	1.00	1.00	2.00	2.00	100.0	oil sample



PERCENT SOLID

Supervisor: Iwona
Analyst: Rubina
Date: 2/1/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 18:20
In Date: 01/31/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:45
Out Date: 02/01/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB129322

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P1312-05	20720	38	1.00	1.00	2.00	2.00	100.0	oil sample
P1313-01	COMP-1	19	1.15	8.79	9.94	8.33	81.7	
P1313-02	COMP-2	20	1.12	8.40	9.52	7.6	77.1	
P1314-05	MOO-23-00071	39	1.00	1.00	2.00	2.00	100.0	debris
P1314-07	MOO-23-00101	43	1.00	1.00	2.00	2.00	100.0	debris
P1314-09	MOO-23-00084	40	1.00	1.00	2.00	2.00	100.0	oily-debris
P1315-01	KZZ386K-END-1-1	31	1.00	1.00	2.00	2.00	100.0	oilc
P1315-02	KZZ386K-END-1-2	32	1.00	1.00	2.00	2.00	100.0	oilc
P1315-03	BC278948LOK-END-1-1	33	1.00	1.00	2.00	2.00	100.0	oilc
P1315-04	BC278948LOK-END-1-2	34	1.00	1.00	2.00	2.00	100.0	oilc
P1315-05	W029020-END-1-1	35	1.00	1.00	2.00	2.00	100.0	oilc
P1315-06	W029020-END-1-2	36	1.00	1.00	2.00	2.00	100.0	oilc
P1316-01	COMP-1	21	1.19	8.39	9.58	7.7	77.6	
P1316-02	COMP-2	22	1.19	8.64	9.83	7.79	76.4	
P1317-01	NB-7-013124	41	1.14	8.46	9.6	8.31	84.8	
P1317-02	NB-7-013124-E2	42	1.11	8.51	9.62	8.39	85.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$