



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 2/7/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 02/06/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 02/07/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB129405

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P1378-01	TP-4	2	1.16	8.70	9.86	8.58	85.3	
P1378-03	TP-5	3	1.13	8.74	9.87	8.9	88.9	
P1379-01	RE137-CARBON-20240205	1	1.18	8.39	9.57	6.94	68.7	
P1390-01	ETGI-357	4	1.00	1.00	2.00	2.00	100.0	stone sample, 100 % solids
P1390-02	ETGI-357	5	1.00	1.00	2.00	2.00	100.0	stone sample, 100 % solids
P1391-01	SP-1	6	1.15	8.80	9.95	9.16	91.0	
P1391-02	SP-1-A	7	1.15	8.81	9.96	9.14	90.7	
P1391-03	SP-1-B	8	1.14	8.83	9.97	9.05	89.6	
P1391-04	SP-1-C	9	1.13	8.77	9.9	9.00	89.7	
P1391-05	SP-1-D	10	1.15	8.45	9.6	8.82	90.8	
P1391-06	SP-1-E	11	1.13	8.55	9.68	8.59	87.3	
P1391-07	SP-1-F	12	1.18	8.59	9.77	8.96	90.6	
P1391-08	SP-1-G	13	1.13	8.68	9.81	8.97	90.3	
P1391-09	SP-1-H	14	1.15	8.82	9.97	9.1	90.1	
P1391-11	SP-2	15	1.12	8.72	9.84	8.91	89.3	
P1391-12	SP-2-A	16	1.19	8.42	9.61	8.72	89.4	
P1391-13	SP-2-B	17	1.18	8.50	9.68	8.8	89.6	
P1391-14	SP-2-C	18	1.15	8.81	9.96	9.06	89.8	
P1391-15	SP-2-D	19	1.11	8.87	9.98	9.24	91.7	
P1391-16	SP-2-E	20	1.15	8.79	9.94	8.75	86.5	
P1391-17	SP-2-F	21	1.18	8.47	9.65	8.71	88.9	
P1391-18	SP-2-G	22	1.19	8.78	9.97	9.16	90.8	
P1391-19	SP-2-H	23	1.15	8.79	9.94	8.8	87.0	
P1391-21	SP-3	24	1.17	8.58	9.75	8.88	89.9	
P1391-22	SP-3-A	25	1.19	8.40	9.59	8.7	89.4	
P1391-23	SP-3-B	26	1.11	8.47	9.58	8.46	86.8	
P1391-24	SP-3-C	27	1.17	8.81	9.98	8.8	86.6	
P1391-25	SP-3-D	28	1.18	8.44	9.62	8.67	88.7	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 2/7/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 02/06/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 02/07/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB129405

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P1391-26	SP-3-E	29	1.13	8.80	9.93	8.9	88.3	
P1391-27	SP-3-F	30	1.18	8.79	9.97	9.31	92.5	
P1391-28	SP-3-G	31	1.19	8.32	9.51	8.63	89.4	
P1391-29	SP-3-H	32	1.11	8.48	9.59	8.61	88.4	
P1391-31	SP-4	33	1.18	8.81	9.99	9.26	91.7	
P1391-32	SP-4-A	34	1.18	8.78	9.96	9.06	89.7	
P1391-33	SP-4-B	35	1.13	8.54	9.67	8.67	88.3	
P1391-34	SP-4-C	36	1.18	8.79	9.97	9.17	90.9	
P1391-35	SP-4-D	37	1.16	8.55	9.71	8.87	90.2	
P1391-36	SP-4-E	38	1.19	8.62	9.81	8.83	88.6	
P1391-37	SP-4-F	39	1.18	8.55	9.73	8.34	83.7	
P1391-38	SP-4-G	40	1.15	8.40	9.55	8.54	88.0	
P1391-39	SP-4-H	41	1.17	8.70	9.87	9.00	90.0	
P1391-41	SP-5	42	1.17	8.60	9.77	8.82	89.0	
P1391-42	SP-5-A	43	1.18	8.74	9.92	9.02	89.7	
P1391-43	SP-5-B	44	1.15	8.60	9.75	8.72	88.0	
P1391-44	SP-5-C	45	1.18	8.73	9.91	9.07	90.4	
P1391-45	SP-5-D	46	1.17	8.80	9.97	8.99	88.9	
P1391-46	SP-5-E	47	1.17	8.81	9.98	9.08	89.8	
P1391-47	SP-5-F	48	1.19	8.64	9.83	8.4	83.4	
P1391-48	SP-5-G	49	1.18	8.60	9.78	8.29	82.7	
P1391-49	SP-5-H	50	1.13	8.45	9.58	8.7	89.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$