



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 2/16/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 02/15/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 02/16/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB129546

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P1435-06	FRENCH-020724-CONC	1	1.00	1.00	2.00	2.00	100.0	conc sample,100 & solids
P1467-04	03A-03B-03C	14	1.00	1.00	2.00	2.00	100.0	caluk
P1505-01	WC-1S	2	1.14	8.47	9.61	8.91	91.7	
P1505-02	WC-1S-EPH	3	1.19	8.61	9.8	9.06	91.4	
P1505-03	WC-1S-VOC	4	1.19	8.61	9.8	9.06	91.4	
P1505-05	WC-1D	5	1.17	8.60	9.77	9.04	91.5	
P1505-06	WC-1D-EPH	6	1.14	8.47	9.61	8.91	91.7	
P1505-07	WC-1D-VOC	7	1.14	8.47	9.61	8.91	91.7	
P1505-09	WC-2S	8	1.15	8.50	9.65	8.88	90.9	
P1505-10	WC-2S-EPH	9	1.18	8.42	9.6	8.83	90.9	
P1505-11	WC-2S-VOC	10	1.18	8.42	9.6	8.83	90.9	
P1505-13	WC-2D	11	1.14	8.56	9.7	8.77	89.1	
P1505-14	WC-2D-EPH	12	1.15	8.81	9.96	9.14	90.7	
P1505-15	WC-2D-VOC	13	1.15	8.81	9.96	9.14	90.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$