



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 3/1/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 16:50
In Date: 02/29/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 03/01/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB129721

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P1620-01	SP-1	1	1.18	8.56	9.74	9.00	91.4	
P1620-02	SP-1A	2	1.13	8.80	9.93	9.16	91.3	
P1620-03	SP-1B	3	1.18	8.38	9.56	8.79	90.8	
P1620-04	SP-1C	4	1.19	8.72	9.91	9.02	89.8	
P1620-05	SP-1D	5	1.19	8.60	9.79	9.02	91.0	
P1620-06	SP-1E	6	1.18	8.59	9.77	8.96	90.6	
P1620-07	SP-1F	7	1.15	8.84	9.99	9.27	91.9	
P1620-08	SP-1G	8	1.18	8.43	9.61	8.75	89.8	
P1620-09	SP-1H	9	1.12	8.72	9.84	8.97	90.0	
P1620-11	SP-2	10	1.15	8.82	9.97	9.33	92.7	
P1620-12	SP-2A	11	1.18	8.70	9.88	9.12	91.3	
P1620-13	SP-2B	12	1.14	8.84	9.98	9.11	90.2	
P1620-14	SP-2C	13	1.15	8.40	9.55	8.89	92.1	
P1620-15	SP-2D	14	1.16	8.40	9.56	9.05	93.9	
P1620-16	SP-2E	15	1.14	8.61	9.75	9.14	92.9	
P1620-17	SP-2F	16	1.15	8.83	9.98	9.26	91.8	
P1620-18	SP-2G	17	1.12	8.70	9.82	8.93	89.8	
P1620-19	SP-2H	18	1.14	8.67	9.81	9.04	91.1	
P1620-21	SP-3	19	1.13	8.76	9.89	9.12	91.2	
P1620-22	SP-3A	20	1.18	8.48	9.66	8.88	90.8	
P1620-23	SP-3B	21	1.18	8.47	9.65	9.12	93.7	
P1620-24	SP-3C	22	1.18	8.47	9.65	8.92	91.4	
P1620-25	SP-3D	23	1.18	8.75	9.93	9.28	92.6	
P1620-26	SP-3E	24	1.18	8.80	9.98	9.19	91.0	
P1620-27	SP-3F	25	1.15	8.78	9.93	9.22	91.9	
P1620-28	SP-3G	26	1.18	8.63	9.81	8.94	89.9	
P1620-29	SP-3H	27	1.18	8.43	9.61	8.93	91.9	
P1620-31	SP-4	28	1.19	8.53	9.72	9.12	93.0	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 3/1/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 16:50
In Date: 02/29/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 03/01/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB129721

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P1620-32	SP-4A	29	1.18	8.43	9.61	8.92	91.8	
P1620-33	SP-4B	30	1.18	8.63	9.81	9.18	92.7	
P1620-34	SP-4C	31	1.16	8.70	9.86	9.24	92.9	
P1620-35	SP-4D	32	1.19	8.51	9.7	9.09	92.8	
P1620-36	SP-4E	33	1.16	8.65	9.81	9.19	92.8	
P1620-37	SP-4F	34	1.18	8.48	9.66	9.12	93.6	
P1620-38	SP-4G	35	1.16	8.54	9.7	9.09	92.9	
P1620-39	SP-4H	36	1.19	8.58	9.77	9.13	92.5	
P1639-01	SP-5	37	1.15	8.40	9.55	9.17	95.5	
P1639-02	SP-5A	38	1.16	8.56	9.72	9.39	96.1	
P1639-03	SP-5B	39	1.18	8.46	9.64	9.24	95.3	
P1639-04	SP-5C	40	1.18	8.69	9.87	9.47	95.4	
P1639-05	SP-5D	41	1.14	8.42	9.56	9.18	95.5	
P1639-06	SP-5E	42	1.16	8.57	9.73	9.33	95.3	
P1639-07	SP-5F	43	1.13	8.46	9.59	9.22	95.6	
P1639-08	SP-5G	44	1.12	8.56	9.68	9.32	95.8	
P1639-09	SP-5H	45	1.19	8.66	9.85	9.5	96.0	
P1639-11	SP-6	46	1.12	8.70	9.82	9.32	94.3	
P1639-12	SP-6A	47	1.18	8.50	9.68	9.26	95.1	
P1639-13	SP-6B	48	1.16	8.63	9.79	9.39	95.4	
P1639-14	SP-6C	49	1.18	8.47	9.65	9.26	95.4	
P1639-15	SP-6D	50	1.16	8.65	9.81	9.36	94.8	
P1639-16	SP-6E	51	1.15	8.60	9.75	9.37	95.6	
P1639-17	SP-6F	52	1.19	8.44	9.63	9.19	94.8	
P1639-18	SP-6G	53	1.15	8.45	9.6	9.19	95.1	
P1639-19	SP-6H	54	1.18	8.50	9.68	9.28	95.3	
P1639-21	SP-7	55	1.19	8.33	9.52	8.66	89.7	
P1639-22	SP-7A	56	1.18	8.36	9.54	8.59	88.6	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 3/1/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 16:50
In Date: 02/29/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 03/01/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB129721

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P1639-23	SP-7B	57	1.18	8.48	9.66	8.78	89.6	
P1639-24	SP-7C	58	1.19	8.63	9.82	8.88	89.1	
P1639-25	SP-7D	59	1.15	8.82	9.97	9.16	90.8	
P1639-26	SP-7E	60	1.19	8.51	9.7	8.85	90.0	
P1639-27	SP-7F	61	1.16	8.43	9.59	8.68	89.2	
P1639-28	SP-7G	62	1.16	8.73	9.89	8.96	89.3	
P1639-29	SP-7H	63	1.19	8.56	9.75	8.87	89.7	
P1639-31	SP-8	64	1.16	8.48	9.64	8.52	86.8	
P1639-32	SP-8A	65	1.14	8.80	9.94	8.83	87.4	
P1639-33	SP-8B	66	1.19	8.41	9.6	8.66	88.8	
P1639-34	SP-8C	67	1.13	8.80	9.93	8.86	87.8	
P1639-35	SP-8D	68	1.19	8.80	9.99	8.89	87.5	
P1639-36	SP-8E	69	1.15	8.41	9.56	8.58	88.3	
P1639-37	SP-8F	70	1.14	8.82	9.96	9.08	90.0	
P1639-38	SP-8G	71	1.15	8.80	9.95	8.91	88.2	
P1639-39	SP-8H	72	1.16	8.59	9.75	8.68	87.5	
P1639-41	SP-8I	73	1.11	8.69	9.8	8.81	88.6	
P1653-01	294	74	1.19	8.65	9.84	8.72	87.1	
P1653-02	294-E2	75	1.18	8.51	9.69	8.63	87.5	
P1654-01	BC136902-1-1	76	1.00	1.00	2.00	2.00	100.0	oilc
P1654-02	BC136902-1-2	77	1.00	1.00	2.00	2.00	100.0	oilc
P1655-01	KZZ007Q-1-1	78	1.00	1.00	2.00	2.00	100.0	oilc
P1655-02	KZZ007Q-1-2	79	1.00	1.00	2.00	2.00	100.0	oilc
P1655-03	KZA510C-1-1	80	1.00	1.00	2.00	2.00	100.0	oilc
P1655-04	KZA510C-1-2	81	1.00	1.00	2.00	2.00	100.0	oilc
P1655-05	HZA550T-1-1	82	1.00	1.00	2.00	2.00	100.0	oilc
P1655-06	HZA550T-1-2	83	1.00	1.00	2.00	2.00	100.0	oilc
P1656-01	BC258874-1-1	84	1.00	1.00	2.00	2.00	100.0	oilc



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 3/1/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 16:50
In Date: 02/29/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 03/01/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB129721

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P1656-02	BC258874-1-2	85	1.00	1.00	2.00	2.00	100.0	pilc
P1656-03	BC258874-2-1	86	1.00	1.00	2.00	2.00	100.0	pilc
P1656-04	BC258874-2-2	87	1.00	1.00	2.00	2.00	100.0	pilc
P1656-05	BC259315-1-1	88	1.00	1.00	2.00	2.00	100.0	pilc
P1656-06	BC259315-1-2	89	1.00	1.00	2.00	2.00	100.0	pilc
P1659-01	TP-16	90	1.15	8.53	9.68	8.56	86.9	
P1659-02	TP-16-EPH	91	1.19	8.50	9.69	8.64	87.6	
P1659-05	TP-7	92	1.17	8.58	9.75	8.58	86.4	
P1659-06	TP-7-EPH	93	1.15	8.80	9.95	8.85	87.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$