



PERCENT SOLID

Supervisor: Sohil
Analyst: jignesh
Date: 4/3/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 04/02/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: m oven-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:27
Out Date: 04/03/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: m sc-4
Thermometer ID: % solids-oven

QC:LB130134

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P1954-01	RBR21147	2	1.14	8.60	9.74	8.08	80.7	
P1954-02	RBR21147-E2	3	1.18	8.42	9.6	7.85	79.2	
P1954-03	RBR21145	4	1.15	8.80	9.95	7.25	69.3	
P1954-04	RBR21145-E2	5	1.16	8.49	9.65	6.22	59.6	
P1955-01	CHRT-20103	6	1.15	8.80	9.95	5.07	44.5	
P1955-02	CHRT-20103-E2	7	1.16	8.79	9.95	7.27	69.5	
P1956-02	231	8	1.17	8.60	9.77	8.39	84.0	
P1956-03	236	9	1.00	1.00	2.00	2.00	100.0	debris
P1956-04	RBR200010	10	1.15	8.39	9.54	8.06	82.4	
P1956-05	RBR200010-E2	11	1.13	8.73	9.86	7.9	77.5	
P1957-01	NM-DRUMS-040224	1	1.15	8.81	9.96	7.4	70.9	sludge sample
P1958-01	VNJ-207	12	1.12	8.64	9.76	8.17	81.6	
P1958-02	VNJ-207-E2	13	1.13	8.85	9.98	8.73	85.9	
P1959-01	BC274768-1-1	14	1.00	1.00	2.00	2.00	100.0	PILC
P1959-02	BC274768-1-2	15	1.00	1.00	2.00	2.00	100.0	PILC
P1959-03	EFB316N-1-1	16	1.00	1.00	2.00	2.00	100.0	PILC
P1959-04	EFB316N-1-2	17	1.00	1.00	2.00	2.00	100.0	PILC
P1960-01	3471	18	1.00	1.00	2.00	2.00	100.0	derbis
P1960-02	20071	19	1.16	8.49	9.65	6.69	65.1	
P1960-03	20071-E2	20	1.12	8.85	9.97	5.72	52.0	
P1961-01	GHA466T-1-1	21	1.00	1.00	2.00	2.00	100.0	PILC
P1961-02	GHA466T-1-1	22	1.00	1.00	2.00	2.00	100.0	PILC
P1961-03	BC274857-1-1	23	1.00	1.00	2.00	2.00	100.0	PILC
P1961-04	BC274857-1-2	24	1.00	1.00	2.00	2.00	100.0	PILC
P1961-05	BC107490-1-1	25	1.00	1.00	2.00	2.00	100.0	PILC
P1961-06	BC107490-1-2	26	1.00	1.00	2.00	2.00	100.0	PILC
P1961-07	BC270896-1-1	27	1.00	1.00	2.00	2.00	100.0	PILC
P1961-08	BC270896-1-2	28	1.00	1.00	2.00	2.00	100.0	PILC



PERCENT SOLID

Supervisor: Sohil
Analyst: jignesh
Date: 4/3/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 04/02/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: m oven-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:27
Out Date: 04/03/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: m sc-4
Thermometer ID: % solids-oven

QC:LB130134

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P1961-09	HEH629P-1-1	29	1.00	1.00	2.00	2.00	100.0	PILC
P1961-10	HEH629P-1-2	30	1.00	1.00	2.00	2.00	100.0	PILC
P1962-01	KMH244Y-1-1	31	1.00	1.00	2.00	2.00	100.0	PILC
P1962-02	KMH244Y-1-2	32	1.00	1.00	2.00	2.00	100.0	PILC
P1962-03	Y2308-003-1-1	33	1.00	1.00	2.00	2.00	100.0	PILC
P1962-04	Y2308-003-1-2	34	1.00	1.00	2.00	2.00	100.0	PILC
P1962-05	Y2310-0287-1-1	35	1.00	1.00	2.00	2.00	100.0	PILC
P1962-06	Y2310-0287-1-2	36	1.00	1.00	2.00	2.00	100.0	PILC
P1962-07	Y2308-0013-1-1	37	1.00	1.00	2.00	2.00	100.0	PILC
P1962-08	Y2308-0013-1-2	38	1.00	1.00	2.00	2.00	100.0	PILC
P1962-09	BC204970-1-1	39	1.00	1.00	2.00	2.00	100.0	PILC
P1962-10	BC204970-1-2	40	1.00	1.00	2.00	2.00	100.0	PILC
P1962-11	DBA291H-1-1	41	1.00	1.00	2.00	2.00	100.0	PILC
P1962-12	DBA291H-1-1	42	1.00	1.00	2.00	2.00	100.0	PILC
P1963-01	391	43	1.00	1.00	2.00	2.00	100.0	derbis

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$