

PERCENT SOLID

Supervisor: Sohil
Analyst: jignesh
Date: 4/5/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 14:35
In Date: 04/04/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 04/05/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB130174

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P1993-01	BH2M3	1	1.13	8.59	9.72	9.34	95.6	
P1993-02	BH2M4	2	1.12	8.47	9.59	8.08	82.2	
P1993-03	BH2M5	3	1.15	8.64	9.79	8.51	85.2	
P1993-05	BH2N9	4	1.16	8.48	9.64	8.97	92.1	
P1993-07	BHA35	5	1.13	8.76	9.89	8.98	89.6	
P1993-08	BHA36	6	1.19	8.58	9.77	9.37	95.3	
P1993-09	BHA37	7	1.19	8.52	9.71	9.22	94.2	
P1993-10	BH2L1	8	1.16	8.53	9.69	9.06	92.6	
P1993-11	BH2L2	9	1.18	8.64	9.82	9.4	95.1	
P1993-12	BH2L3	10	1.15	8.81	9.96	8.18	79.8	
P1993-13	BH2L4	11	1.15	8.38	9.53	9.06	94.4	
P1993-14	BH2L5	12	1.13	8.44	9.57	9.41	98.1	
P1993-15	BH2L6	13	1.16	8.44	9.6	8.46	86.5	
P1993-21	VHBLK001	14	1.00	1.00	2.00	2.00	100.0	vhblk

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$