

Analysis Method: 9034
 Parameter: Reactive Sulfide
 Run Number: LB130190

ANALYST: Nikita
 SUPERVISOR REVIEW BY: jignesh
 Constant: 16000
 Normality1: 0.025
 Normality2: 0.025

Reagent/Standard	Lot/Log #
SODIUM THIOSULFATE, 0.025N, 4LITRE	W2976
IODINE SOLUTION .025N 1L	W3050
Starch Solution, 4L	W2977

Seq	Lab ID	True Value (mg/l)	DF	Initial Weight (g)	Final Volume (ml)	T1 (ml)	T2 Initial	T2 Final	T2 Diff. (ml)	T1 - T2 Diff (mL)	Value Corrected With Blank	Result (ppm)	Anal Date	Anal Time
1	PB160039BL		1	5.00	50	2.00	0.00	1.94	1.94	0.06	0.00	0.00	04/05/2024	11:00
2	P1974-03		1	5.04	50	2.00	0.00	1.92	1.92	0.08	0.02	1.59	04/05/2024	11:03
3	P1974-03DUP		1	5.03	50	2.00	0.00	1.92	1.92	0.08	0.02	1.59	04/05/2024	11:06
4	P1974-06		1	5.02	50	2.00	0.00	1.86	1.86	0.14	0.08	6.37	04/05/2024	11:09
5	P1974-09		1	5.01	50	2.00	0.00	1.90	1.90	0.10	0.04	3.19	04/05/2024	11:12
6	P1975-01		1	5.04	50	2.00	0.00	1.86	1.86	0.14	0.08	6.35	04/05/2024	11:15
7	P1995-02		1	5.05	50	2.00	0.00	1.90	1.90	0.10	0.04	3.17	04/05/2024	11:18
8	P1995-04		1	5.02	50	2.00	0.00	1.88	1.88	0.12	0.06	4.78	04/05/2024	11:21
9	P1995-06		1	5.04	50	2.00	0.00	1.90	1.90	0.10	0.04	3.17	04/05/2024	11:24
10	P1996-03		1	5.02	50	2.00	0.00	1.88	1.88	0.12	0.06	4.78	04/05/2024	11:27
11	P1996-06		1	5.03	50	2.00	0.00	1.92	1.92	0.08	0.02	1.59	04/05/2024	11:30
12	P1996-09		1	5.02	50	2.00	0.00	1.88	1.88	0.12	0.06	4.78	04/05/2024	11:33

T1 = Titrant1

T2 = Titrant2

T2 Diff = T2 Final - T2 Initial

Value Corrected With Blank = ((T1 - T2 Diff) - Blank Correction(BL))

Result = ((T1 * Normality1) - ((T1 - Value Corrected With Blank) * Normality2)) * Constant / Initial Volume