

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 4/11/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 04/10/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 04/11/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB130264

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P2061-01	KMA818K-1-1	1	1.00	1.00	2.00	2.00	100.0	oilc
P2061-02	KMA818K-1-2	2	1.00	1.00	2.00	2.00	100.0	oilc
P2062-01	31550	3	1.00	1.00	2.00	2.00	100.0	wipe sample
P2062-02	31551	4	1.00	1.00	2.00	2.00	100.0	wipe sample
P2063-01	BUR-1289	5	1.13	8.57	9.7	9.53	98.0	
P2068-01	AC01-1	13	1.00	1.00	2.00	2.00	100.0	oilc
P2068-02	AC01-2	14	1.00	1.00	2.00	2.00	100.0	oilc
P2068-03	AC01-3	15	1.00	1.00	2.00	2.00	100.0	oilc
P2068-04	AC01-4	16	1.00	1.00	2.00	2.00	100.0	oilc
P2068-05	AC01-5	17	1.00	1.00	2.00	2.00	100.0	oilc
P2068-06	AC02-1	18	1.00	1.00	2.00	2.00	100.0	oilc
P2068-07	AC02-2	19	1.00	1.00	2.00	2.00	100.0	oilc
P2068-08	AC02-3	20	1.00	1.00	2.00	2.00	100.0	oilc
P2068-09	AC02-4	21	1.00	1.00	2.00	2.00	100.0	oilc
P2068-10	AC02-5	22	1.00	1.00	2.00	2.00	100.0	oilc
P2068-11	AC03-1	23	1.00	1.00	2.00	2.00	100.0	oilc
P2068-12	AC03-2	24	1.00	1.00	2.00	2.00	100.0	oilc
P2068-13	AC03-3	25	1.00	1.00	2.00	2.00	100.0	oilc
P2068-14	AC03-4	26	1.00	1.00	2.00	2.00	100.0	oilc
P2068-15	AC03-5	27	1.00	1.00	2.00	2.00	100.0	oilc
P2068-16	AC04-1	28	1.00	1.00	2.00	2.00	100.0	oilc
P2068-17	AC04-2	29	1.00	1.00	2.00	2.00	100.0	oilc
P2068-18	AC04-3	30	1.00	1.00	2.00	2.00	100.0	oilc
P2068-19	AC04-4	31	1.00	1.00	2.00	2.00	100.0	oilc
P2068-20	AC04-5	32	1.00	1.00	2.00	2.00	100.0	oilc
P2068-21	AC05-1	33	1.00	1.00	2.00	2.00	100.0	oilc
P2068-22	AC05-2	34	1.00	1.00	2.00	2.00	100.0	oilc
P2068-23	AC05-3	35	1.00	1.00	2.00	2.00	100.0	oilc

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 4/11/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 04/10/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 04/11/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB130264

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P2068-24	AC05-4	36	1.00	1.00	2.00	2.00	100.0	oilc
P2068-25	AC05-5	37	1.00	1.00	2.00	2.00	100.0	oilc
P2068-26	AC06-1	38	1.00	1.00	2.00	2.00	100.0	oilc
P2068-27	AC06-2	39	1.00	1.00	2.00	2.00	100.0	oilc
P2068-28	AC06-3	40	1.00	1.00	2.00	2.00	100.0	oilc
P2068-29	AC06-4	41	1.00	1.00	2.00	2.00	100.0	oilc
P2068-30	AC06-5	42	1.00	1.00	2.00	2.00	100.0	oilc
P2069-01	AC07-1	43	1.00	1.00	2.00	2.00	100.0	oilc
P2069-02	AC07-2	44	1.00	1.00	2.00	2.00	100.0	oilc
P2069-03	AC07-3	45	1.00	1.00	2.00	2.00	100.0	oilc
P2069-04	AC07-4	46	1.00	1.00	2.00	2.00	100.0	oilc
P2069-05	AC07-5	47	1.00	1.00	2.00	2.00	100.0	oilc
P2069-06	AC08-1	48	1.00	1.00	2.00	2.00	100.0	oilc
P2069-07	AC08-2	49	1.00	1.00	2.00	2.00	100.0	oilc
P2069-08	AC08-3	50	1.00	1.00	2.00	2.00	100.0	oilc
P2069-09	AC08-4	51	1.00	1.00	2.00	2.00	100.0	oilc
P2069-10	AC08-5	52	1.00	1.00	2.00	2.00	100.0	oilc
P2069-11	AC09-1	53	1.00	1.00	2.00	2.00	100.0	oilc
P2069-12	AC09-2	54	1.00	1.00	2.00	2.00	100.0	oilc
P2069-13	AC09-3	55	1.00	1.00	2.00	2.00	100.0	oilc
P2069-14	AC09-4	56	1.00	1.00	2.00	2.00	100.0	oilc
P2069-15	AC09-5	57	1.00	1.00	2.00	2.00	100.0	oilc
P2070-01	TP-1	6	1.18	8.64	9.82	8.99	90.4	
P2070-03	TP-2	7	1.18	8.49	9.67	8.76	89.3	
P2070-05	TP-3	8	1.17	8.60	9.77	9.16	92.9	
P2070-07	TP-4	9	1.11	8.75	9.86	8.61	85.7	
P2070-09	TP-5	10	1.15	8.68	9.83	9.3	93.9	
P2071-01	PL-01-041024	11	1.19	8.38	9.57	8.79	90.7	

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 4/11/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 04/10/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 04/11/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB130264

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P2071-02	PL-01-041024-E2	12	1.15	8.84	9.99	9.06	89.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$