

**PERCENT SOLID**

Supervisor: Iwona  
Analyst: jignesh  
Date: 4/18/2024

OVENTEMP IN Celsius(°C): 107  
Time IN: 17:05  
In Date: 04/17/2024  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103  
Time OUT: 08:10  
Out Date: 04/18/2024  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
BalanceID: M SC-4  
Thermometer ID: % SOLID- OVEN

QC:LB130367

| Lab ID   | Client SampleID | Dish # | Dish Wt (g) (A) | Sample Wt (g) | Dish + Sample Wt (g) (B) | Dish+Dry Sample Wt (g) (C) | % Solid | Comments |
|----------|-----------------|--------|-----------------|---------------|--------------------------|----------------------------|---------|----------|
| P2169-01 | YCTD7           | 1      | 1.19            | 8.55          | 9.74                     | 8.36                       | 83.9    |          |
| P2169-02 | YCTD8           | 2      | 1.16            | 8.50          | 9.66                     | 7.8                        | 78.1    |          |
| P2169-03 | YCTD9           | 3      | 1.19            | 8.67          | 9.86                     | 9.04                       | 90.5    |          |
| P2169-04 | YCTE0           | 4      | 1.15            | 8.77          | 9.92                     | 9.05                       | 90.1    |          |
| P2169-05 | YCTE1           | 5      | 1.19            | 8.49          | 9.68                     | 8.56                       | 86.8    |          |
| P2169-06 | YCTE2           | 6      | 1.15            | 8.67          | 9.82                     | 8.79                       | 88.1    |          |
| P2169-07 | YCTE3           | 7      | 1.15            | 8.38          | 9.53                     | 8.35                       | 85.9    |          |
| P2169-08 | YCTE3D          | 8      | 1.15            | 8.38          | 9.53                     | 8.35                       | 85.9    |          |
| P2169-09 | YCTE3S          | 9      | 1.15            | 8.38          | 9.53                     | 8.35                       | 85.9    |          |
| P2169-10 | YCTE4           | 10     | 1.16            | 8.42          | 9.58                     | 8.84                       | 91.2    |          |
| P2169-11 | YCMR6           | 11     | 1.16            | 8.50          | 9.66                     | 9.09                       | 93.3    |          |
| P2169-12 | YCMR8           | 12     | 1.15            | 8.40          | 9.55                     | 8.85                       | 91.7    |          |
| P2169-13 | YCMR9           | 13     | 1.14            | 8.45          | 9.59                     | 9.15                       | 94.8    |          |
| P2169-14 | YCMS0           | 14     | 1.18            | 8.60          | 9.78                     | 9.22                       | 93.5    |          |
| P2169-15 | YCTE5           | 15     | 1.16            | 8.51          | 9.67                     | 8.77                       | 89.4    |          |
| P2169-16 | YCTE6           | 16     | 1.15            | 8.66          | 9.81                     | 8.57                       | 85.7    |          |
| P2169-17 | YCTE7           | 17     | 1.18            | 8.67          | 9.85                     | 8.73                       | 87.1    |          |
| P2169-18 | YCTE8           | 18     | 1.16            | 8.44          | 9.6                      | 8.4                        | 85.8    |          |
| P2169-19 | YCTF5           | 19     | 1.17            | 8.58          | 9.75                     | 9.4                        | 95.9    |          |

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$