

**PERCENT SOLID**

Supervisor: Iwona  
Analyst: jignesh  
Date: 4/30/2024

OVENTEMP IN Celsius(°C): 107  
Time IN: 17:00  
In Date: 04/29/2024  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103  
Time OUT: 08:11  
Out Date: 04/30/2024  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
BalanceID: M SC-4  
Thermometer ID: % SOLID- OVEN

QC:LB130545

| Lab ID   | Client SampleID | Dish # | Dish Wt (g) (A) | Sample Wt (g) | Dish + Sample Wt (g) (B) | Dish+Dry Sample Wt (g) (C) | % Solid | Comments                 |
|----------|-----------------|--------|-----------------|---------------|--------------------------|----------------------------|---------|--------------------------|
| P2307-01 | WC-1-VOC        | 9      | 1.13            | 8.62          | 9.75                     | 8.58                       | 86.4    |                          |
| P2307-02 | WC-1            | 10     | 1.12            | 8.46          | 9.58                     | 8.62                       | 88.7    |                          |
| P2307-04 | WC-2-VOC        | 11     | 1.14            | 8.83          | 9.97                     | 8.5                        | 83.4    |                          |
| P2307-05 | WC-2            | 12     | 1.15            | 8.40          | 9.55                     | 8.54                       | 88.0    |                          |
| P2307-07 | WC-3-VOC        | 13     | 1.16            | 8.69          | 9.85                     | 8.8                        | 87.9    |                          |
| P2307-08 | WC-3            | 14     | 1.13            | 8.51          | 9.64                     | 8.73                       | 89.3    |                          |
| P2308-01 | 346             | 1      | 1.14            | 8.57          | 9.71                     | 9.36                       | 95.9    |                          |
| P2308-03 | 72-12016        | 2      | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   | conc sample,100 & solids |
| P2309-01 | TR-05-042924    | 3      | 1.14            | 8.80          | 9.94                     | 9.46                       | 94.5    |                          |
| P2309-02 | TR-05-042924-E2 | 4      | 1.18            | 8.45          | 9.63                     | 8.5                        | 86.6    |                          |
| P2310-01 | SU-02-042924    | 5      | 1.18            | 8.51          | 9.69                     | 9.39                       | 96.5    |                          |
| P2310-02 | SU-02-042924-E2 | 6      | 1.16            | 8.54          | 9.7                      | 9.00                       | 91.8    |                          |
| P2310-03 | SU-03-042924    | 7      | 1.15            | 8.40          | 9.55                     | 9.22                       | 96.1    |                          |
| P2310-04 | SU-03-042924-E2 | 8      | 1.11            | 8.74          | 9.85                     | 9.33                       | 94.1    |                          |

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$