

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 5/6/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 05/03/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 05/04/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB130610

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P2241-05	344	1	1.00	1.00	2.00	2.00	100.0	conc sample,100 & solids
P2362-01	INSIDE	2	1.15	8.40	9.55	9.14	95.1	
P2362-02	OUTSIDE-FOOTINGS	3	1.19	8.52	9.71	9.00	91.7	
P2366-01	AT024P-SD02-043024-00	14	1.15	8.81	9.96	7.56	72.8	
P2366-02	AT024P-SD01-043024-00	15	1.15	8.81	9.96	5.91	54.0	
P2366-03	AT024P-SD05-043024-00	16	1.19	8.51	9.7	7.63	75.7	
P2366-04	AT024P-SD03-043024-00	17	1.14	8.83	9.97	7.35	70.3	
P2366-05	AT024P-SD04-050124-00	18	1.16	8.47	9.63	7.45	74.3	
P2366-06	SS051P-SD04-050124-00	19	1.15	8.35	9.5	7.6	77.2	
P2367-01	SS051P-SD01-050224-00	20	1.18	8.80	9.98	8.07	78.3	
P2367-02	SS051P-SD03-050224-00	21	1.15	8.83	9.98	8.16	79.4	
P2367-03	SS051P-SD02-050224-00	22	1.19	8.73	9.92	8.07	78.8	
P2367-04	SS051P-SD12-050224-00	23	1.13	8.84	9.97	8.12	79.1	
P2367-05	SS051P-SD11-050224-00	24	1.1	8.75	9.85	7.95	78.3	
P2377-05	SVOC-GPC-BLANK	4	1.00	1.00	2.00	2.00	100.0	
P2377-06	PEST-GPC-BLANK	5	1.00	1.00	2.00	2.00	100.0	
P2377-07	PEST-GPC-BLANK-SPIKE	6	1.00	1.00	2.00	2.00	100.0	
P2377-08	PCB-GPC-BLANK	7	1.00	1.00	2.00	2.00	100.0	
P2377-09	PCB-GPC-BLANK-SPIKE	8	1.00	1.00	2.00	2.00	100.0	
P2377-10	SVOC-GPC2-BLANK	9	1.00	1.00	2.00	2.00	100.0	
P2377-11	PEST-GPC2-BLANK	10	1.00	1.00	2.00	2.00	100.0	
P2377-12	PEST-GPC2-BLANK-SPIKE	11	1.00	1.00	2.00	2.00	100.0	
P2377-13	PCB-GPC2-BLANK	12	1.00	1.00	2.00	2.00	100.0	
P2377-14	PCB-GCP2-BLANK-SPIKE	13	1.00	1.00	2.00	2.00	100.0	
P2378-01	OR-02-05032024	25	1.16	8.50	9.66	8.8	89.9	
P2378-02	OR-02-05032024	26	1.12	8.70	9.82	8.74	87.6	
P2384-01	CABLE-1	27	1.00	1.00	2.00	2.00	100.0	pile
P2384-02	CABLE-2	28	1.00	1.00	2.00	2.00	100.0	pile

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 5/6/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 05/03/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 05/04/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB130610

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P2384-03	CABLE-3	29	1.00	1.00	2.00	2.00	100.0	oilc
P2384-04	CABLE-4	30	1.00	1.00	2.00	2.00	100.0	oilc
P2384-05	CABLE-5	31	1.00	1.00	2.00	2.00	100.0	oilc
P2384-06	CABLE-6	32	1.00	1.00	2.00	2.00	100.0	oilc
P2384-07	CABLE-7	33	1.00	1.00	2.00	2.00	100.0	oilc
P2384-08	CABLE-8	34	1.00	1.00	2.00	2.00	100.0	oilc
P2384-09	CABLE-9	35	1.00	1.00	2.00	2.00	100.0	oilc
P2384-10	CABLE-10	36	1.00	1.00	2.00	2.00	100.0	oilc
P2384-11	CABLE-11	37	1.00	1.00	2.00	2.00	100.0	oilc
P2384-12	CABLE-12	38	1.00	1.00	2.00	2.00	100.0	oilc
P2384-13	CABLE-13	39	1.00	1.00	2.00	2.00	100.0	oilc
P2384-14	CABLE-14	40	1.00	1.00	2.00	2.00	100.0	oilc
P2384-15	CABLE-15	41	1.00	1.00	2.00	2.00	100.0	oilc
P2384-16	CABLE-16	42	1.00	1.00	2.00	2.00	100.0	oilc
P2384-17	CABLE-17	43	1.00	1.00	2.00	2.00	100.0	oilc
P2384-18	CABLE-18	44	1.00	1.00	2.00	2.00	100.0	oilc
P2384-19	CABLE-19	45	1.00	1.00	2.00	2.00	100.0	oilc
P2384-20	CABLE-20	46	1.00	1.00	2.00	2.00	100.0	oilc
P2384-21	CABLE-21	47	1.00	1.00	2.00	2.00	100.0	oilc
P2384-22	CABLE-22	48	1.00	1.00	2.00	2.00	100.0	oilc
P2384-23	CABLE-23	49	1.00	1.00	2.00	2.00	100.0	oilc
P2384-24	CABLE-24	50	1.00	1.00	2.00	2.00	100.0	oilc
P2384-25	CABLE-25	51	1.00	1.00	2.00	2.00	100.0	oilc
P2385-01	WALL-02	52	1.00	1.00	2.00	2.00	100.0	wipe sample
P2385-02	WALL-03	53	1.00	1.00	2.00	2.00	100.0	wipe sample
P2385-03	WALL-04	54	1.00	1.00	2.00	2.00	100.0	wipe sample
P2385-04	WALL-05	55	1.00	1.00	2.00	2.00	100.0	wipe sample
P2385-05	CEILING-02	56	1.00	1.00	2.00	2.00	100.0	wipe sample

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 5/6/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 05/03/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 05/04/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB130610

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P2385-06	CEILING-03	57	1.00	1.00	2.00	2.00	100.0	wipe sample
P2385-07	FLOOR-02	58	1.00	1.00	2.00	2.00	100.0	wipe sample
P2385-08	FLOOR-03	59	1.00	1.00	2.00	2.00	100.0	wipe sample
P2385-09	PIPE-01	60	1.00	1.00	2.00	2.00	100.0	wipe sample
P2385-10	PIPE-02	61	1.00	1.00	2.00	2.00	100.0	wipe sample
P2385-11	STAIRS-02	62	1.00	1.00	2.00	2.00	100.0	wipe sample
P2385-12	BEAM-02	63	1.00	1.00	2.00	2.00	100.0	wipe sample
P2386-01	EO-01-050324-A	64	1.18	8.38	9.56	8.57	88.2	
P2386-02	EO-01-050324-A-E2	65	1.15	8.36	9.51	8.53	88.3	
P2386-03	EO-01-050324-B	66	1.17	8.54	9.71	8.97	91.3	
P2386-04	EO-01-050324-B-E2	67	1.16	8.67	9.83	9.06	91.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$