

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 5/9/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 05/08/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 05/09/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB130669

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P2411-01	WCS-TPA	1	1.16	8.80	9.96	9.4	93.6	
P2411-02	WCS-TPA-EPH	2	1.12	8.40	9.52	8.87	92.3	
P2411-03	WCS-TPA-VOC	3	1.12	8.50	9.62	9.00	92.7	
P2412-01	SS-1-COMP	4	1.18	8.45	9.63	8.8	90.2	
P2412-02	SS-1-GRAB	5	1.13	8.64	9.77	8.96	90.6	
P2413-01	WC-CC	6	1.11	8.72	9.83	7.58	74.2	
P2413-02	WC-CC-E2	7	1.13	8.55	9.68	7.47	74.2	
P2413-03	WC-CC-VOC	8	1.18	8.77	9.95	8.98	88.9	
P2413-05	WC-EE	9	1.14	8.40	9.54	7.64	77.4	
P2413-06	WC-EE-E2	10	1.17	8.60	9.77	8.24	82.2	
P2413-07	WC-EE-VOC	11	1.15	8.84	9.99	9.36	92.9	
P2414-01	1011UMR-SS004-0212-01	12	1.19	8.50	9.69	8.21	82.6	
P2415-01	1011UMR-SS005-0212-01	13	1.19	8.53	9.72	8.11	81.1	
P2421-01	CABLE-1A	14	1.00	1.00	2.00	2.00	100.0	oilc
P2421-02	CABLE-2A	15	1.00	1.00	2.00	2.00	100.0	oilc
P2421-03	CABLE-3A	16	1.00	1.00	2.00	2.00	100.0	oilc
P2421-04	CABLE-4A	17	1.00	1.00	2.00	2.00	100.0	oilc
P2421-05	CABLE-5A	18	1.00	1.00	2.00	2.00	100.0	oilc
P2421-06	CABLE-6A	19	1.00	1.00	2.00	2.00	100.0	oilc
P2421-07	CABLE-7A	20	1.00	1.00	2.00	2.00	100.0	oilc
P2421-08	CABLE-8A	21	1.00	1.00	2.00	2.00	100.0	oilc
P2421-09	CABLE-9A	22	1.00	1.00	2.00	2.00	100.0	oilc
P2421-10	CABLE-10A	23	1.00	1.00	2.00	2.00	100.0	oilc
P2421-11	CABLE-11A	24	1.00	1.00	2.00	2.00	100.0	oilc
P2421-12	CABLE-12A	25	1.00	1.00	2.00	2.00	100.0	oilc
P2421-13	CABLE-13A	26	1.00	1.00	2.00	2.00	100.0	oilc
P2421-14	CABLE-17A	27	1.00	1.00	2.00	2.00	100.0	oilc
P2421-15	CABLE-21A	28	1.00	1.00	2.00	2.00	100.0	oilc

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 5/9/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 05/08/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 05/09/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB130669

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P2421-16	CABLE-22A	29	1.00	1.00	2.00	2.00	100.0	oilc
P2421-17	CABLE-24A	30	1.00	1.00	2.00	2.00	100.0	oilc
P2421-18	CABLE25-A	31	1.00	1.00	2.00	2.00	100.0	oilc
P2421-19	CABLE-26	32	1.00	1.00	2.00	2.00	100.0	oilc
P2421-20	CABLE-26-A	33	1.00	1.00	2.00	2.00	100.0	oilc
P2421-21	CABLE-27	34	1.00	1.00	2.00	2.00	100.0	oilc
P2421-22	CABLE-27A	35	1.00	1.00	2.00	2.00	100.0	oilc
P2421-23	CABLE-28	36	1.00	1.00	2.00	2.00	100.0	oilc
P2421-24	CABLE-28A	37	1.00	1.00	2.00	2.00	100.0	oilc
P2421-25	CABLE-29	38	1.00	1.00	2.00	2.00	100.0	oilc
P2421-26	CABLE-29A	39	1.00	1.00	2.00	2.00	100.0	oilc
P2421-27	CABLE-30	40	1.00	1.00	2.00	2.00	100.0	oilc
P2421-28	CABLE-30A	41	1.00	1.00	2.00	2.00	100.0	oilc

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$