

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 5/10/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 05/09/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:08
Out Date: 05/10/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB130686

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P2422-01	WCS-TP-D	1	1.18	8.63	9.81	8.75	87.7	
P2422-02	WCS-TP-D-E2	2	1.15	8.81	9.96	8.81	86.9	
P2422-03	WCS-TP-D-VOC	3	1.12	8.68	9.8	8.66	86.9	
P2422-05	WCS-TP-F	4	1.14	8.40	9.54	8.15	83.5	
P2422-06	WCS-TP-F-E2	5	1.18	8.60	9.78	8.88	89.5	
P2422-07	WCS-TP-F-VOC	6	1.19	8.79	9.98	9.36	92.9	
P2423-01	WC-GG	7	1.12	8.70	9.82	8.37	83.3	
P2423-02	WC-GG-E2	8	1.16	8.68	9.84	8.06	79.5	
P2423-03	WC-GG-VOC	9	1.17	8.81	9.98	8.79	86.5	
P2423-05	WC-HH	10	1.16	8.75	9.91	8.81	87.4	
P2423-06	WC-HH-E2	11	1.11	8.77	9.88	8.67	86.2	
P2423-07	WC-HH-VOC	12	1.18	8.81	9.99	8.89	87.5	
P2430-01	1011UMR-SS001-1224-01	13	1.19	8.52	9.71	7.68	76.2	
P2430-03	1011UMR-SS002-1224-01	14	1.16	8.56	9.72	7.62	75.5	
P2430-05	1011UMR-SS003-1224-01	15	1.19	8.67	9.86	7.98	78.3	
P2431-01	AT029P-SD06-050624-00	17	1.15	8.45	9.6	5.73	54.2	
P2431-02	SS052P-SD04-050624-00	18	1.15	8.81	9.96	8.45	82.9	
P2431-03	SS030P-SD04-050624-00	19	1.15	8.39	9.54	7.58	76.6	
P2431-04	FT013P-SD04-050624-00	20	1.16	8.48	9.64	7.66	76.7	
P2431-05	AT054P-SD04-050624-00	21	1.12	8.82	9.94	8.19	80.2	
P2431-07	SS050P-SD01-050624-00	22	1.18	8.53	9.71	8.02	80.2	
P2431-08	SS050P-SD04-050624-00	23	1.18	8.76	9.94	7.95	77.3	
P2431-09	DA053P-SD02-050624-00	24	1.16	8.81	9.97	8.05	78.2	
P2432-01	B-02-011	25	1.13	8.62	9.75	8.48	85.3	
P2432-02	B-01-011	26	1.19	8.50	9.69	8.68	88.1	
P2432-03	SW-01-011	27	1.18	8.50	9.68	9.11	93.3	
P2432-04	SW-02-011	28	1.11	8.47	9.58	8.99	93.0	
P2432-05	SW-03-011	29	1.18	8.47	9.65	8.79	89.8	

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 5/10/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 05/09/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:08
Out Date: 05/10/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB130686

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P2432-06	SW-04-011	30	1.11	8.47	9.58	9.04	93.6	
P2432-07	P-01-011	31	1.16	8.40	9.56	9.12	94.8	
P2433-01	72-12013	16	1.13	8.58	9.71	9.36	95.9	
P2434-01	SS050P-SD07-050724-00	32	1.16	8.81	9.97	8.76	86.3	
P2434-02	SS050P-SD08-050724-00	33	1.18	8.80	9.98	7.74	74.5	
P2435-01	SS055P-SD05-050724-00	34	1.19	8.69	9.88	7.77	75.7	
P2435-02	SS055P-SD06-050724-00	35	1.15	8.80	9.95	7.9	76.7	
P2435-03	SS057P-SD04-050724-00	36	1.17	8.81	9.98	8.07	78.3	
P2437-01	359	37	1.12	8.46	9.58	8.75	90.2	
P2437-02	359-E2	38	1.13	8.69	9.82	8.76	87.8	
P2437-03	SOM-44	39	1.00	1.00	2.00	2.00	100.0	debris
P2437-04	COMP-1	40	1.00	1.00	2.00	2.00	100.0	debris
P2437-05	WAREHOUSE-TOTE-OIL	41	1.00	1.00	2.00	2.00	100.0	oil sample
P2438-01	50316	42	1.00	1.00	2.00	2.00	100.0	WIPE Sample
P2438-02	50317	43	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
P2438-03	50358	44	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$