


Sample data

Ident CCV
Sample type Check standard 1
Determination start 2024-05-08 10:17:11 UTC-4
Method anions
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.32 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.440	0.2792	2.285	2.055	Fluoride
2	6.788	0.2822	1.830	3.075	Chloride
3	9.007	0.6168	3.156	3.011	Nitrite
4	11.153	0.4208	1.860	10.347	Bromide
5	12.925	0.5878	2.262	2.571	Nitrate
6	14.277	0.4049	1.094	4.967	Phosphate
7	17.225	1.0658	2.793	15.356	Sulfate


Sample data

Ident CCB
Sample type Sample
Determination start 2024-05-08 10:38:41 UTC-4
Method anions
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.15 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions




Sample data

Ident LB130688BLW
Sample type Sample
Determination start 2024-05-08 11:00:11 UTC-4
Method anions
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.04 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions




Sample data

Ident LB130688BSW
Sample type Check standard 1
Determination start 2024-05-08 11:21:41 UTC-4
Method anions
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 10.98 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.442	0.2790	2.310	2.053	Fluoride
2	6.765	0.2828	1.846	3.081	Chloride
3	8.965	0.6164	3.168	3.010	Nitrite
4	11.083	0.4226	1.883	10.391	Bromide
5	12.847	0.5907	2.291	2.583	Nitrate
6	14.295	0.4165	1.140	5.109	Phosphate
7	17.270	1.0704	2.814	15.421	Sulfate


Sample data

Ident P2403-08
Sample type Sample
Determination start 2024-05-08 11:43:13 UTC-4
Method anions
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.20 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.438	0.0561	0.459	0.431	Fluoride
2	6.753	0.0564	0.370	0.623	Chloride
3	8.940	0.1226	0.629	0.630	Nitrite
4	11.055	0.0866	0.380	2.186	Bromide
5	12.817	0.1197	0.460	0.548	Nitrate
6	14.280	0.0839	0.225	1.052	Phosphate
7	17.255	0.2163	0.563	3.378	Sulfate


Sample data

Ident P2403-07
Sample type Sample
Determination start 2024-05-08 12:04:45 UTC-4
Method anions
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 10.92 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.435	0.0278	0.228	0.224	Fluoride
2	6.745	0.0281	0.183	0.316	Chloride
3	8.930	0.0604	0.311	0.330	Nitrite
4	11.040	0.0429	0.189	1.120	Bromide
5	12.802	0.0592	0.228	0.287	Nitrate
6	14.272	0.0414	0.111	0.534	Phosphate
7	17.253	0.1080	0.281	1.852	Sulfate


Sample data

Ident P2403-07RE 0.25ML
Sample type Sample
Determination start 2024-05-08 12:26:18 UTC-4
Method anions
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.04 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.428	0.0141	0.118	0.125	Fluoride
2	6.727	0.0143	0.095	0.167	Chloride
3	8.902	0.0315	0.161	0.191	Nitrite
4	11.007	0.0224	0.099	0.619	Bromide
5	12.763	0.0312	0.119	0.166	Nitrate
6	14.238	0.0219	0.059	0.297	Phosphate
7	17.217	0.0560	0.147	1.119	Sulfate


Sample data

Ident CCV
Sample type Check standard 1
Determination start 2024-05-08 14:13:54 UTC-4
Method anions
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.26 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.435	0.2775	2.301	2.042	Fluoride
2	6.725	0.2812	1.847	3.064	Chloride
3	8.897	0.6137	3.182	2.997	Nitrite
4	10.983	0.4205	1.894	10.339	Bromide
5	12.733	0.5883	2.305	2.573	Nitrate
6	14.263	0.4126	1.127	5.061	Phosphate
7	17.270	1.0644	2.814	15.337	Sulfate

Sample data

Ident CCB
Sample type Sample
Determination start 2024-05-08 14:42:22 UTC-4
Method anions
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.22 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions