

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 5/16/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:25
In Date: 05/15/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:27
Out Date: 05/16/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB130770

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P2480-01	WC-MH-JJ	21	1.13	8.47	9.6	8.61	88.3	
P2480-02	EPH-MH-JJ	1	1.13	8.58	9.71	8.73	88.6	
P2480-03	VOC-MH-JJ	2	1.12	8.76	9.88	8.84	88.1	
P2481-01	INK	3	1.00	1.00	2.00	2.00	100.0	ink sample
P2482-01	B-01-229	4	1.13	8.54	9.67	8.52	86.5	
P2482-02	SW-01-229	5	1.18	8.44	9.62	8.64	88.4	
P2482-03	SW-02-229	6	1.14	8.68	9.82	8.8	88.2	
P2482-04	SW-03-229	7	1.15	8.42	9.57	8.47	86.9	
P2482-05	B-02-229	8	1.18	8.59	9.77	9.14	92.7	
P2482-06	SW-04-229	9	1.16	8.52	9.68	8.38	84.7	
P2482-07	P-01-229	10	1.18	8.74	9.92	8.84	87.6	
P2486-01	WCS-TPK	11	1.14	8.40	9.54	8.48	87.4	
P2486-02	WCS-TPK-EPH	12	1.18	8.80	9.98	8.85	87.2	
P2486-03	WCS-TPK-VOC	13	1.14	8.83	9.97	8.9	87.9	
P2486-05	WCS-TPL	14	1.16	8.81	9.97	9.14	90.6	
P2486-06	WCS-TPL-EPH	15	1.11	8.50	9.61	8.8	90.5	
P2486-07	WCS-TPL-VOC	16	1.13	8.64	9.77	8.99	91.0	
P2487-01	WC-A-COMP	17	1.17	8.71	9.88	8.86	88.3	
P2487-03	WC-A-GRAB	18	1.19	8.68	9.87	8.53	84.6	
P2488-01	SP-G-COMP	19	1.15	8.82	9.97	9.06	89.7	
P2488-03	SP-G-GRAB	20	1.11	8.46	9.57	8.5	87.4	
P2489-01	BC226734-1-1	22	1.00	1.00	2.00	2.00	100.0	pile
P2489-02	BC226734-1-2	23	1.00	1.00	2.00	2.00	100.0	pile
P2489-03	BC226734-2-1	24	1.00	1.00	2.00	2.00	100.0	pile
P2489-04	BC226734-2-2	25	1.00	1.00	2.00	2.00	100.0	pile
P2490-01	72-11978	26	1.15	8.40	9.55	8.2	83.9	
P2490-02	72-11978-E2	27	1.13	8.70	9.83	8.38	83.3	
P2491-01	VNJ-277	28	1.18	8.65	9.83	9.26	93.4	

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 5/16/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:25
In Date: 05/15/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:27
Out Date: 05/16/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB130770

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P2491-02	VNJ-277-E2	29	1.13	8.74	9.87	9.38	94.4	
P2492-01	BC247601-1-1	30	1.00	1.00	2.00	2.00	100.0	oilc
P2492-02	BC247601-1-2	31	1.00	1.00	2.00	2.00	100.0	oilc
P2492-03	BC247601-2-1	32	1.00	1.00	2.00	2.00	100.0	oilc
P2492-04	BC247601-2-2	33	1.00	1.00	2.00	2.00	100.0	oilc
P2493-01	D243	34	1.00	1.00	2.00	2.00	100.0	debris
P2493-02	D241	35	1.00	1.00	2.00	2.00	100.0	oil sample
P2493-03	D242	36	1.00	1.00	2.00	2.00	100.0	oil sample
P2496-01	CHRT-24376	37	1.16	8.80	9.96	8.72	85.9	
P2496-02	CHRT-24376-E2	38	1.15	8.81	9.96	9.22	91.6	
P2496-03	CHRT-28452	39	1.18	8.56	9.74	7.54	74.3	
P2496-04	CHRT-28452-E2	40	1.19	8.42	9.61	7.97	80.5	
P2497-01	280856	41	1.16	8.50	9.66	9.43	97.3	
P2497-03	RT3823	42	1.11	8.56	9.67	9.41	97.0	
P2497-05	RT3997	43	1.18	8.66	9.84	9.21	92.7	
P2497-07	280646	44	1.15	8.40	9.55	9.42	98.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$