

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 5/17/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 05/16/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 05/17/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB130788

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P2494-01	1	1	1.14	8.56	9.7	9.63	99.2	
P2494-02	2	2	1.15	8.50	9.65	9.54	98.7	
P2494-03	3	3	1.18	8.39	9.57	9.48	98.9	
P2494-04	4	4	1.14	8.59	9.73	9.5	97.3	
P2494-05	5	5	1.18	8.51	9.69	9.51	97.9	
P2494-06	6	6	1.15	8.40	9.55	9.36	97.7	
P2494-07	7	7	1.18	8.44	9.62	9.41	97.5	
P2494-08	8	8	1.17	8.43	9.6	9.5	98.8	
P2494-09	9	9	1.17	8.35	9.52	9.44	99.0	
P2494-10	10	10	1.16	8.38	9.54	9.45	98.9	
P2494-11	11	11	1.19	8.68	9.87	9.72	98.3	
P2494-12	12	12	1.13	8.55	9.68	9.54	98.4	
P2494-13	13	13	1.14	8.52	9.66	9.51	98.2	
P2494-14	14	14	1.14	8.40	9.54	9.43	98.7	
P2494-15	15	15	1.19	8.43	9.62	9.48	98.3	
P2494-16	16	16	1.12	8.76	9.88	9.2	92.2	
P2494-17	17	17	1.13	8.43	9.56	9.54	99.8	
P2494-18	18	18	1.17	8.60	9.77	9.68	99.0	
P2494-20	19	19	1.19	8.65	9.84	9.7	98.4	
P2495-01	RT-3463	20	1.18	8.52	9.7	9.04	92.3	
P2495-02	RT-3463	21	1.17	8.61	9.78	9.26	94.0	
P2498-01	OK-01-051524	22	1.19	8.60	9.79	8.52	85.2	
P2498-02	OK-01-051524-E2	23	1.15	8.79	9.94	9.11	90.6	
P2500-01	B-112-SB01	24	1.18	8.47	9.65	8.62	87.8	
P2500-02	B-112-SB02	25	1.13	8.47	9.6	7.92	80.2	
P2501-01	33351	26	1.00	1.00	2.00	2.00	100.0	oil sample
P2511-01	510-C	27	1.17	8.60	9.77	7.17	69.8	
P2511-02	510-D	28	1.15	8.78	9.93	6.04	55.7	

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 5/17/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 05/16/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 05/17/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB130788

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P2511-03	510-CD	29	1.18	8.79	9.97	7.19	68.4	
P2511-04	665	30	1.14	8.49	9.63	9.38	97.1	
P2512-01	WCS-TPJ	31	1.11	8.78	9.89	8.91	88.8	
P2512-02	WCS-TPJ-EPH	32	1.19	8.64	9.83	8.83	88.4	
P2512-03	WCS-TPJ-VOC	33	1.12	8.77	9.89	8.92	88.9	
P2512-05	WCS-TPI	34	1.17	8.58	9.75	8.98	91.0	
P2512-06	WCS-TPI-EPH	35	1.17	8.62	9.79	8.94	90.1	
P2512-07	WCS-TPI-VOC	36	1.15	8.81	9.96	8.97	88.8	
P2513-01	TP-13	37	1.12	8.86	9.98	9.22	91.4	
P2513-02	TP-13-E2	38	1.14	8.67	9.81	9.06	91.3	
P2513-04	TP-14	39	1.13	8.64	9.77	9.09	92.1	
P2513-05	TP-14-E2	40	1.17	8.57	9.74	8.9	90.2	
P2513-07	TP-15	41	1.15	8.38	9.53	8.91	92.6	
P2513-08	TP-15-E2	42	1.19	8.43	9.62	7.49	74.7	
P2513-10	TP-16	43	1.16	8.81	9.97	9.3	92.4	
P2513-11	TP-16-E2	44	1.13	8.80	9.93	9.12	90.8	
P2513-13	TP-17	45	1.18	8.48	9.66	8.89	90.9	
P2513-14	TP-17-E2	46	1.15	8.37	9.52	8.81	91.5	
P2519-01	B2COMP (0-8)	47	1.12	8.61	9.73	8.42	84.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$