

Analysis Method: 9034
 Parameter: Reactive Sulfide
 Run Number: LB130794

ANALYST: Rubina
 SUPERVISOR REVIEW BY: Iwona
 Constant: 16000
 Normality1: 0.025
 Normality2: 0.025

Reagent/Standard	Lot/Log #
SODIUM THIOSULFATE, 0.025N, 4LITRE	W3105
IODINE SOLUTION .025N 1L	W3050
Starch Solution, 4L	W2977

Seq	Lab ID	True Value (mg/l)	DF	Initial Weight (g)	Final Volume (ml)	T1 (ml)	T2 Initial	T2 Final	T2 Diff. (ml)	T1 - T2 Diff (mL)	Value Corrected With Blank	Result (ppm)	Anal Date	Anal Time
1	PB160942BL		1	5.00	50	2.00	0.00	1.94	1.94	0.06	0.00	0.00	05/16/2024	14:40
2	P2476-02		1	5.06	50	2.00	0.00	1.92	1.92	0.08	0.02	1.58	05/16/2024	14:43
3	P2476-02DUP		1	5.05	50	2.00	0.00	1.92	1.92	0.08	0.02	1.58	05/16/2024	14:45
4	P2477-03		1	5.02	50	2.00	0.00	1.84	1.84	0.16	0.10	7.97	05/16/2024	14:48
5	P2477-04		1	5.03	50	2.00	0.00	1.88	1.88	0.12	0.06	4.77	05/16/2024	14:50
6	P2480-04		1	5.05	50	2.00	0.00	1.90	1.90	0.10	0.04	3.17	05/16/2024	14:52
7	P2486-04		1	5.03	50	2.00	0.00	1.84	1.84	0.16	0.10	7.95	05/16/2024	14:55
8	P2486-08		1	5.02	50	2.00	0.00	1.92	1.92	0.08	0.02	1.59	05/16/2024	14:58
9	P2487-02		1	5.05	50	2.00	0.00	1.86	1.86	0.14	0.08	6.34	05/16/2024	15:00
10	P2488-02		1	5.04	50	2.00	0.00	1.88	1.88	0.12	0.06	4.76	05/16/2024	15:02
11	P2493-01		1	5.02	50	2.00	0.00	1.90	1.90	0.10	0.04	3.19	05/16/2024	15:05
12	P2493-02		1	5.06	50	2.00	0.00	1.86	1.86	0.14	0.08	6.32	05/16/2024	15:08

T1 = Titrant1

T2 = Titrant2

T2 Diff = T2 Final - T2 Initial

Value Corrected With Blank = ((T1 - T2 Diff) - Blank Correction(BL))

Result = ((T1 * Normality1) - ((T1 - Value Corrected With Blank) * Normality2)) * Constant / Initial Volume

Analysis Method: 9034
 Parameter: Reactive Sulfide
 Run Number: LB130794

ANALYST: Rubina
 SUPERVISOR REVIEW BY: Iwona
 Constant: 16000
 Normality1: 0.025
 Normality2: 0.025

Reagent/Standard	Lot/Log #
SODIUM THIOSULFATE, 0.025N, 4LITRE	W3105
IODINE SOLUTION .025N 1L	W3050
Starch Solution, 4L	W2977

Seq	Lab ID	True Value (mg/l)	DF	Initial Weight (g)	Final Volume (ml)	T1 (ml)	T2 Initial	T2 Final	T2 Diff. (ml)	T1 - T2 Diff (mL)	Value Corrected With Blank	Result (ppm)	Anal Date	Anal Time
13	P2493-03		1	5.07	50	2.00	0.00	1.90	1.90	0.10	0.04	3.16	05/16/2024	15:10
14	P2497-02		1	5.05	50	2.00	0.00	1.86	1.86	0.14	0.08	6.34	05/16/2024	15:13
15	P2497-04		1	5.05	50	2.00	0.00	1.90	1.90	0.10	0.04	3.17	05/16/2024	15:15
16	P2497-06		1	5.05	50	2.00	0.00	1.86	1.86	0.14	0.08	6.34	05/16/2024	15:18
17	P2497-08		1	5.02	50	2.00	0.00	1.92	1.92	0.08	0.02	1.59	05/16/2024	15:20

T1 = Titrant1

T2 = Titrant2

T2 Diff = T2 Final - T2 Initial

Value Corrected With Blank = ((T1 - T2 Diff) - Blank Correction(BL))

Result = ((T1 * Normality1) - ((T1 - Value Corrected With Blank) * Normality2)) * Constant / Initial Volume