

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 5/29/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 05/28/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 05/29/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB130943

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P2608-01	522	1	1.00	1.00	2.00	2.00	100.0	wipe sample
P2608-02	523	2	1.00	1.00	2.00	2.00	100.0	wipe sample
P2608-03	524	3	1.00	1.00	2.00	2.00	100.0	wipe sample
P2610-01	HR-01-05242024	32	1.19	8.51	9.7	8.94	91.1	
P2610-02	HR-01-05242024	33	1.15	8.81	9.96	9.14	90.7	
P2610-03	HR-02-05242024	34	1.16	8.81	9.97	9.4	93.5	
P2610-04	HR-02-05242024	35	1.17	8.79	9.96	9.26	92.0	
P2614-01	52324	4	1.00	1.00	2.00	2.00	100.0	wipe sample
P2622-01	S01-01A	5	1.15	8.80	9.95	8.72	86.0	
P2622-02	S01-01B	6	1.16	8.58	9.74	8.23	82.4	
P2622-03	S01-01C	7	1.19	8.44	9.63	7.87	79.1	
P2622-04	S01-02A	8	1.17	8.44	9.61	8.47	86.5	
P2622-05	S01-02B	9	1.12	8.50	9.62	8.48	86.6	
P2622-06	S01-02C	10	1.12	8.60	9.72	8.46	85.3	
P2622-07	S01-03A	11	1.14	8.83	9.97	8.81	86.9	
P2622-08	S01-03B	12	1.15	8.37	9.52	8.33	85.8	
P2622-09	S01-03C	13	1.16	8.52	9.68	8.25	83.2	
P2622-10	S02-01A	14	1.14	8.82	9.96	8.01	77.9	
P2622-11	S02-01B	15	1.14	8.84	9.98	8.33	81.3	
P2622-12	S02-01C	16	1.15	8.79	9.94	8.44	82.9	
P2622-13	S02-02A	17	1.13	8.80	9.93	8.07	78.9	
P2622-14	S02-02B	18	1.17	8.74	9.91	8.27	81.2	
P2622-15	S02-02C	19	1.16	8.59	9.75	8.23	82.3	
P2622-16	S02-03A	20	1.15	8.40	9.55	8.07	82.4	
P2622-17	S01-03B	21	1.12	8.85	9.97	8.31	81.2	
P2622-18	S01-03C	22	1.14	8.44	9.58	7.41	74.3	
P2622-19	S02-04A	23	1.16	8.68	9.84	8.45	84.0	
P2622-20	S02-04B	24	1.11	8.86	9.97	8.44	82.7	

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 5/29/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 05/28/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 05/29/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB130943

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P2622-21	S02-04C	25	1.15	8.83	9.98	8.99	88.8	
P2622-22	S02-05A	26	1.16	8.79	9.95	8.4	82.4	
P2622-23	S02-05B	27	1.16	8.60	9.76	8.11	80.8	
P2622-24	S02-05C	28	1.13	8.58	9.71	8.15	81.8	
P2622-25	S02-06A	29	1.19	8.45	9.64	7.98	80.4	
P2622-26	S02-06B	30	1.18	8.79	9.97	9.14	90.6	
P2622-27	S02-06C	31	1.2	8.79	9.99	8.99	88.6	
P2629-01	SS057P-SD05-0522424-00	37	1.15	8.76	9.91	5.4	48.5	
P2629-02	SS029P-SD05-0522424-00	38	1.13	8.74	9.87	2.74	18.4	
P2629-03	FT011P-SD06-0522424-00	39	1.14	8.82	9.96	2.82	19.0	
P2629-04	FT011P-SD05-0522424-00	40	1.13	8.71	9.84	2.9	20.3	
P2629-05	AT024P-SD05-0522424-00	41	1.15	8.83	9.98	6.7	62.9	
P2630-01	FT0133P-SD05-052424-00	42	1.16	8.83	9.99	5.01	43.6	
P2631-01	34795	36	1.00	1.00	2.00	2.00	100.0	wipe sample
P2637-01	WCS-TPT	43	1.18	8.50	9.68	8.65	87.9	
P2637-02	WCS-TPT-EPH	44	1.17	8.57	9.74	8.68	87.6	
P2637-03	WCS-TPT-VOC	45	1.15	8.82	9.97	8.79	86.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$